

MOBILITY





TRANSPORT CHAIR WITH FLIP-UP DESK ARMS ITEM # 353



THE FEATURES

- 1 *One-Grip* locking handbrakes
- 2 Lightweight aluminum frame
- 3 Flip-up desk arms
- 4 Non-skid foot plates & heel loops
- 5 Foldable
- 6 8" front wheels & 12" locking rear wheels
- 7 Anti-tippers

Available in:

-  353B - Blue
-  353R - Red
-  353BK - Black
-  353P - Purple



NEW
One-Grip Handbrakes
(see back page for details)

We love the 353 Transport Chair because it's lightweight, yet versatile and all-terrain, with *One-Grip* locking handbrakes, flip-up desk arms, rugged rear 12" wheels and anti-tippers.

NOVA Retail Price

\$409.95 ← *Budget Friendly with lower MAP and MSRP!*

Item #	Weight Capacity	Overall Weight (w/out footrests)	Width btw Arms	Seat Height	Seat Dimensions	Overall Dimensions	Overall Folded Dimensions
353	300 lb	27 lb	19"	19"	18.5"w x 15.5"d	24"w x 30"d x 39.5"h	11.5"w x 29"d x 30.5"h

Los Angeles - Chicago - Atlanta

©2024 NOVA



CARE • QUALITY • INNOVATION



THE NOVA PROMISE

novajoy.com

CONTACT US!
We're here to help

 info@novajoy.com
 1 800 557 6682

One-Grip Braking System

It's as Easy as 1, 2, 3!



- 1 Squeeze lever partially (about fifty percent engaged) to **control the speed**
- 2 Squeeze lever further to **stop completely**.
Brake is locked when you hear a click
- 3 To **release brake**, engage the Release Tab

