



Yes, you can.®

TheraPure™ 3752G & 3802G Recessed Side-Entry Tubs



3802G Shown

NEW

Featuring SANIJET®
Pipeless™ technology

For environments requiring a more “residential” looking tub, the 3752G and 3802G tubs provide the easy access of a side-opening door in a design suitable for a variety of installations, including being recessed in a standard tub/shower alcove.

Featuring 4 easily cleaned whirlpool jets powered by SANIJET Pipeless technology, these Invacare tubs combine resident ease of use with caregiver ease of disinfection and peace of mind.

TheraPure™ 3752G & 3802G Recessed Side-Entry Tubs



Integrated console (3752G shown)

Provides the resident with easy access to the auto-fill/ auto shut-off and anti-scald mixing valve, thermometer and whirlpool jet controls.



Innovation at hand

The shower wand features an integrated flow control lever to deliver the desired water pressure, while the spray setting can be adjusted by rotating the showerhead.



Indulge residents with a soothing spa experience

An Invacare bath is quieter, cleaner and stays warm longer than air bubble and piped systems for the ultimate in relaxation.



Pipeless™ Purity

Without pipes, pumps or air channels that can provide a breeding ground for bacteria, **four jets allow caregivers to quickly remove up to 99.9%* of bacteria after each bath.**

Specifications

	3802G	3752G
Outside length:	59.5"	59.5"
Outside width:	31.5"	31.5"
Inside length:	50"	50"
Inside width:	22.5"	22.5"
Rim height:	41.5"	41.5"
Seat height:	19"	19"
Fill rate:	17 GPM	17 GPM
Hydromassage:	72-80 GPM flow rate	72-80 GPM flow rate

Features

- All wetted components are easily accessible for quick and thorough disinfection
- Pipeless technology conserves water, chemicals and time
- Heavy-duty hinges suspended on a spring-assisted raising and lowering mechanism support the door
- Door is easily locked and unlocked with an easy-to-use padded handle
- Fast-fill mixing anti-scald valve dispenses 17 gallons per minute
- Fiberglass end panels are available for corner and one wall installations
- Door seal is guaranteed for life

*SANIJET web site, retrieved from http://www.sanijet.com/WhyPipeless_Genuine.aspx

Invacare Continuing Care, Inc.

www.invacare-cc.com

info_icc@invacare.com

USA

One Invacare Way
Elyria, Ohio 44035
USA
800-668-2337

Canada

994 Hargrieve Road
London, Ontario
N6E 1P5 Canada
800-668-2337

Trademarks are identified by the symbol™ and ®. All trademarks are owned or licensed to Invacare Corporation or its subsidiaries unless otherwise noted. SANIJET® and Pipeless™ are trademarks of Maple Industries, LLC, d.b.a. SANIJET.

© 2010 Invacare Corporation. Form 04.17.5/10



Yes, you can: