

Protekt™ Aire 4000

Low Air Loss/Alternating Pressure Mattress System



 **Proactive**
medical products

Medicare
Code
EO277

270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com

The Protekt™ Aire 4000 is the 1st defense for facilities to treat pressure ulcers stages I-IV. The high volume user friendly pump combines low air loss and alternating pressure in one unit. The air flow and pressure range can be customized according to each individual patients weight. The “static” mode provides firm support for easy patient transfer. The mattress is comprised of (20) 8” air bladders that allow for deep immersion and excellent pressure redistribution. The top cover is quilted for added patient comfort.

FEATURES AND BENEFITS

MATTRESS & COVER

- 18 laser air holes to keep patient cool and dry.
- Removable, water resistant, non-shear, anti-microbial nylon cover is quilted for extra comfort.
- CPR valve provides rapid deflation of mattress.
- Non-skid cover prevents the mattress from slipping.
- Dimensions: 35”x80”x8”.
- Weight: 9 lbs.
- Alternating cell depth: 8”.

PROTEKT™ AIRE 4000 PUMP

- High flow pump with 8 liters of air per minute.
- Adjustable comfort setting allows pressure to be customized according to each patients therapeutic requirements.
- Low pressure alarm indicates when the pumps airflow may be compromised.
- Low pressure alarm mute button.
- Static function for safe and easy patient transfer.
- UL/CSA approved power cord.
- Weight Capacity: 350 lbs.
- Dimension: 11”x8”x4”.
- Weight: 6.2 lbs.
- 18 month warranty.

80040	Protekt™ Aire 4000	8” Alternating/Low Air Loss Mattress System
80041	Protekt™ Aire 4000 Pump	Pump only for Protekt™ Aire 4000
80042	Protekt™ Aire 4000 Mattress	Mattress only for Protekt™ Aire 4000
80043	Protekt™ Aire 4000 Cover	Cover only for Protekt™ Aire 4000
80044	Protekt™ Aire 4000 Bladder	Individual Bladder for Protekt™ Aire 4000



270 Washington Street Mount Vernon, NY 10553

Tel: 855-BE-PROACTIVE (855-237-7622) • Fax: 845-205-6005

www.proactivemedical.com • sales@proactivemedical.com