



**2 SIZES  
AVAILABLE**



Facebook.com/  
LeggeroUSA



14



12



# TRAK™

- ✔ Free Sunshade & Storage Bag
- ✔ Lightweight (29 lbs/13 kg)
- ✔ 30° Hip Angle Adjustment
- ✔ 45° Tilting Frame
- ✔ Configurable Seating
- ✔ Folds Compactly
- ✔ Rear Suspension



## SPECIFICATIONS & ADJUSTMENT RANGE

|                            | TRAK™ 12                             | TRAK™ 14                             |
|----------------------------|--------------------------------------|--------------------------------------|
| Weight Limit               | 80 lbs / 36 kg                       | 80 lbs / 36 kg                       |
| Seat Depth                 | 6.5-15 in / 17-38 cm                 | 6.5-15 in / 17-38 cm                 |
| Seat Width                 | 6-12 in / 15-30 cm                   | 8-14 in / 20-36 cm                   |
| Back Height                | 12-19 in / 30-48 cm                  | 13-19 in / 33-48 cm                  |
| Seat to Footplate          | 4-17 in / 10-43 cm                   | 4-17 in / 10-43 cm                   |
| Tilt Range                 | 45°                                  | 45°                                  |
| Hip Angle Positions        | 90°, 105°, 120°                      | 90°, 105°, 120°                      |
| Maximum Safe Slope         | 5°                                   | 5°                                   |
| Overall Width              | 22.5 in / 57 cm                      | 24.5 in / 62 cm                      |
| Push Handle Height         | 40-50 in / 101-127 cm                | 40-50 in / 101-127 cm                |
| Overall Length             | 35 in / 89 cm                        | 35 in / 89 cm                        |
| Caster Diameter            | 7 in / 18 cm                         | 7 in / 18 cm                         |
| Rear Wheel Diameter        | 16 in / 41 cm                        | 16 in / 41 cm                        |
| Weight of Chair            | 29 lbs / 13 kg                       | 29 lbs / 13 kg                       |
| Folded without Rear Wheels | 35 x 17 x 22.5 in<br>89 x 43 x 57 cm | 35 x 17 x 24.5 in<br>89 x 43 x 62 cm |
| Folded with Wheels         | 36 x 21 x 22.5 in<br>92 x 51 x 57 cm | 36 x 21 x 24.5 in<br>92 x 51 x 62 cm |

30°  
Hip  
Angle  
Range

45° Tilt Range

Easy folding with  
seating components  
in place.

## FEATURES & OPTIONS

- Easy to operate foot tilt lever
- Accessories include Oxygen Tank Holder, IV Pole, LTV Vent Bracket & Ventilator/Battery Tray
- Numerous Positioning Components
- Height Adjustable Push Handle
- Foot Activated Wheel Locks
- Polyurethane Wheels & Casters
- Height Adjustable Footplate
- Rear Anti-tips

## RESPOND™ SEATING

- Contoured Back Cushion
- Contoured Seat Cushion
- Adjustable Lateral Thoracic Supports
- Adjustable Hip Guides

