

Caster and Wheel Options

Select one wheel / caster / STF combination from the options below. Then select related options. STF heights are approximate and measured with the seat rails level.

Wheel Size	Caster Size	Seat To Floor Height Options													
↓	↓	Item number builds from selection: 3101-WWW-CC-STF where W = Wheel Selection, C = Caster Selection, STF = Seat to Floor Selection													
		13.5"/34cm	14"/35cm	14.5"/37cm	15"/38cm	15.5"/39cm	16"/40cm	16.5"/42cm	17"/43cm	17.5"/44cm	18"/46cm	18.5"/47cm	19"/48cm	19.5"/49cm	20"/51cm

12" or 16" Poly Mag Leggero Smooth Roll Wheels and Casters [n/c]

12"/30cm 3101-2010	5"/13cm -05	13.5"/34cm -135 □ * ‡	14"/35cm -140 □		15"/38cm -150 □	15.5"/39cm -155 □	16"/40cm -160 □	16.5"/42cm -165 □	17"/43cm -170 □						
	6"/15cm -06				15"/38cm -150 □	15.5"/39cm -155 □	16"/40cm -160 □	16.5"/42cm -165 □	17"/43cm -170 □						
	7"/18cm -07					15.5"/39cm -155 □	16"/40cm -160 □	16.5"/42cm -165 □	17"/43cm -170 □						
16"/40cm 3101-2011	5"/13cm -05		14"/35cm -140 □ ‡	14.5"/37cm -145 □ ‡	15"/38cm -150 □ ‡	15.5"/39cm -155 □ ‡	16"/40cm -160 □		17"/43cm -170 □	17.5"/44cm -175 □	18"/46cm -180 □				
	6"/15cm -06				15"/38cm -150 □ ‡	15.5"/39cm -155 □ ‡	16"/40cm -160 □		17"/43cm -170 □	17.5"/44cm -175 □	18"/46cm -180 □	18.5"/47cm -185 □			
	7"/18cm -07					15.5"/39cm -155 □ ‡	16"/40cm -160 □		17"/43cm -170 □	17.5"/44cm -175 □	18"/46cm -180 □	18.5"/47cm -185 □	19"/48cm -190 □		

12" or 16" Hub Lock (HL) Poly Mag Leggero Smooth Roll Wheels and Casters [n/c]

HL12"/30cm 3101-2012	5"/13cm -05				15"/38cm -150 □ ‡	15.5"/39cm -155 □ ‡	16"/40cm -160 □ ‡	16.5"/42cm -165 □	17"/43cm -170 □						
	6"/15cm -06				15"/38cm -150 □ ‡	15.5"/39cm -155 □ ‡	16"/40cm -160 □ ‡	16.5"/42cm -165 □	17"/43cm -170 □						
	7"/18cm -07					15.5"/39cm -155 □ ‡	16"/40cm -160 □ ‡	16.5"/42cm -165 □	17"/43cm -170 □						
HL16"/40cm 3101-2013	5"/13cm -05								17"/43cm -170 □ ‡	17.5"/44cm -175 □ ‡	18"/46cm -180 □ ‡				
	6"/15cm -06								17"/43cm -170 □ ‡	17.5"/44cm -175 □ ‡	18"/46cm -180 □ ‡	18.5"/47cm -185 □			
	7"/18cm -07								17"/43cm -170 □ ‡	17.5"/44cm -175 □ ‡	18"/46cm -180 □ ‡	18.5"/47cm -185 □	19"/48cm -190 □		

16" Drum Brake (DB) Compatible Poly Mag Wheels and Leggero Smooth Roll Casters [n/c]

DB16"/40cm 3101-2014	5"/13cm -05		14"/35cm -140 □ ‡	14.5"/37cm -145 □ ‡	15"/38cm -150 □ ‡	15.5"/39cm -155 □ ‡	16"/40cm -160 □		17"/43cm -170 □	17.5"/44cm -175 □	18"/46cm -180 □				
	6"/15cm -06				15"/38cm -150 □ ‡	15.5"/39cm -155 □ ‡	16"/40cm -160 □		17"/43cm -170 □	17.5"/44cm -175 □	18"/46cm -180 □	18.5"/47cm -185 □			
	7"/18cm -07					15.5"/39cm -155 □ ‡	16"/40cm -160 □		17"/43cm -170 □	17.5"/44cm -175 □	18"/46cm -180 □	18.5"/47cm -185 □	19"/48cm -190 □		

20" Mag Wheels and Leggero Smooth Roll Casters [n/c]

20"/51cm 3101-2015	5"/13cm -05					15.5"/39cm -155 □ ‡	16"/40cm -160 □ ‡	16.5"/42cm -165 □ ‡	17"/43cm -170 □ ‡	17.5"/44cm -175 □					
	6"/15cm -06					15.5"/39cm -155 □ ‡	16"/40cm -160 □ ‡	16.5"/42cm -165 □ ‡	17"/43cm -170 □ ‡	17.5"/44cm -175 □		18.5"/47cm -185 □			
	7"/18cm -07					15.5"/39cm -155 □ ‡	16"/40cm -160 □ ‡	16.5"/42cm -165 □ ‡	17"/43cm -170 □ ‡	17.5"/44cm -175 □		18.5"/47cm -185 □	19"/48cm -190 □	19.5"/49cm -195 □	20"/51cm -200 □

Tire Options for 20" Mag Wheels

☐	3101-2027	20" Poly Tire	n/c
☐	3101-2017	E2213 20" Rubber Tire with Airless Insert	\$90.00

Handrim Options for 20" Mag Wheels

☐	3101-2016	20" Anodized Aluminum Handrim	\$119.00
--------------------------	-----------	-------------------------------	----------

22" Spoke Wheels with Anodized Aluminum Handrims [+ \$204] and Leggero Smooth Roll Casters

22"/56cm 3101-2018	5"/13cm -05							16.5"/42cm -165 □ ‡	17"/43cm -170 □ ‡	17.5"/44cm -175 □ ‡	18"/46cm -180 □ ‡				
	6"/15cm -06							16.5"/42cm -165 □ ‡	17"/43cm -170 □ ‡	17.5"/44cm -175 □ ‡	18"/46cm -180 □ ‡	18.5"/47cm -185 □			
	7"/18cm -07							16.5"/42cm -165 □ ‡	17"/43cm -170 □ ‡	17.5"/44cm -175 □ ‡	18"/46cm -180 □ ‡	18.5"/47cm -185 □	19.5"/49cm -195 □	20"/51cm -200 □	

Tire Options for 22" Spoke Wheels

☐	3101-2023	22" Poly Tire	n/c
☐	3101-2024	22" Pneumatic Tire	n/c
☐	3101-2019	E2213 22" Airless Insert	+ \$90.00

Handrim Options for 22" Spoke Wheels

☐	3101-2099	Omit Anodized Alum Handrim	n/c
--------------------------	-----------	----------------------------	-----

24" Spoke Wheels with Anodized Aluminum Handrims [+ \$204] and Leggero Smooth Roll Casters

24"/61cm 3101-2021	5"/13cm -05									17.5"/44cm -175 □ ‡	18"/46cm -180 □ ‡				
	6"/15cm -06									17.5"/44cm -175 □ ‡	18"/46cm -180 □ ‡	18.5"/47cm -185 □ ‡			
	7"/18cm -07									17.5"/44cm -175 □ ‡	18"/46cm -180 □ ‡	18.5"/47cm -185 □ ‡	19"/48cm -190 □ ‡	19.5"/49cm -195 □	

Tire Options for 24" Spoke Wheels

☐	3101-2025	24" Poly Tire	n/c
☐	3101-2026	24" Pneumatic Tire	n/c
☐	3101-2022	E2213 24" Airless Insert	+ \$90.00

Handrim Options for 24" Spoke Wheels

☐	3101-2199	Omit Anodized Alum Handrim	n/c
--------------------------	-----------	----------------------------	-----

* Anterior tilt is not available with 13.5"/34cm STF option

‡ Indicates wheel position is not compatible with the extended rear seat rail option (3101-3110)

Customer Service: 1-844-503-5437 | Fax: 1-406-551-9370 | Email: sales@leggero.us

Prices and specifications are subject to change without notice. ©2020 Leggero LLC.