



ENZO



From the first concept to the final system, every aspect of the Leggero Enzo was designed to maximize in aesthetics, performance, and safety.

The Enzo utilizes lightweight, weld-free cast aluminum, and precise machining to provide unparalleled rigidity in folding and rigid configurations.

WEIGHT CAPACITY 185lbs (84kg)

- SEAT DEPTH 11-20" (28-51cm)
- SEAT WIDTH 12-18" (30-46cm)
- SEAT TO FLOOR 13.5-20" (34-51cm)
- TILT ANGLE 8-14" (20-35cm)
- HIP ANGLE -5° to 45°
- FOLDED WITH WHEELS 90° to 110°

The Leggero Enzo is available in a range of metallic colors.



Safe



Customizable



Adjustable



Compact



www.leggero.us

☎ 1-844-503-5437

✉ sales@leggero.us

Stealth EPiC Seating
is exclusively available on the Leggero Enzo
directly out of the box!

What is EPiC Seating?

EPiC Seating™ is a Postural Management System that introduces movement into a wheelchair seating environment.

The system provides the ability to change positions of function while maintaining control of the pelvis, mitigating the experience of shearing.

EPiC Seating is now coded K0108!



Why combine Leggero's Enzo with Stealth's EPiC Seating?

- Does not void the Enzo's factory warranty
- Allows a client to reposition independently for comfort or tool-free adjustment by a caregiver
- Provides multiple resistance levels to accommodate a wide range of clients



Leggero