



REACH™

Crafted from premium cast aluminum, the Reach effortlessly combines durability with a featherlight feel.

With its sleek, modern design, the Leggero Reach is the epitome of style, functionality, and adjustability.


MAX WEIGHT CAPACITY 150lbs (68kg)

- AVAILABLE IN 3 SIZES
- LIGHTWEIGHT
- ADJUSTABLE FIXED TILT
- KNEE ANGLE ADJUSTABILITY
- QUICK HIP ANGLE ADJUSTABILITY
- FREE SUNSHADE & STORAGE BAG

The Leggero Reach is available in 4 bold fabric colors.




Safe



Adjustable



Customizable



Compact




www.leggero.us

☎ 1-844-503-5437

✉ sales@leggero.us

What is NEW on the Leggero Reach?



Magenta Fabric Option



Shadow Frame Finish



Perforated Leather Handle



Black Wheels

	REACH ¹²	REACH ¹⁴	REACH ¹⁶
WEIGHT LIMIT	80 lbs (36 kg)	110 lbs (50 kg)	150 lbs (68 kg)
SEAT WIDTH	8-12 in (20-30 cm)	10-14 in (25-35 cm)	12-16 in (30-40 cm)
SEAT DEPTH	8-13 in (20-33 cm)	11-16 in (28-40 cm)	13-18 in (33-46 cm)
BACK HEIGHT	11-18 in (28-46 cm)	13-21 in (33-53 cm)	16-24 in (40-61 cm)
SEAT TO FOOTPLATE	8-13 in (20-33 cm)	11-16 in (28-40 cm)	13-18 in (33-46 cm)
ADJUSTABLE FIXED TILT	10°, 20°, or 30°	10°, 20°, or 30°	10°, 20°, or 30°
SEAT TO BACK ANGLE RANGE (10° FIXED TILT)	85°, 90°, 95°, 100°, 105°, 110°, 115°	85°, 90°, 95°, 100°, 105°, 110°, 115°	85°, 90°, 95°, 100°, 105°, 110°, 115°
SEAT TO BACK ANGLE RANGE (20° FIXED TILT)	85°, 90°, 95°, 100°, 105°	85°, 90°, 95°, 100°, 105°	85°, 90°, 95°, 100°, 105°
SEAT TO BACK ANGLE RANGE (30° FIXED TILT)	85°, 90°, 95°	85°, 90°, 95°	85°, 90°, 95°
CHAIR WEIGHT	27 lbs (12 kg)	29 lbs (13 kg)	31 lbs (14 kg)
MAXIMUM SAFE SLOPE	5°	5°	5°
REAR WHEEL DIAMETER	12 in (30 cm)	12 in (30 cm)	12 in (30 cm)
FRONT CASTER DIAMETER	7 in (18 cm)	7 in (18 cm)	7 in (18 cm)
CHAIR WIDTH	21.5 in (55 cm)	23.5 in (60 cm)	25.5 in (65 cm)
FOLDED DIMENSIONS	31 x 21.5 x 18 in (79 x 55 x 46 cm)	31 x 23.5 x 18 in (79 x 60 x 46 cm)	33 x 25.5 x 18 in (84 x 65 x 46 cm)

