

Live the Longer
Life You Love®

7000 MILLENNIAL CRUTCH

Finally, a crutch that works with you, not against you



Features & Benefits

- **Spring Assist/Shock Absorbing System:**
 - Absorbs the impact when walking
 - Returns the energy positively in forward motion
 - Provides comfort even in the uneven places
- **Anatomical handle**
 - Keeps wrist and hand in their natural position
 - Reduces damage to the carpal tunnel in the wrist & hand
- **Other Features:**
 - The sleek design and extra cushion eliminates under arm jamming
 - Convenient Folding system allows you to fold and store
 - Comes in a variety of different colors

Specifications

#7001 Junior Pair: for people 4' 7" to 5' 7", adjusts from 42" to 50" in height, weight: 2 lb per crutch

#7000 Standard Pair: for people 5' 7" to 6' 7" adjusts from 49" to 62" in height, weight: 2.5 lb per crutch

Cradle dimensions: 11" long, 1.5" diameter

Weight Capacity: 400 lbs

