

Live the Longer
Life You Love[®]

8000 EZ ADJUSTABLE BED RAIL

Extendable Bed Rail Pivots Away

EXTENDS 26" - 42" →

↻
PIVOTS DOWN 180°



EXTENDS 26" - 42"



EASILY GET OUT OF BED



EASILY PIVOTS DOWN

Features & Benefits

- **Bed Rail/Hand Rail** - Adjustable length allows for use as both a bed rail to prevent falling out of bed or as a hand rail to assist in getting out of bed.
- **Adjustable** - Bed rail easily adjusts from 26" to 34" to 42" after it has been placed on the bed.
- **Collapsible** - Rail folds down to side of bed to allow user to get out of bed or provide space for making bed.
- **Reversible** - Can be placed on either side of the bed
- **Dual Safety Strap** - Secures bed rail to any home or hospital bed. Accommodates up to King size beds.

Specifications

Width of Rail: Locks in place at 26", 34" & 42" after installed on bed
 Height of Rail: 23" high from the base of mattress
 Gap Between Rails: 3"
 Bed Type: Attaches to any home or hospital bed with included Safety Strap
 Weight of Product: 14 lbs
 Weight Capacity: 400 lbs
 Assembly: Installs in minutes with 4 bolts and Allen Wrench (Included)
 Package Dimensions: 28" x 21" x 1.5"
 Case Dimensions: 4/case, 28" x 21" x 6" 60 lbs

How Do We Compare

	EZ Adjust Bed Rail	Other Bed Rails
Adjustable Length Before and After Installation	✓	✗
Pivots Out of the Way	✓	✗
Versatility (Bed Rail and Bed Handle)	✓	✗
Accommodates All Bed Types	✓	✗
Durable Powder Coat Finish	✓	✓
Prevents Falling from the Bed	✓	✓

Easy Assembly (No Additional Tools Required)



Attach L Shaped Tubes



Slide on Strap



Attach Rear Crossbar



Place On Bed