

Bath Chair

Foldable

9406

We love the all new Foldable Bath Chair because it's so easy to fold and store when not in use. It features a rust-resistant aluminum frame, adjustable seat height, comfortable foam seat and back padding, and slip-resistant rubber tips for added stability. Best of all, it's super light-weight at only 8.75 lb!



MSRP: \$109.95
UMAP: \$99.95

Item #	Weight Cap.	Weight	Seat Height	Seat Dim.	Folded Dim.	Overall Dim.
9406	275 lb	8.75 lb	14" - 18.5"	16.25"w x 15"d	40.75"l x 21.25"w x 6"d	21.25"w x 26.5"d x 35.25"h