



# Uplift Walker

Get up and on the move more easily



The Carex Uplift Walker offers assistance while transitioning to a standing position and walking . The walker is height adjustable from 30-37 inches. The single button design allows the walker to fold for easy transport and storage. Features soft hand grips and a matte silver finish with blue crossbar.

- Assists in standing and walking
- Easily folds for storage and transport
- Height adjustable 28-33 inches
- Tri-lingual hang tag

**Product #:** FGA849-C0  
**UPC:** 023601018496  
**Product:** 30-37”h x 20”w x 19.5”d  
**Weight Capacity:** 250 lbs.  
**MSRP:** \$71.00

**Carex Health Brands**  
customerservice@carex.com  
Phone 800.526.8051  
www.carex.com

