

## AL225 Axis II Lite Powered Inside Lift



### FEATURES

- **A light duty and economical, yet effective solution**
- **Powered rotation and lifting**
- **Adjustable base with telescoping boom strap**
- **Weather resistant hand control**
- **Curved, offset post allows mobility device to be delivered to passenger's or driver's side**

### SPECIFICATIONS

Lifting Capacity:	250 lbs (113 kgs)
Operation:	Power lifting & rotation
Adj. Boom Height:	30" - 40"
Adj. Boom Length:	22" - 34"
Adj. Boom Pitch:	0 - 30 degrees
Warranty:	3-year transferable

### DESCRIPTION

Harmar's AL225 is an innovative and feature filled boom style lift. Without sacrificing quality and safety, the AL225 is value engineered for use with lighter power chairs and scooters. Power lifting combined with power rotation make loading and unloading your mobility device an effortless task.