



P.O. number:*

Company Name:*

Date:

Contact name:* Phone:*

Email:

Marked for:

Ship to:*

Name:

Address:


City: State:

Zip code:

Phone:

* Required

Taxes not included

MODEL										USD		
1	APWC11	Ultralight Aluminum APEX Rigid Wheelchair (250 lb weight limit)						Tripple-butted aluminum frame / Aluminum back canes and back rigidizer / Aluminum seat rigidizer / Ultralight carbon fiber camber tube and rear frame rigidizer		\$2 595	☐	
2	APWC12	Apex Carbon fiber frame upgrade						Ultralight carbon fiber frame, back canes, back rigidizer, seat rigidizer, camber tube and rear frame rigidizer.		\$1 700	☐	
SEAT WIDTH												
3	APSW	Seat width:	12" ☐	13" ☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐	NC	☐	
4	APSW	Seat width:	19" ☐	20" ☐							\$495	☐
SEAT DEPTH												
5	APSD	Seat depth:	12" ☐	13" ☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐	19" ☐	NC	☐
6	APSD	Seat depth:	20" ☐								\$475	☐
FRAME LENGTH												
7	APSP01	Standard frame length configuration - The frame bend starts 0.5" in front of the seat sling								STD	☐	
8	APSP02	Frame length 1" longer than standard configuration								NC	☐	
9	APSP03	Frame length 2" longer than standard configuration								NC	☐	
10	APSP04	Frame length 1" shorter than standard configuration								NC	☐	
11	APSP05	Frame length 2" shorter than standard configuration								NC	☐	

FRONT Seat-to-floor height

Caster sizes	Height available +/- 1/4"	Specify desired height: (inches)
3" Caster	14" to 21"	
4" Caster	14.5" to 21"	
5" Caster	15.5" to 21"	
6" Caster	16.5" to 21" (N/A with 90° front frame)	



REAR WHEELS									
51	APRW01	Mag (aluminum or plastic handrim only)	20" ☐	22" ☐	24" ☐			\$50	☐
52	APRW02	Newton One - spoke wheel	20" ☐	22" ☐	24" ☐		New	STD	☐
53	APRW03	Newton Gravity - ultralight wheel	20" ☐	22" ☐	24" ☐	25" ☐	26" ☐	\$300	☐
54	APRW04	Spinergy Spox		22" ☐	24" ☐	25" ☐	26" ☐	\$795	☐
55		Spinergy Spox spokes colors*	Blk ☐	STD	Wht ☐	Blue ☐	Red ☐	Yel ☐	NC
56	APRW05	Spinergy LX		22" ☐	24" ☐	25" ☐	26" ☐	\$995	☐
57		Spinergy LX spokes colors *	Blk ☐	STD	Wht ☐	Blue ☐	Red ☐	Yel ☐	NC

* Black spokes standard lead-time / Colour spokes up to 2 weeks extra lead-time

TIRES						Available			
58	APTIO1	Soft Urethane 1 1/4 "				20", 22" 24"		NC	☐
59	APTIO2	Soft Urethane knobby 1 3/8 "				22", 24", 25", 26"	New	\$95	☐
60	APTIO4	Soft Urethane - Shox 1"				22", 24", 25", 26"	New	\$115	☐
61	APTIO5	Pneumatic 1 3/8 "				20", 22", 24", 25", 26"		STD	☐
62	APTIO8	Pneumatic 1 3/8" with Airless insert				22", 24"	New	\$80	☐
63	APTIO6	Pneumatic HP - Schwalbe Right Run 1" *				20", 22", 24", 25", 26"		\$95	☐
64	APTIO7	Pneumatic HP - Schwalbe Marathon Plus 1" *				24", 25", 26"		\$125	☐

* Not available with Mag Wheels

HANDRIMS						Assembly			
65	APHR01	Aluminum anodized				Regular STD ☐ Narrow* ☐		STD	☐
66	APHR07	Newton Air Grip * Black						\$95	☐
67	APHR	Newton Air Grip * (08) Sunk. Orange ☐ (09) Burgundy ☐ (10) Blue Sapphire ☐ (11) Monster Green ☐ (12) Charcoal ☐ (13) Fuchsia ☐ (14) Steel Blue ☐ (15) Ferrari Red ☐ (16) Acid Green ☐ (17) White ☐ (18) Sky Blue ☐						\$240	☐
68	APHR02	Plastic coated						\$95	☐
69	APHR03	Q-Grip High friction Neoprene coated*						\$295	☐
70	APHR05	Surge LT* ¹ (oval aluminum handrim w/ a rubber high friction strip)						\$380	☐
71	APHR04	Natural Fit* ² (oval aluminum handrim w/thumb piece)						\$375	☐
72	APHR06	Bbraver* ³ (ultra-performance aluminum handrim, adjustable)				Orange ☐ White ☐ Black ☐		See Price list	☐
73	APHR99	Omit handrim						NC	☐

* Not available with Mag wheels 1 - Available in 22", 24", 25" 2 - Available in 22", 24", 25" 26" 3 - Available in 24", 25" 26"

AXLES									
74	APAX01	Quick release axle						STD	☐
75	APAX03	Quad release axle						\$125	☐

REAR WHEEL CAMBER									
76	APCM00	0° camber						STD	☐
77	APCM02	2° camber						NC	☐
78	APCM04	4° camber						NC	☐
79	APCM06	6° camber						NC	☐
80	APCM08	8° camber						NC	☐

CENTER OF GRAVITY

Choose position of rear wheel axle plate:

Measure distance from front of back post to center of axle receiver.

0 - 6" range

Set At: (inches)

If not indicated, default setting is at 2.5"

WHEEL LOCK									
81	APWL08	Newton wheel lock - push to lock						STD	☐
82	APWL14	6" removable extension handle for Newton push to lock						\$75	☐
83	APWL04	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)						\$150	☐
84	APWL13	6" removable extension handle for pull to lock and grade aid						\$75	☐
85	APWL05	Scissor lock						\$240	☐
86	APWL99	Omit wheel locks						NC	☐

[illegible]


BioFIT HCPCS coding E2607

- Cool feeling Visco Gel top layer
- Multi-layer design with SPS technology
- Anatomically Contoured Cushion design
- Moderate to high risk users



NuFIT HCPCS coding E2607

- Multi-layer design with SPS technology for comfort and positioning
- Anatomically Contoured Cushion design
- Low to medium risk users



FITS WHEELCHAIR WIDTH / DEPTH

10" / 10"
11", 12" / 12"
13", 14" / 14"
13", 14" / 16"
15", 16" / 16"
15", 16" / 18"
15", 16" / 20"
17", 18" / 16"
17", 18" / 18"
17", 18" / 20"
19", 20" / 16"
19", 20" / 18"
19", 20" / 20"
21", 22" / 16"
21", 22" / 18"
21", 22" / 20"
23", 24" / 16"
23", 24" / 18"
23", 24" / 20"

BIOFIT CUSHION MODEL NO.

○ NXCABF1010
○ NXCABF1212
○ NXCABF1414
○ NXCABF1416
○ NXCABF1616
○ NXCABF1618
○ NXCABF1620
○ NXCABF1816
○ NXCABF1818
○ NXCABF1820
○ NXCABF2016
○ NXCABF2018
○ NXCABF2020
○ NXCABF2216
○ NXCABF2218
○ NXCABF2220
○ NXCABF2416
○ NXCABF2418
○ NXCABF2420

RETAIL

395
395
395
395
395
395
395
395
395
395
395
395
395
425
425
425
425
425
425

NUFIT CUSHION MODEL NO.

—
—
—
—
○ NXCANF1616
○ NXCANF1618
○ NXCANF1620
○ NXCANF1816
○ NXCANF1818
○ NXCANF1820
○ NXCANF2016
○ NXCANF2018
○ NXCANF2020
○ NXCANF2216
○ NXCANF2218
○ NXCANF2220

RETAIL

365
365
365
365
365
365
365
365
365
365
365
395
395
395

OPTIMA HA^{BACK} HEIGHT ADJUSTABLE HCPCS coding E2620/E2621

- Stepless back height adjustment from 16" to 20" (KC models vary)
- Supports pelvis and allows independent positioning of the laterals and posterior support
- Pelvic cut out for improved positioning, comfort & support

AC^{BACK}

Mild 3.5" lateral contour



PC^{BACK}

Depth adjustable 5.5 - 7" lateral contour



KC^{BACK}

HCPCS coding E2611/E2612
Gentle 1" lateral contour



FITS WHEELCHAIR WIDTH / BACK HEIGHT

10" / 8" - 10.5"
12" / 9" - 11.5"
14" / 10" - 11.5"
16", 17", 18" ** / 16" - 20"
18", 19", 20" ** / 16" - 20"
20", 21", 22" ** / 16" - 20"
22", 23", 24" ** / 16" - 20"
24", 25", 26" ** / 16" - 20"

HA AC BACK MODEL NO.

—
—
—
○ NHAACB1616/20
○ NHAACB1816/20
○ NHAACB2016/20
○ NHAACB2216/20
○ NHAACB2416/20

RETAIL

625
625
625
675
675

HA PC Back MODEL NO.

—
—
—
○ NHAPCB1616/20
○ NHAPCB1816/20
○ NHAPCB2016/20
○ NHAPCB2216/20
○ NHAPCB2416/20

RETAIL

625
625
625
675
675

HA KC BACK MODEL NO.

○ NHAKCB1010/10
○ NHAKCB1212/11
○ NHAKCB1414/11

RETAIL

625
625
625

OPTIMA FHP^{BACK} FIXED HEIGHT PELVIC CONTOUR HCPCS coding E2620/E2621

- Fixed heights available from 12" - 20"
- Scapula cut out shape to improve range of motion
- Pelvic cut out for improved positioning, comfort & support

AC^{BACK}

Mild 3" lateral contour



PC^{BACK}

Depth adjustable 5.5 - 7" lateral contour



FITS WHEELCHAIR WIDTH / BACK HEIGHT

16", 17", 18" ** / 16"
16", 17", 18" ** / 18"
16", 17", 18" ** / 20"
18", 19", 20" ** / 16"
18", 19", 20" ** / 18"
18", 19", 20" ** / 20"
20", 21", 22" ** / 16"
20", 21", 22" ** / 18"
20", 21", 22" ** / 20"
22", 23", 24" ** / 16"
22", 23", 24" ** / 18"
22", 23", 24" ** / 20"

FHP AC BACK MODEL NO.

○ NACB1616
○ NACB1618
○ NACB1620
○ NACB1816
○ NACB1818
○ NACB1820
○ NACB2016
○ NACB2018
○ NACB2020
○ NACB2216
○ NACB2218
○ NACB2220

RETAIL

625
625
625
625
625
625
625
625
625
675
675
675


FHP PC Back MODEL NO.

○ NPCB1616
○ NPCB1618
○ NPCB1620
○ NPCB1816
○ NPCB1818
○ NPCB1820
○ NPCB2016
○ NPCB2018
○ NPCB2020
○ NPCB2216
○ NPCB2218
○ NPCB2220

RETAIL

625
625
625
625
625
625
625
625
625
675
675
675

Mild 3" lateral contour



MULTI-ADJUSTABLE POSTURAL SUPPORT

- Upper segment adjusts 40° anterior & 20° posterior
- Mounting attachments can be adjusted in width symmetrically or asymmetrically
- Optional width adjustable lateral/pelvic supports
- Lower segment adjusts 30° anterior & 30° posterior



14°, 15°, 16° / 9°
14°, 15°, 16° / 13°
14°, 15°, 16° / 16°
15°, 16°, 17° / 9°
15°, 16°, 17° / 13°
15°, 16°, 17° / 16°
16°, 17°, 18° / 9°
16°, 17°, 18° / 13°
16°, 17°, 18° / 16°
16°, 17°, 18° / 19°
17°, 18°, 19° / 9°
17°, 18°, 19° / 13°
17°, 18°, 19° / 16°
17°, 18°, 19° / 19°
18°, 19°, 20° / 9°
18°, 19°, 20° / 13°
18°, 19°, 20° / 16°
18°, 19°, 20° / 19°
19°, 20°, 21° / 9°
20°, 21°, 22° / 9°
20°, 21°, 22° / 13°
20°, 21°, 22° / 16°
20°, 21°, 22° / 19°
21°, 22°, 23° / 16°
21°, 22°, 23° / 19°
22°, 23°, 24° / 13°
22°, 23°, 24° / 16°
22°, 23°, 24° / 19°

- ☐ NOACB1409
- ☐ NOACB1413
- ☐ NOACB1416
- ☐ NOACB1509
- ☐ NOACB1513
- ☐ NOACB1516
- ☐ NOACB1609
- ☐ NOACB1613
- ☐ NOACB1616
- ☐ NOACB1619
- ☐ NOACB1709
- ☐ NOACB1713
- ☐ NOACB1716
- ☐ NOACB1719
- ☐ NOACB1809
- ☐ NOACB1813
- ☐ NOACB1816
- ☐ NOACB1819
- ☐ NOACB1916
- ☐ NOACB1919
- ☐ NOACB2009
- ☐ NOACB2013
- ☐ NOACB2016
- ☐ NOACB2019
- ☐ NOACB2116
- ☐ NOACB2119
- ☐ NOACB2213
- ☐ NOACB2216
- ☐ NOACB2219

[illegible]

- ☐ NOPCB1509
- ☐ NOPCB1513
- ☐ NOPCB1516
- ☐ NOPCB1609
- ☐ NOPCB1613
- ☐ NOPCB1616
- ☐ NOPCB1619
- ☐ NOPCB1709
- ☐ NOPCB1773
- ☐ NOPCB1716
- ☐ NOPCB1719
- ☐ NOPCB1809
- ☐ NOPCB1813
- ☐ NOPCB1816
- ☐ NOPCB1819
- ☐ NOPCB1916
- ☐ NOPCB1919
- ☐ NOPCB2009
- ☐ NOPCB2013
- ☐ NOPCB2016
- ☐ NOPCB2019
- ☐ NOPCB2116
- ☐ NOPCB2119
- ☐ NOPCB2213
- ☐ NOPCB2216
- ☐ NOPCB2219

[illegible]

14" / 16"
16" / 16"
18" / 18"
20" / 18"
22" / 18"
24" / 18"

- NMAPSB1
- NMAPSB1
- NMAPSB1
- NMAPSB2
- NMAPSB2
- NMAPSB2

795
795
795
795
945
945

HCPCS coding E0955

HEAD SUPPORT

- The MultiFit head support is available in a large or contoured pad design that can be shaped to fit



HEAD SUPPORT STYLE

CONTOURED PAD
LARGE RECTANGULAR PAD
EXTRA LARGE
RECTANGULAR PAD

MULTIFIT MODEL NO.	RETAIL
☐ NIB-0785	225
☐ NIB-0784	225
☐ NIB-0818	225

HCPCS
coding N/A

- Mild 3" lateral contour
- Height adjustment
- 5" of depth
- 25° angle adjustment

- Removable QuickFit
mounting hardware

OR

Check **F** for Fixed QuickFit mounting hardware

- Removable QuickFit
mounting hardware

OR
Check **F** for Fixed QuickFit
mounting hardware

14", 15", 16" / 9"
14", 15", 16" / 13"
14", 15", 16" / 16"
15", 16", 17" / 9"
15", 16", 17" / 13"
15", 16", 17" / 16"
16", 17", 18" / 9"
16", 17", 18" / 13"
16", 17", 18" / 16"
17", 18", 19" / 9"
17", 18", 19" / 13"
17", 18", 19" / 16"
18", 19", 20" / 9"
18", 19", 20" / 13"
18", 19", 20" / 16"
20", 21", 22" / 9"
20", 21", 22" / 13"
20", 21", 22" / 16"
Rel handle 14"-18" backs
Rel handle 18"-22" backs

- ☐ NOCFACB1409
- ☐ NOCFACB1413
- ☐ NOCFACB1416
- ☐ NOCFACB1419
- ☐ NOCFACB1513
- ☐ NOCFACB1516
- ☐ NOCFACB1609
- ☐ NOCFACB1613
- ☐ NOCFACB1616
- ☐ NOCFACB1709
- ☐ NOCFACB1713
- ☐ NOCFACB1716
- ☐ NOCFACB1809
- ☐ NOCFACB1813
- ☐ NOCFACB1816
- ☐ NOCFACB2009
- ☐ NOCFACB2013
- ☐ NOCFACB2016
- ☐ NIB-1174
- ☐ NIB-1052

835
865
895
835
865
895
835
865
895
835
865
895
835
865
895
C Option
C Option

- ☐ NOCFPCB1509
- ☐ NOCFPCB1513
- ☐ NOCFPCB1516
- ☐ NOCFPCB1609
- ☐ NOCFPCB1613
- ☐ NOCFPCB1616
- ☐ NOCFPCB1709
- ☐ NOCFPCB1713
- ☐ NOCFPCB1716
- ☐ NOCFPCB1809
- ☐ NOCFPCB1813
- ☐ NOCFPCB1816
- ☐ NOCFPCB2009
- ☐ NOCFPCB2013
- ☐ NOCFPCB2016

835
865
895
835
865
895
835
865
895
835
865
895
835
865
895

HCPCS coding E1028

MULTI-AXIS
HARDWARE & BACK
MOUNTING BRACKET
HEAD SUPPORT/
FLIP-BACK KIT

- NIB-0747 145
- NIB-1210 95

