



2915 Oglethorpe RD # 2270
Newark, DE 19713
T. +1-866-650-6555 F. +1-888-966-6555
orders@motioncomposites.com
www.motioncomposites.com

4.7



APEX

ORDER FORM - USA



P.O. number: * _____ Date: _____

Company Name: * _____ Phone: * _____

Ordered by: * _____

Clinician/ATP/RTS: _____

Email: _____

Marked for: _____

Ship to: * _____

Name: _____

Address: _____

City: _____ State: _____

Postal Code: _____

Phone: _____

* Required

MODEL				USD		
1	APWC11	Ultralight Aluminum APEX Rigid Wheelchair (265 lb weight limit)	Tripple-buttet aluminum frame / Aluminum back canes and back rigidizer / Aluminum seat rigidizer / Ultralight carbon fiber camber tube and rear frame rigidizer	K0005	\$2 595	☐
2	APWC12	Apex Carbon fiber frame upgrade	Ultralight carbon fiber frame, back canes, back rigidizer, seat rigidizer, camber tube and rear frame rigidizer.		\$1 700	☐


SEAT WIDTH												
3	APSW	Seat width:	12" ☐	13" ☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐		NC	☐
4	APSW	Seat width:	19" ☐	20" ☐						E2201	\$495	☐


SEAT DEPTH												
5	APSD	Seat depth:	12" ☐	13" ☐	14" ☐	15" ☐	16" ☐	17" ☐	18" ☐	19" ☐	NC	☐
6	APSD	Seat depth:	20" ☐						E2203	\$475	☐	

FRAME LENGTH					
7	APSP01	Standard frame length configuration		STD	☐
8	APSP02	Frame length 1" longer than standard configuration		NC	☐
9	APSP03	Frame length 2" longer than standard configuration		NC	☐
10	APSP04	Frame length 1" shorter than standard configuration		NC	☐
11	APSP05	Frame length 2" shorter than standard configuration		NC	☐


FRONT Seat-to-floor height

Caster sizes	Height available +- 1/4"	Specify desired height: _____ (inches)
3" Caster	14" to 21"	
4" Caster	14.5" to 21"	
5" Caster	15.5" to 21"	
6" Caster	16.5" to 21" (N/A with 90° front frame)	





FRONT FRAME ANGLE							
20	APFA70	70° Front frame angle			NC	☐	
21	APFA75	75° Front frame angle			NC	☐	
22	APFA80	80° Front frame angle			NC	☐	
23	APFA85	85° Front frame angle			NC	☐	
24	APFA90	90° Front frame angle			NC	☐	



Footrest width / Taper selection (circle selection)				
Seat Width	Straight	1" Taper	2" Taper	
Outside Measurement	2 ½" Narrower than Seat Width	4 ½" Narrower than Seat Width	6 ½" Narrower than Seat Width	
12"	9.5"	7.5"		
13"	10.5"	8.5"		
14"	11.5"	9.5"		
15"	12.5"	10.5"		
16"	13.5"	11.5"	9.5"	
17"	14.5"	12.5"	10.5"	
18"	15.5"	13.5"	11.5"	
19"	16.5"	14.5"	12.5"	
20"	17.5"	15.5"	13.5"	

Note: higher than 20" Front STFH will reduce the footrest width up to 1" with a 2" taper

FOOTPLATES					
25	APFP12	Tubular open loop ¹		STD	☐
26	APFP13	Tubular with ABS cover	K0040	\$40	☐
27	APFP31	Newton adjustable angle - Aluminum (one piece)	K0040	\$95	☐
28	APFP32	Newton adjustable angle - Carbon fiber <i>(one piece with rubber inserts)</i>	K0040	\$175	☐
29	APFP14	Flip back, adjustable angle aluminum	K0040	\$195	☐
30	APFP15	Flip back, adjustable angle carbon fiber	K0040	\$275	☐
31	APFP16	2" riser <i>(only available with aluminum angle adjustable footplate)</i>	K0040	\$140	☐


1: Not available on all non tapered frames; Only available with 1 inch taper with seat width 17" and greater; Not available with 2 inch taper.

CALF STRAPS						
32	APFP09	Calf strap		K0038	\$35	☐
33	APFP30	Bodypoint Aeromesh Padded calf strap	New	SEFND0010	\$115	☐

FOLD DOWN BACKREST					
34	APBC20	Locking backrest		STD	☐
35	APBC21	Ultralight non locking backrest (not locking in upright and down position)		\$95	☐

Back canes fixed height Without push handle											
ADJUSTABLE ANGLE FOLD DOWN BACKREST											
36	APBC	Fixed height without push handle		(30) 9"□	(31) 9.5"□	(32) 10"□	(33) 10.5"□	(34) 11"□		NC	□
		(35) 11.5"□	(36) 12"□	(37) 12.5"□	(38) 13"□	(39) 13.5"□	(40) 14"□	(41) 14.5"□			
		(42) 15"□	(43) 15.5"□	(44) 16"□	(45) 16.5"□	(46) 17"□	(47) 17.5"□	(48) 18"□			
		(49) 18.5"□	(50) 19"□	(51) 19.5"□	(52) 20"□	(53) 20.5"□	(54) 21"□				

Back canes, height and angle adjustable, with option									
ADJUSTABLE ANGLE BACK									
37	APBC	Straight back canes		(22) 9 - 12" □	(23) 12 - 15" □	(24) 15 - 18" □	(25) 18 - 21" □	NC	☐
Set height to: *		(in increments of 1/2", measured from seat sling to top of back upholstery) * If not indicated, default height setting is at the lowest position							


38		Select back cane angle: ¹	101°□	99°□	96°□	93° ¹ □	90°□	87°□	83°□	80°□		STD	☐
----	--	--------------------------------------	-------	------	------	--------------------	------	------	------	------	---	-----	--------------------------

1) If not indicated default angle setting is at 93° * If not indicated, default height setting is at the lowest position

			PUSH HANDLE OPTION						
39	APPH01	Push handles						STD	☐
40	APPH05	Fold down push handles ¹					New	\$120	☐
41	APPH06	Clamp on height adjustable push handle ²		9 - 20" ☐	STD	9-23" ☐	New	\$300	☐
42	APPH02	Without push handles						NC	☐

1- 265 pounds weight limit

2 - Not compatible with swing away armrests ; 265 pounds weight limit on carbon fiber back cane

BACKRESTS OPTIONS					
43	APBC57	Standard 3 inch deep rigidizer bar		STD	☐
44	APBC55	Four inch deep Aluminum rigidizer bar		NC	☐
45	APBC56	Depth adjustable back (<i>provides 3" of depth adjustment</i>)		\$180	☐
46	APBC88	Custom height rigidizer bar set at _____ (<i>inches</i>) <i>Measured from the top of the seat sling to the center of the rigidizer - standard height is 5.5 "</i> <i>Lowest position is 4.25"</i>		\$125	☐

BACK UPHOLSTERY					
47	APBH03	Slip-on Syntech-Suede back upholstery		STD	☐
48	APBH04	Adjustable tension Syntech-Suede back upholstery	E2611	\$300	☐
49	APBH99	Omit Back upholstery		NC	☐

SEAT CUSHION						
50	APSC01	2" Cushion		E2601	\$165	☐
51	APSC02	3" Cushion		E2601	\$185	☐

REAR WHEELS													
52	APRW01	Mag ¹	20" ☐	22" ☐	24" ☐							\$50	☐
53	APRW02	Newton One - spoke wheel ²	20" ☐	22" ☐	24" ☐							STD	☐
54	APRW03	Newton Gravity - ultralight wheel	20" ☐	22" ☐	24" ☐	25" ☐	26" ☐					\$300	☐
55	APRW04	Spinergy Spox		22" ☐	24" ☐	25" ☐	26" ☐					\$795	☐
56		Spinergy Spox spokes colors* Blk ☐ <i>STD</i>		Wht ☐	Blue ☐	Red ☐	Yel ☐					NC	☐
57	APRW05	Spinergy LX		22" ☐	24" ☐	25" ☐	26" ☐					\$995	☐
58		Spinergy LX spokes colors * Blk ☐ <i>STD</i>	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	Pink ☐	Grn ☐	Org ☐		NC	☐
59	APRW07	Spinergy CLX Carbon Blade			24" ☐	25" ☐					New	\$1 480	☐
60		Spinergy CLX spokes colors * Blk ☐ <i>STD</i>	Pur ☐	Wht ☐	Blue ☐	Red ☐	Yel ☐	Pink ☐	Grn ☐	Org ☐		NC	☐

* Colour spokes require an additional 10 days lead-time

1 - HP tires are not compatible with Mag Wheels -use only aluminum or plastic coated handrim 2- Not compatible with Shox tires

TIRES				Available			
61	APTIO1	Soft Urethane 1 1/4 "		20", 22" 24"	E2218	NC	☐
62	APTIO2	Soft Urethane knobby 1 3/8"		22", 24", 25", 26"	E2218	\$95	☐
63	APTIO4	Soft Urethane - Shox 1"*		22", 24", 25", 26"	E2218	\$115	☐
64	APTIO5	Pneumatic 1 3/8 "		20", 22", 24", 25", 26"	E2211/E2212	STD	☐
65	APTIO8	Pneumatic with Airless insert 1 3/8"		22", 24"	E2211/E2213	\$80	☐
66	APTIO6	Pneumatic HP - 1" **		20", 22", 24", 25", 26"	E2211/E2212	\$95	☐
67	APTIO7	Pneumatic HP - Schwalbe Marathon Plus 1" **		24", 25", 26"	E2211/E2212	\$125	☐

* Not compatible with Newton One rear wheels ** Not compatible with Mag Wheels

HANDRIMS						
68	APHR01	Aluminum anodized				STD ☐
69	APHR07	Newton Air Grip* <i>Black</i> ☐				\$95 ☐
70	APHR	Newton Air Grip * <i>(08) Sunkissed Orange</i> ☐ <i>(09) Burgundy</i> ☐ <i>(10) Sapphire Blue</i> <i>(11) Monster Green</i> ☐ <i>(12) Charcoal</i> ☐ <i>(13) Fuchsia</i> ☐ <i>(14) Steel Blue</i> ☐ <i>(15) Ferrari Red</i> ☐ <i>(16) Acid Green</i> ☐ <i>(17) White</i> ☐ <i>(18) Sky Blue</i> ☐				\$240 ☐
Assembly position for aluminum and Newton AirGrip handrims						
71	APHR29	Standard Handrim position				STD ☐
72	APHR30	Narrow Handrim position - <i>Short tab (not available in 22")</i>				\$25 ☐
73	APHR31	Super Narrow Handrim - <i>Short tab (not available in 22")</i>				\$25 ☐

HANDRIM						
74	APHR02	Plastic coated ²				\$95 ☐
75	APHR03	Q-Grip High friction Neoprene coated* ^{1,2}				\$350 ☐
76	APHR28	Surge * ^{1,2} <i>(1-1/2" Oval handrim w/ Gripton high friction strip)</i>				\$375 ☐
77	APHR05	Surge LT* ^{1,2} <i>(1-1/8" Oval handrim w/ Gripton high friction strip)</i>				\$405 ☐
78	APHR04	Natural Fit* Standard Grip ^{1,2} <i>(1-1/2" Oval handrim w/standard thumb piece)</i>				\$395 ☐
79	APHR25	Natural Fit* Super Grip ^{1,2} <i>(1-1/2" Oval handrim w/super grip thumb piece)</i>				\$450 ☐
80	APHR99	Omit handrim				NC ☐

* Not available with Mag wheels 1 - Not available in 20" 2- Only standard handrim position available

AXLES						
81	APAX01	Quick release axle				STD ☐
82	APAX03	Quad release axle		K0108	\$125	☐

REAR WHEEL CAMBER						
83	APCM00	0° camber				STD ☐
84	APCM02	2° camber				NC ☐
85	APCM04	4° camber				NC ☐
86	APCM06	6° camber				NC ☐
87	APCM08	8° camber				NC ☐

CENTER OF GRAVITY

Specify the position of rear wheel axle plate:


Measure distance from front of back post to center of axle receiver.

0 - 6" range

Set At: _____ (inches)

If not indicated, default setting is at 2.5"

If a configuration issue arises due to the choice of the Center Of Gravity (COG) Motion Composites reserves the right to reduce or increase the COG to the closest position.



WHEEL LOCK				
88	APWL08	Newton wheel lock - push to lock		STD
89	APWL04	Newton grade aid push to lock (suggested with pneumatic or high pressure tire)	E0974	\$150
90	APWL13 APWL14	6" removable extension handle	E0961	\$75
91	APWL05	Scissor lock		\$240
92	APWL99	Omit wheel locks		NC

"T" ARMREST*		Height Adjustment		
93	APAR03	10" Desk Length	8 - 12" ☐ STD 10 - 14" ☐	E0973 \$325 ☐
94	APAR04	14" Full Length	8 - 12" ☐ STD 10 - 14" ☐	E0973 \$325 ☐

* May have a 1/2 " adjustment restriction on both limits of the range

SWING AWAY ARMREST**		Height Adjustment		
95	APAR05	Tubular arm pad	8 - 11"	E0973 \$225 ☐
96	APAR06	10" Desk length arm pad	8 - 11"	E0973 \$225 ☐

**not available with the Ultralight non-locking backrest

SIDE GUARDS (not available with "T" Armrest)				
97	APSG01	Plastic side guards (removable)	K0108	\$195 ☐
98	APSG03	Carbon fiber side guards (removable)	K0108	\$290 ☐
99	APSG06	Plastic fender side guards (removable)* *Not compatible with 20" and 22" rear wheels	WAMV	\$270 ☐
100	APSG04	Carbon fiber fender side guards (removable)	K0108	\$540 ☐

POSITIONING BELT				
101	APSB01	Velcro adjustable	Standard 35" ☐ Long 60" ☐	E0978 \$45 ☐
102	APSB02	Auto buckle	Standard 35" ☐ Long 60" ☐	E0978 \$85 ☐
103	APSB07	Hardware attachment for 2 point seatbelt (1 pair)		NC ☐
104	APSB03	Additional hardware for 2 extra attachment points (needed for the 4 point seat belt) (1 pair)		\$15 ☐
105	APSB04	BodyPoint Hip Belt - 2 point non padded	M ☐ L ☐	E0978 \$80 ☐
106	APSB05	BodyPoint Hip Belt - 2 point padded	S ☐ M ☐ L ☐	E0978 \$170 ☐
107	APSB06	BodyPoint Hip Belt - 4 point padded	S ☐ M ☐ L ☐	E0978 \$190 ☐
108	APSB08	BodyPoint EVOFLEX® Positioning Belt Rehab latch and Band Clamps	M ☐ L ☐	E0978 \$285 ☐
109	APSB99	Omit positioning belt		STD ☐

ACCESSORIES				
110	APAC03	Anti-tippers	E0971	\$140 ☐
111	APAC05	WC-19 Transit Tie Down w/Q'straint belt* *Not available with the swing away armrest ; Not available with HD option	Short belt 20-41" ☐ Long belt 20-46" ☐	\$320 ☐
112	APAC08	Spoke guards** Available: 20", 22", 24", 26"	K0065	\$115 ☐
113	APAC10	Seat pouch		\$75 ☐
114	APAC26	BodyPoint Mobility Bag	New	\$125 ☐
115	APAC23	Newton Backpack		\$120 ☐
116	APAC13	Tool kit		\$45 ☐
117	APAC16	Motion Composites Polo Shirt W/S ☐ W/M ☐ W/L ☐ M/S ☐ M/M ☐ M/L ☐ M/XL ☐ M/XXL ☐ M/XXXL ☐ M/XXXXL ☐		\$32 ☐
118	APAC15	Military Patches Army ☐ Navy ☐ Air Force ☐ Marines ☐ Coast Guard ☐		\$90 ☐
119	APAC15	Combat Wounded Purple Heart Patch		\$90 ☐

** Not available with Mag wheels

COLOR									
120	APFC	Highlight color:	(10) Ferrari Red ☐	(11) White ☐	(12) Black ☐	(13) Sunkissed Orange ☐		NC	☐
			(14) Sky Blue ☐	(15) Fuchsia ☐	(16) Acid Green ☐				
							Total		

Notes :