

ActiveCare Air Series




Low Air Loss APM Mattress

- 3 in 1 powered air support surface with three therapy modes, including Alternating, Pulsation, and Static (Constant low pressure) plus additional therapy with the true low air loss advanced coverlet.
- Advanced true low air loss coverlet blows air directly into the cover which provides microclimate management of skin temperature and moisture.
- Auto Boost helps max inflate the surface during rotation which is intended to provide maximum pressure to the seat section to help with entering and exiting the ActiveCare bed.
- Alternating pressure cycles of 5, 10, 15, 20, 25 & 30 minutes allows clinicians to customize pressure redistribution and skin surface interface.
- Innovative design allows the surface to work with either a right or left hand turn bed as the hose routing can be configured at anytime in the field.
- Mattress dimensions: 36"x80"x8" and 42"x80"x8" with included heel pillow with heel slope.

Mattress Solutions

Deluxe


Multiple-layer pressure redistribution foam provides sufficient comfort.



Prevention


StatMatt

Convertible surface reacts automatically to each patient's weight/body profile.



Reactive

Comfort



Elite

Copper comfort gel infused memory foam provides the ultimate in comfort.

Treatment



ActiveCare Air

Powered air surface for the prevention and treatment of pressure injuries.

Available on:




Deluxe

Premium Aesthetics with Powered Height with OneButtonSafeTurn, Night Light, and Battery BackUp



Standard

Powered Height with OneButtonSafeTurn



Fixed Height

Fixed Height with OneButtonSafeTurn

is at the forefront of innovating medical technology for healthcare. We are focused on partnering with you to help address healthcare challenges in the following key areas—




FALLS PREVENTION



BARIATRIC CARE



MOBILITY



REHABILITATION



WOUND CARE