



Afikim Ultra-Lite Wheelchair



Ultra-Lite Dims

- Dims: 41.9" x 26.6"x 36.8"
- Folded: 14.2"W x 31.3"H
- Seat Width: 16.5"
- Seat Depth: 16.1"
- Seat to Floor: 20.3"
- Backrest: 16.1"
- Weight Limit: 275 Lbs

Weight

- High strength and Light Weight
- 16.2 lbs without Wheels
 - 26.7 lbs overall

Accessories

- Elevating Footrests
- Seatbelts

Warranty

- Lifetime Limited Warranty on wheelchair Frame
- 1 Year on Non-wearable parts

Ultra-Lite

This wheelchair is expertly crafted from a Magnesium-Aluminum alloy, a material commonly used in aviation. Its high strength and lightweight design ensure convenient usability for the user.

Features

- Wheel locks
- Caregiver Braking System
- Swing-Away/Removable Legs
- Solid Flat-Free Tires
- Anti -Tip Wheels
- Raisable Armrests
- Rear Storage Bag & Drink Holder
- Quick release rear wheels

AFIKIM MOBILITY

1011 8th St unit 310 Daytona Beach, FL

Get in Touch

 1-844-313-7116

 orders@afiscooters.com

 www.afiscooters.com