

P O W E R C H A I R



Seat Colors
Charcoal Gray

Body Colors
Candy Apple Red

Roll

New!

^{GOLDEN} *Alanté* **SPORT**

The Alanté name continues to mean quality and value! The first thing that will catch your eye with the all-new Alanté Sport is the unique black and red shroud design that is stylish and functional. Featuring a full size, captain's seat, the Alanté Sport is a front wheel drive power base with our black, non-marking drive tires on high strength steel rims and shiny hubcaps.



YOUR LIFE IN MOTION

Alanté Sport GP-205 Specifications

Model Number	GP205CC
Medicare Code	K0823
Weight Capacity	300 lbs
Drive Wheels	Front
Maximum Speed	3.5 mph
Operating Range ¹	16.5 miles
Ground Clearance	1.75"
Turning Radius	25"
Type Batteries	2-U1
Colors	Red
Freewheel Mode	Yes
Electronic Speed Control	Yes
Electro-Mechanical Brakes	Yes
Charger	Off Board
Controller	PG VR-2 50 Amp Remote
Dimensions:	
Length	40.5"
Width	22.25"
Height (ground to top of back of seat)	37" - 39"
Height (ground to top of back of headrest)	44.5" - 48"
Ground to Top of Seat	20.5" - 22.5"
Number of Seat Height Adjustments	3
Size of Increments	1"
Weight of Unit: (assembled)	
Base	67 lbs
Seat w/Arms	46 lbs
Batteries	22 lbs ea.
Tires:	
Front	9"
Casters	6"
Front Anti-Tip Wheels	4"
Standard Seat:	
Back Height - no headrest	19"
Back Height with headrest	26.5" - 30"
Width x Depth	18" X 18"
Color	Gray Vinyl
Optional Seat:	
Back Height - no headrest	21"
Back Height with headrest	28.5" - 32.5"
Width x Depth	20" x 18"
Color	Coffee/Sand, Charcoal/Gray
Footrest	
Height adjustable	Yes
Angle adjustable	Yes
Warranty:	
Frame	5 Years
Drive Train	13 Months
Electronics ²	13 Months

¹. Battery range will vary due to rider weight, drive surface and terrain.

². Electronics warranty excludes batteries.

Features

- 25" turning radius
- 300 lb weight capacity
- 9" driving wheels with alloy rims
- Enlarged 4" anti-tipper wheels

5 Year Warranty:
Power Chair Frame

13 Month Limited Warranty:
Electronic Controller
Drive Train Components

Power Chair Accessories and Options



Pack-N-Go



Quad Cane Holder




Crutch Holder



Rear Basket



Oxygen Holder



Walker Holder



Cup Holder



Disassembles easily into 6 pieces!

Literature is current at time of printing. Golden Technologies reserves the right to make changes to the product or literature at any time.