

Specifications

“Low” Deck Height Range	3.625" to 19.375"	Weight with Standard Boards	295 lbs.
“High” Deck Height Range	7.375" to 23.125"	Maximum Head Angle	64° ± 3°
Safe Working Load	600 lbs.	Maximum Knee Angle	37° ± 3°
Overall Bed Widths	35.5", 38.5", 41.5"	Power Requirements	120 VAC, 60 Hz
Overall Bed Length	94.00"	Maximum Current Draw	2.7 amps
Overall Bed Weight	270 lbs.	Controller Output Voltage	24 VDC

Assembly AllCare C

NOTE:

Disregard the connections for Step 7 (page 11) in the AllCare User Manual and use the connections shown below for Step 7 of the assembly process for AllCare C

