



THE MOST
COMFORTABLE
RECLINER
IN THE WORLD!

THE 4-MOTOR RECLINER

WITH
COUNTLESS
POSITION
COMBINATIONS



Shown in Badlands
Steel Fabric

ATLAS



The Atlas Collection by **VivaLift®** delivers comfort and relaxation that is fully customizable to meet your needs. A convenient remote with soft-touch buttons adjusts each function independently. Enjoy a beverage with the accessible cup holder. The Atlas collection offers an array of elegant fabrics with contrast stitching to complement any living room or family room.

Pride Mobility® (US) 800-800-8586 | (CAN) 888-570-1113
pridemobility.com | prideliftchairs.com |    

ATLAS

COLLECTION	ATLAS v.2
POSITION	Infinite Lay Flat Relaxed Reclined Upward Lift
WEIGHT CAPACITY	375 lbs.
SEAT-TO-FLOOR	19.5"
SEAT DEPTH	21"
SEAT WIDTH	22"
TOP OF BACK TO SEAT	26"
FULLY PADDED CHAISE	Yes
FOOTREST EXTENSION	Standard
SEATING TYPE	Total Comfort
CHAIR HEIGHT	42"
DISTANCE FROM WALL	16"
OVERALL WIDTH	39"
CHAIR WEIGHT	161 lbs.
WARRANTY	7-Yr. Pro-rated


FUNCTIONS AND FEATURES:

-  POWER HEADREST
-  POWER LUMBAR
-  USB REMOTE
-  LITHIUM BATTERY BACKUP
-  FOOTREST EXTENSION
-  POWER BACKREST

POSITION GUIDE:

-  RELAXED
-  RECLINED
-  UPWARD LIFT
-  INFINITE LAY FLAT

STANDARD FABRIC: BADLANDS



MUSHROOM



WALNUT



STEEL

Note: Measurements may vary +/- 1". Specifications subject to change without prior notice. The information contained herein is correct at the time of publication; we reserve the right to alter specifications without prior notice.



= Indicates shortest available lead time

Recommended HCPCS: E0627-Seat lift mechanism, electric, any type DME providers are responsible for determining the appropriate billing codes when submitting for insurance reimbursement. HCPCS codes provided should not be considered as legal advice and do not guarantee reimbursement. Payer coding, coverage and bundling guidelines may apply.



Pride Mobility® (US) 800-800-8586 | (CAN) 888-570-1113

pridemobility.com | prideliftchairs.com |    