



zeen®

Mobility • Fall Prevention • Independence • Joy

What is the Zeen?

The Zeen is a revolutionary multi-purpose mobility device. It is designed to empower you to stand up from sitting, walk safely without the fear of falling, and sit and rest comfortably when you want.

Zeen's motor-free technology lifts 75% of your body weight, raising you from a seated position to standing. Fully supported through your hips and buttocks, you can safely walk upright with the option of having your hands free. When resting, you can sit at any height, including barstool height, to comfortably interact with family and friends at eye-level. Zeen offers the quality of life and dignity you've been waiting for!

Is the Zeen Right for Me?

- ✓ Are you motivated and determined to regain or maintain your mobility or movement?
- ✓ Do you have moderate grip strength in at least one hand?
- ✓ Do you find it difficult to stand up from a seated position?
- ✓ Have you limited walking/standing due to fear of falling, strength, balance, or endurance challenges?
- ✓ Are you able to learn a new device on your own, or do you have someone who could assist you?

Health Benefits Associated with Standing and Weight Bearing

- Improves bone & joint integrity, circulation, leg strength, and range of motion
- Strengthens cardiovascular and digestive systems
- Decreases risk of heart disease, diabetes, lower back pain, weight gain, and pressure sores
- Improves cognition, alertness, mood, sleep quality, and overall mental wellness
- Increases energy level, balance, and coordination



sit. stand. go

Zeen by Exokinetics, Inc.

19 Hagerty Blvd. Unit 15, West Chester, PA 19382



**Zeen is approved for
FSA/HSA purchase.**

**Watch our
demo video** →





"It has been amazing -- **I am able to walk and stand up straight without the fear of falling.** I know that whenever I needed a few minutes to rest, I have a seat right behind me." -- **Amy G.**

The Zeen Replaces the Need for Multiple Pieces of Costly Equipment

- ✓ **Fall Prevention**
Full support at hips through the seat and seatbelt
 - ✓ **Assisted Sit-to-Stand**
Patented non-motorized lift assist technology
 - ✓ **Stander**
Safely promotes weight bearing
 - ✓ **Walker**
Posterior support for upright and hands-free option
 - ✓ **Transport Chair**
Smooth maneuverability, lightweight, and compact
 - ✓ **Exercise Aid**
Improves stride, balance, endurance, and strength
- In addition, staff-to-patient ratios can be reduced by eliminating the need for two-person gait training therapy sessions

Sit down when you need to,
Stand tall when you want to,
and **Go** wherever life takes you!



Usability & Clearance:

Doorways: Minimum 27"
Turning Diameter: 45"
Width of Seat: 21"
Weight: ~42 Pounds

Folds for Transport

How to Select Your Zeen

The Zeen is designed to be customized for an individual's body and mobility needs. When selecting your Zeen, you must consider three things.

- 1) **Hand Control Style** - Standard (2-Hand Squeeze)
Deluxe (1-Hand Squeeze)
- 2) **Stand Over Height** - Standard 23" - 36" (58-91 cm)
Tall 23" - 38" (58 - 97 cm)
- 3) **Gas Springs for Lift** - Light (50-175 lbs.)
Medium (100-250 lbs.)
Heavy (150-325 lbs.)

Once you know which unit fits your needs, your dealer will determine your height and weight settings so your Zeen will be ready to use right out of the box. No assembly required!

We have various Accessories for the Zeen to compliment your lifestyle: seat cushion, trunk lift straps, cup holder, and more!

DEALER INFO



*Coverage Zone and
Demonstration Types*

Contact Name

P: (833)867-5309

E: email@example.com

Business Address

www.Example.com

