




SIZING GUIDE



Axiom

BACKS

Wheelchair Frame Back Widths




Axiom Back Width	Fits Back Frame Widths	
	Inner Slots	Outer Slots
12" (30.5cm)	11" (28.0cm)	12"-14" (30.5cm-35.5cm)
14" (35.5cm)	12"-13" (30.5cm-33.0cm)	14"-16" (35.5cm-40.5cm)
16" (40.5cm)	14"-15" (35.5cm-38.0cm)	16"-18" (40.5cm-45.5cm)
18" (45.5cm)	16"-17" (40.5cm-43.0cm)	18"-20" (45.5cm-50.5cm)
20" (51.0cm)	18"-19" (46.0cm-48.5cm)	20"-22" (50.5cm-55.5cm)
22" (56.0cm)	20"-21" (51.0cm-53.5cm)	22"-24" (56.0cm-61.0cm)

AXIOM

BACKS

Posterior Back (AXPB) Sizing Guide



A - Height of Foam
B - Apex of Lateral

D - Height of Shell
E - Width of Shell

F - Width of Foam
G - Depth of Shell

H - Depth of Foam

POSTERIOR BACK (AXPB)							
Low							
Assy	A	B	D	E	F	G	H
AXPB1208	10" (25.5cm)	5.5" (14cm)	8" (20.5cm)	9" (23cm)	9.75" (25cm)	2" (5cm)	2.25" (5.5cm)
AXPB1210	12" (30.5cm)	7.5" (19cm)	10" (25.5cm)	9" (23cm)	9.75" (25cm)	2" (5cm)	2.25" (5.5cm)
AXPB1212	14" (35.5cm)	9.5" (24cm)	12" (30.5cm)	9" (23cm)	9.75" (25cm)	2" (5cm)	2.25" (5.5cm)
AXPB1408	10" (25.5cm)	5.5" (14cm)	8" (20.5cm)	11" (28cm)	12" (30.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB1410	12" (30.5cm)	7.5" (19cm)	10" (25.5cm)	11" (28cm)	12" (30.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB1412	14" (35.5cm)	9.5" (24cm)	12" (30.5cm)	11" (28cm)	12" (30.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB1608	10" (25.5cm)	5.5" (14cm)	8" (20.5cm)	13" (33cm)	14" (35.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB1610	12" (30.5cm)	7.5" (19cm)	10" (25.5cm)	13" (33cm)	14" (35.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB1612	14" (35.5cm)	9.5" (24cm)	12" (30.5cm)	13" (33cm)	14" (35.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB1808	10" (25.5cm)	5.5" (14cm)	8" (20.5cm)	15" (38cm)	16" (40.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB1810	12" (30.5cm)	7.5" (19cm)	10" (25.5cm)	15" (38cm)	16" (40.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB1812	14" (35.5cm)	9.5" (24cm)	12" (30.5cm)	15" (38cm)	16" (40.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB2008	10" (25.5cm)	5.5" (14cm)	8" (20.5cm)	17" (43cm)	18" (45.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB2010	12" (30.5cm)	7.5" (19cm)	10" (25.5cm)	17" (43cm)	18" (45.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB2012	14" (35.5cm)	9.5" (24cm)	12" (30.5cm)	17" (43cm)	18" (45.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB2208	10" (25.5cm)	5.5" (14cm)	8" (20.5cm)	19" (48.5cm)	20" (51cm)	2" (5cm)	2.25" (5.5cm)
AXPB2210	12" (30.5cm)	7.5" (19cm)	10" (25.5cm)	19" (48.5cm)	20" (51cm)	2" (5cm)	2.25" (5.5cm)
AXPB2212	14" (35.5cm)	9.5" (24cm)	12" (30.5cm)	19" (48.5cm)	20" (51cm)	2" (5cm)	2.25" (5.5cm)

AXION


BACKS

(CONTINUED) POSTERIOR BACK (AXPB)							
Mid							
Assy	A	B	D	E	F	G	H
AXPB1214	16" (40.5cm)	7.5" (19cm)	14" (35.5cm)	9" (23cm)	9.75" (25cm)	2" (5cm)	2.25" (5.5cm)
AXPB1216	18" (45.5cm)	8.5" (21.5cm)	16" (40.5cm)	9" (23cm)	9.75" (25cm)	2" (5cm)	2.25" (5.5cm)
AXPB1218	20" (51cm)	9.5" (24cm)	18" (45.5cm)	9" (23cm)	9.75" (25cm)	2" (5cm)	2.25" (5.5cm)
AXPB1414	16" (40.5cm)	7.5" (19cm)	14" (35.5cm)	11" (28cm)	12" (30.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB1416	18" (45.5cm)	8.5" (21.5cm)	16" (40.5cm)	11" (28cm)	12" (30.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB1418	20" (51cm)	9.5" (24cm)	18" (45.5cm)	11" (28cm)	12" (30.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB1614	16" (40.5cm)	7.5" (19cm)	14" (35.5cm)	13" (33cm)	14" (35.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB1616	18" (45.5cm)	8.5" (21.5cm)	16" (40.5cm)	13" (33cm)	14" (35.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB1618	20" (51cm)	9.5" (24cm)	18" (45.5cm)	13" (33cm)	14" (35.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB1814	16" (40.5cm)	7.5" (19cm)	14" (35.5cm)	15" (38cm)	16" (40.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB1816	18" (45.5cm)	8.5" (21.5cm)	16" (40.5cm)	15" (38cm)	16" (40.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB1818	20" (51cm)	9.5" (24cm)	18" (45.5cm)	15" (38cm)	16" (40.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB2014	16" (40.5cm)	7.5" (19cm)	14" (35.5cm)	17" (43cm)	18" (45.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB2016	18" (45.5cm)	8.5" (21.5cm)	16" (40.5cm)	17" (43cm)	18" (45.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB2018	20" (51cm)	9.5" (24cm)	18" (45.5cm)	17" (43cm)	18" (45.5cm)	2" (5cm)	2.25" (5.5cm)
AXPB2214	16" (40.5cm)	7.5" (19cm)	14" (35.5cm)	19" (48.5cm)	20" (51cm)	2" (5cm)	2.25" (5.5cm)
AXPB2216	18" (45.5cm)	8.5" (21.5cm)	16" (40.5cm)	19" (48.5cm)	20" (51cm)	2" (5cm)	2.25" (5.5cm)
AXPB2218	20" (51cm)	9.5" (24cm)	18" (45.5cm)	19" (48.5cm)	20" (51cm)	2" (5cm)	2.25" (5.5cm)

AXION

BACKS

Posterior Lateral and Lateral (AXPL & AXLB) Sizing Guide



A - Height of Foam
B - Apex of Lateral

D - Height of Shell
E - Width of Shell

F - Width of Foam
G - Depth of Shell

H - Depth of Foam

POSTERIOR LATERAL BACK (AXPL) & LATERAL BACK (AXLB)

Low

Assy	A	B	D	E	F	G	H
AXPL/LB1208	10" (25.5cm)	5.5" (14cm)	8" (20.5cm)	10" (25.5cm)	9" (23cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1210	12" (30.5cm)	7.5" (19cm)	10" (25.5cm)	10" (25.5cm)	9" (23cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1212	14" (35.5cm)	9.5" (24cm)	12" (30.5cm)	10" (25.5cm)	9" (23cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1408	10" (25.5cm)	5.5" (14cm)	8" (20.5cm)	12" (30.5cm)	11" (28cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1410	12" (30.5cm)	7.5" (19cm)	10" (25.5cm)	12" (30.5cm)	11" (28cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1412	14" (35.5cm)	9.5" (24cm)	12" (30.5cm)	12" (30.5cm)	11" (28cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1608	10" (25.5cm)	5.5" (14cm)	8" (20.5cm)	14" (35.5cm)	13" (33cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1610	12" (30.5cm)	7.5" (19cm)	10" (25.5cm)	14" (35.5cm)	13" (33cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1612	14" (35.5cm)	9.5" (24cm)	12" (30.5cm)	14" (35.5cm)	13" (33cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1808	10" (25.5cm)	5.5" (14cm)	8" (20.5cm)	16" (40.5cm)	15" (38cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1810	12" (30.5cm)	7.5" (19cm)	10" (25.5cm)	16" (40.5cm)	15" (38cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1812	14" (35.5cm)	9.5" (24cm)	12" (30.5cm)	16" (40.5cm)	15" (38cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB2008	10" (25.5cm)	5.5" (14cm)	8" (20.5cm)	18" (45.5cm)	17" (43cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB2010	12" (30.5cm)	7.5" (19cm)	10" (25.5cm)	18" (45.5cm)	17" (43cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB2012	14" (35.5cm)	9.5" (24cm)	12" (30.5cm)	18" (45.5cm)	17" (43cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB2208	10" (25.5cm)	5.5" (14cm)	8" (20.5cm)	20" (51cm)	19" (48.5cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB2210	12" (30.5cm)	7.5" (19cm)	10" (25.5cm)	20" (51cm)	19" (48.5cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB2212	14" (35.5cm)	9.5" (24cm)	12" (30.5cm)	20" (51cm)	19" (48.5cm)	3.5" (9cm)	3.75" (9.5cm)

AXION

BACKS


(CONTINUED) POSTERIOR LATERAL BACK (AXPL) & LATERAL BACK (AXLB)

Mid							
Assy	A	B	D	E	F	G	H
AXPL/LB1214	16" (40.5cm)	7.5" (19cm)	14" (35.5cm)	10" (25.5cm)	9" (23cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1216	18" (45.5cm)	8.5" (21.5cm)	16" (40.5cm)	10" (25.5cm)	9" (23cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1218	20" (51cm)	9.5" (24cm)	18" (45.5cm)	10" (25.5cm)	9" (23cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1414	16" (40.5cm)	7.5" (19cm)	14" (35.5cm)	12" (30.5cm)	11" (28cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1416	18" (45.5cm)	8.5" (21.5cm)	16" (40.5cm)	12" (30.5cm)	11" (28cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1418	20" (51cm)	9.5" (24cm)	18" (45.5cm)	12" (30.5cm)	11" (28cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1614	16" (40.5cm)	7.5" (19cm)	14" (35.5cm)	14" (35.5cm)	13" (33cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1616	18" (45.5cm)	8.5" (21.5cm)	16" (40.5cm)	14" (35.5cm)	13" (33cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1618	20" (51cm)	9.5" (24cm)	18" (45.5cm)	14" (35.5cm)	13" (33cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1814	16" (40.5cm)	7.5" (19cm)	14" (35.5cm)	16" (40.5cm)	15" (38cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1816	18" (45.5cm)	8.5" (21.5cm)	16" (40.5cm)	16" (40.5cm)	15" (38cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1818	20" (51cm)	9.5" (24cm)	18" (45.5cm)	16" (40.5cm)	15" (38cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB2014	16" (40.5cm)	7.5" (19cm)	14" (35.5cm)	18" (45.5cm)	17" (43cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB2016	18" (45.5cm)	8.5" (21.5cm)	16" (40.5cm)	18" (45.5cm)	17" (43cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB2018	20" (51cm)	9.5" (24cm)	18" (45.5cm)	18" (45.5cm)	17" (43cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB2214	16" (40.5cm)	7.5" (19cm)	14" (35.5cm)	20" (51cm)	19" (48.5cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB2216	18" (45.5cm)	8.5" (21.5cm)	16" (40.5cm)	20" (51cm)	19" (48.5cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB2218	20" (51cm)	9.5" (24cm)	18" (45.5cm)	20" (51cm)	19" (48.5cm)	3.5" (9cm)	3.75" (9.5cm)
Tall							
Assy	A	B	D	E	F	G	H
AXPL/LB1220	22" (56cm)	10.5" (26.5cm)	20" (51cm)	10" (25.5cm)	9" (23cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1222	24" (61cm)	12.5" (32cm)	22" (56cm)	10" (25.5cm)	9" (23cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1420	22" (56cm)	10.5" (26.5cm)	20" (51cm)	12" (30.5cm)	11" (28cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1422	24" (61cm)	11.5" (29cm)	22" (56cm)	12" (30.5cm)	11" (28cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1620	22" (56cm)	10.5" (26.5cm)	20" (51cm)	14" (35.5cm)	13" (33cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1622	24" (61cm)	11.5" (29cm)	22" (56cm)	14" (35.5cm)	13" (33cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1820	22" (56cm)	10.5" (26.5cm)	20" (51cm)	16" (40.5cm)	15" (38cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB1822	24" (61cm)	11.5" (29cm)	22" (56cm)	16" (40.5cm)	15" (38cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB2020	22" (56cm)	10.5" (26.5cm)	20" (51cm)	18" (45.5cm)	17" (43cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB2022	24" (61cm)	11.5" (29cm)	22" (56cm)	18" (45.5cm)	17" (43cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB2220	22" (56cm)	10.5" (26.5cm)	20" (51cm)	20" (51cm)	19" (48.5cm)	3.5" (9cm)	3.75" (9.5cm)
AXPL/LB2222	24" (61cm)	11.5" (29cm)	22" (56cm)	20" (51cm)	19" (48.5cm)	3.5" (9cm)	3.75" (9.5cm)

AXIOM

BACKS

Deep Lateral (AXDL) Sizing Guide



A - Height of Foam
B - Apex of Lateral

D - Height of Shell
E - Width of Shell

F - Width of Foam
G - Depth of Shell

H - Depth of Foam

DEEP LATERAL (AXDL)							
Mid							
Assy	A	B	D	E	F	G	H
AXDL1214	16" (40.5cm)	6" (15cm)	14" (35.5cm)	10.5" (26.5cm)	8.75" (22cm)	5.75" (14.5cm)	6" (15cm)
AXDL1216	18" (45.5cm)	7.25" (18.5cm)	16" (40.5cm)	10.5" (26.5cm)	8.75" (22cm)	5.75" (14.5cm)	6" (15cm)
AXDL1218	20" (51cm)	8.5" (21.5cm)	18" (45.5cm)	10.5" (26.5cm)	8.75" (22cm)	5.75" (14.5cm)	6" (15cm)
AXDL1414	16" (40.5cm)	6.25" (16cm)	14" (35.5cm)	12.5" (32cm)	10.75" (27.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL1416	18" (45.5cm)	7.5" (19cm)	16" (40.5cm)	12.5" (32cm)	10.75" (27.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL1418	20" (51cm)	8.75" (22cm)	18" (45.5cm)	12.5" (32cm)	10.75" (27.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL1614	16" (40.5cm)	6.25" (16cm)	14" (35.5cm)	14.5" (37cm)	12.75" (32.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL1616	18" (45.5cm)	7.5" (19cm)	16" (40.5cm)	14.5" (37cm)	12.75" (32.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL1618	20" (51cm)	8.75" (22cm)	18" (45.5cm)	14.5" (37cm)	12.75" (32.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL1814	16" (40.5cm)	6.5" (16.5cm)	14" (35.5cm)	16.5" (42cm)	14.75" (37.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL1816	18" (45.5cm)	7.75" (19.5cm)	16" (40.5cm)	16.5" (42cm)	14.75" (37.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL1818	20" (51cm)	9" (23cm)	18" (45.5cm)	16.5" (42cm)	14.75" (37.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL2014	16" (40.5cm)	6.5" (16.5cm)	14" (35.5cm)	18.5" (47cm)	16.75" (42.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL2016	18" (45.5cm)	7.75" (19.5cm)	16" (40.5cm)	18.5" (47cm)	16.75" (42.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL2018	20" (51cm)	9" (23cm)	18" (45.5cm)	18.5" (47cm)	16.75" (42.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL2214	16" (40.5cm)	6.5" (16.5cm)	14" (35.5cm)	20.5" (52cm)	18.75" (47.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL2216	18" (45.5cm)	7.75" (19.5cm)	16" (40.5cm)	20.5" (52cm)	18.75" (47.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL2218	20" (51cm)	9" (23cm)	18" (45.5cm)	20.5" (52cm)	18.75" (47.5cm)	5.75" (14.5cm)	6" (15cm)

AXION

BACKS

(CONTINUED) DEEP LATERAL (AXDL)

Tall

Assy	A	B	D	E	F	G	H
AXDL1220	22" (56cm)	8.5" (21.5cm)	20" (51cm)	10.5" (26.5cm)	8.75" (22cm)	5.75" (14.5cm)	6" (15cm)
AXDL1222	24" (61cm)	9.75" (25cm)	22" (56cm)	10.5" (26.5cm)	8.75" (22cm)	5.75" (14.5cm)	6" (15cm)
AXDL1420	22" (56cm)	8.75" (22cm)	20" (51cm)	12.5" (32cm)	10.75" (27.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL1422	24" (61cm)	10" (25.5cm)	22" (56cm)	12.5" (32cm)	10.75" (27.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL1620	22" (56cm)	9" (23cm)	20" (51cm)	14.5" (37cm)	12.75" (32.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL1622	24" (61cm)	10.25" (26cm)	22" (56cm)	14.5" (37cm)	12.75" (32.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL1820	22" (56cm)	9.25" (23.5cm)	20" (51cm)	16.5" (42cm)	14.75" (37.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL1822	24" (61cm)	10.5" (26.5cm)	22" (56cm)	16.5" (42cm)	14.75" (37.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL2020	22" (56cm)	9.25" (23.5cm)	20" (51cm)	18.5" (47cm)	16.75" (42.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL2022	24" (61cm)	10.5" (26.5cm)	22" (56cm)	18.5" (47cm)	16.75" (42.5cm)	5.75" (14.5cm)	6" (15cm)
AXDL2222	24" (61cm)	10.75" (27.5cm)	22" (56cm)	20.5" (52cm)	18.75" (47.5cm)	5.75" (14.5cm)	6" (15cm)