



AXIOM ENDEAVOUR PUSH CHAIR



Color : Size:

<input type="checkbox"/> Red	<input type="checkbox"/> 2
<input type="checkbox"/> Navy	<input type="checkbox"/> 3
	<input type="checkbox"/> 4

Accessories :

<input type="checkbox"/> Head Rest	<input type="checkbox"/> Anti-Tippers
<input type="checkbox"/> Seat Abductor	<input type="checkbox"/> Bug Canopy
<input type="checkbox"/> Hip Abductor - Single	<input type="checkbox"/> Seat Bottom Insert Pad
<input type="checkbox"/> Hip Abductor - Double	<input type="checkbox"/> Seat Back Insert Pad
<input type="checkbox"/> 2pc Lateral Support	<input type="checkbox"/> Swivel Wheel Axis Kit
<input type="checkbox"/> Positioning Wedges (Set)	<input type="checkbox"/> Travel Bag
<input type="checkbox"/> Harness Extension	<input type="checkbox"/> Rain Canopy Size 2,3
<input type="checkbox"/> Sheerling Insert	<input type="checkbox"/> Rain Canopy Size 4
<input type="checkbox"/> Flashing Light	

Customizations:

**Always require a
measurement chart.
ALL CUSTOMS ARE
FINAL SALE**

<input type="checkbox"/> Seatback Height Increase to _____" (max height 35")	<input type="checkbox"/> Beach Wheels
<input type="checkbox"/> Seat Depth Increase to _____" (max height 18")	<input type="checkbox"/> Adjustable Lower Leg Lift
	<input type="checkbox"/> 8" Swivel Wheel Axis Kit

Axiom Endeavour Order Form

AXIOM ENDEAVOUR PUSH CHAIR FOLDED W/OUT WHEELS



Specifications :

Size	2	3	4
Weight Limit	150 LBS	200 LBS	250 LBS
Seatback Height	28"	31"	33"
Seat Depth	15"	15"	18"
Seat Width	15"	18"	21"
Lower Leg Length	15"-22"	15"-24"	17"-24"
Assembled Size	61x24x43	61x27x43	69X 29 1/4 X53
Folded w/ Wheels	52x24x19	52x27x19	57X291/4X20
Folded w/o Wheels	45x19x17.5	45x22x18	51x24x19
Unit Weight	34.2 LBS	36 LBS	38.8 LBS

Included with Chair

*Seat	*Zipper Pouch
*5 Point Harness	*16" Quick release Wheels
*Sun Canopy	*Hand Brake
*Basket underneath chair	*Rear Footbrake
*Adjustable Foot Buckets	*2 Water bottle Holder
*Fender	*Adjustable Push Handle
*Saftey Wrist Strap	

***Multi positional patient system with integrated seat, operated by the caregiver.**

PO Box 369
Selah WA 98942
Toll Free 800-278-9626
Fax 509-453-8326
Adaptivestar.com

Dealer Name: _____

Company Name: _____

Phone: _____

E-mail: _____

Ship to Address _____

Reference Name _____