



AXIOM IMPROV  
PUSH CHAIR



AXIOM IMPROV PUSH CHAIR  
FOLDED W/OUT WHEELS



PO Box 369  
Selah WA 98942  
Toll Free 800-278-9626  
Fax 509-453-8326  
Adaptivestar.com

Dealer Name: \_\_\_\_\_  
Company Name: \_\_\_\_\_  
Phone: \_\_\_\_\_  
E-mail: \_\_\_\_\_  
Ship to Address \_\_\_\_\_  
Reference Name \_\_\_\_\_

### Specifications :

Size	2	3
Weight Limit	150 LBS	200 LBS
Seatback Height	28"	31"
Seat Depth	15"	15"
Seat Width	15"	18"
Lower Leg Length	17"	17"
Assembled Size	63.5x24x46	63.5x27x45
Folded w/ Wheels	51.5x24x19	51.5x27x19
Folded w/o Wheels	43.5x19x12.25	43.5x22x12.25
Unit Weight	32.8 LBS	35.2 LBS

Color : Size:

☐ Red ☐ 2  
☐ Navy ☐ 3

### Accessories :

<input type="checkbox"/> Head Rest	<input type="checkbox"/> Anti-Tippers
<input type="checkbox"/> Seat Abductor	<input type="checkbox"/> Bug Canopy
<input type="checkbox"/> Hip Abductor - Single	<input type="checkbox"/> Seat Bottom Insert Pad
<input type="checkbox"/> Hip Abductor - Double	<input type="checkbox"/> Seat Back Insert Pad
<input type="checkbox"/> 2pc Lateral Support	<input type="checkbox"/> Swivel Wheel Axis Kit
<input type="checkbox"/> Positioning Wedges (Set)	<input type="checkbox"/> Travel Bag
<input type="checkbox"/> Harness Extension	<input type="checkbox"/> Rain Canopy
<input type="checkbox"/> Sheerling Insert	<input type="checkbox"/> Fender
<input type="checkbox"/> Flashing Light	

### Customizations:

☐ Seatback Height Increase to \_\_\_\_\_" (max height 35")  
☐ Seat Depth Increase to \_\_\_\_\_" (max height 18")  
☐ Beach Wheels  
☐ 8" Swivel Wheel Axis Kit

Always require a  
measurement chart.  
ALL CUSTOMS ARE  
FINAL SALE

### Included with Chair

*Seat	*16" Single Fixed Front Wheel
*5 Point Harness	*Two 16" Quick Release Rear Wheels
*Sun Canopy	*Hand Brake
*Basket underneath chair	*Rear Footbrake
*Zipper Pouch	*2 Water bottle Pocket Holder
*Adjustable Push Handle	*Safety Wrist Strap

\*Multi positional patient system with integrated seat, operated by the caregiver.