



Axiom LASSEN Order Form

PO Box 369
Selah WA 98942
Toll Free 800-278-9626
Fax 509-453-8326
Adaptivestar.com

AXIOM LASSEN
PUSH CHAIR



AXIOM LASSEN PUSH CHAIR
FOLDED W/OUT WHEELS



Dealer Name: _____
Company Name: _____
Phone: _____
E-mail: _____
Ship to Address: _____
Reference Name: _____

Specifications :

Size	2	3	4
Weight Limit	150 LBS	200 LBS	250 LBS
Seatback Height	28"	31"	33"
Seat Depth	15"	15"	18"
Seat Width	15"	18"	21"
Lower Leg Length	17"	17"	18"
Assembled Size	47x24x44	47x27x44	47x30x44
Folded w/ Wheels	49x24x16	49x27x16	49x29x16
Folded w/o Wheels	43x19x11	43x22x11	43x25x11
Unit Weight	31.7 LBS	33.2 LBS	36.2 LBS

Color :

☐

Red

Size:

☐

2

☐

Navy

☐

4

Accessories:

☐

Head Rest

☐

Seat Abductor

☐

Hip Abductor - Single

☐

Hip Abductor - Double

☐

2pc Lateral Support

☐

Positioning Wedges (Set)

☐

Harness Extension

☐

Sheerling Insert

☐

Flashing Light

☐

Anti-Tippers

☐

Seat Bottom Insert Pad

☐

Seat Back Insert Pad

☐

Bug Canopy

☐

Swivel Wheel Axis Kit

☐

Travel Bag

☐

Rain Canopy Size 2,2

☐

Rain Canopy Size 4

Customizations:

Always require a
measurement chart.
ALL CUSTOMS ARE FINAL
SALE

☐

Seatback Height Increase to
_____ " (max height 35")

☐

Seat Depth Increase to
_____ " (max height 18")

☐

(2) 8" swivel wheels

Included with Chair

*Seat	*Safety Wrist Strap
*5 Point Harness	*Rear Footbrake
*Sun Canopy	*2 Water bottle Holder
*Basket underneath chair	*16" Rear Quick Release Wheels
*Zipper Pouch	* (2) 6" Front Swivel Wheels
*Adjustable Push Handle	