

# Kid Kart™



## **About Kid Kart**

Kid Kart offers eye-catching contemporary styling and clinically sophisticated yet user-friendly components to give you the confidence that Kid Kart is the right product choice for you.

The positioning shell and seating system on the Xpress and TLC offer a 65° tilt range, extensive adjustment capabilities to accommodate growth, and options for recline and transit.

The Mighty Lite is an ultra lightweight portable chair that's perfect for people on the go.

Kid Kart Xpress™


Kid Kart TLC™

Kid Kart Mighty Lite™





#### **Kid Kart Xpress**



**Kid Kart TLC** 

Canopy shown in Cross Over Stripe style

Lightweight high strength frame is easy to fold and lift for transport

Height adjustable push handle.

> Convenient storage compartment carries up to 15 lbs.

Multiple quick release rear wheel and front caster options make travel easy over any terrain

> Canopy shown in Center Stripe style

rearward facing

Ability to position shell forward or

Removable storage pack

Comfortable foam padded

varying caregiver needs.

height-adjustable handlebar for



Portable: Folds with or without the seating system.



Quick-release rear wheels

Kid Kart TLC in reverse seat configuration

#### **Color Options**

All canopies available in Cross Over or Center Stripe style. See opposite page for pictures. Cross Over option shown below.









Purple

Orange

Soft Pink

Navy Blue



Green



Butterflies (with Purple Highlights)



Red Race Flags (with Red Highlights)

## **Xpress and TLC Base Accessories**



Height Adjustable Booster Base





Rain Cape

I.V. Pole/O2 Tank Holder Option and Vent Tray




### **Mighty Lite Accessories**



3-Point Positioning Belt: Increases pelvic stabalization.



Canopy: Height adjustable canopy provides cover over the headrest and includes a flapcovered observation window.



Foam Wedge: Modify positioning and increase stability with this wedge insert.



H-Harness & Padded Covers: Imitation lamb's wool increases comfort while the H-shape of the harness provides additional support.




Utliity Bag: Keep essentials items close at hand with this roomy, zippered bag.



Headrest Extension: Fully adjustable in both width and height.

- Black nylon fabric is standard
  Navy or Raspberry Cordora is available as a no charge option.



## **Specifications**

	Xpress	TLC
HCPCS Code	E1232	E1232
	E1234 (without seating)	E1234 (without seating)
Tilt Range	20° anterior	0° to 30°
	45° posterior	
Recline Range	0° to 45°	85° to 130°
Footbed Adjustment	70° to 170°	90° to 160°
User Weight Capacity	55 lbs.	75 lbs.
Weight	34 lbs	43 lbs
Shell only	16 lbs	22 lbs
Base only	18 lbs	21 lbs
Width	24"	26.5"
Handle bar height	36" to 41"	37" to 41"
Folded Dimensions	24" × 28" × 20"	
Shell only	16" x 19" x 14"	12.5" x 18" x 28"
Base only	24" × 28" × 15"	26.5" x 14" x 28"
Booster base height range	10.5" - 24"	N/A on TLC
Feeder base height	24"	24"

Might	y Lite		
E1236			
Black			
Navy or Raspberry			
15° and 22°			
85°, 90°, 95°			
From 4" up to 26" in <sup>3</sup> / <sub>4</sub> " increments			
1 <i>5</i> 0 lbs	s		
12	14	16	18
11.75"	13.75"	15.5"	17.5"
11.5"	13"	13.5"	15"
13.5"	15"	15.5"	17"
22.25"	21.25"	23.75"	23.5"
7"	7"	7"	7"
	E1236 Black Navy o 15° an 85°, 90 From 4 150 lbs 12 11.75" 11.5" 13.5" 22.25"	Black Navy or Raspber 15° and 22° 85°, 90°, 95° From 4" up to 26 150 lbs 12 14 11.75" 13.75" 11.5" 13" 13.5" 15" 22.25" 21.25"	E1236 Black Navy or Raspberry 15° and 22° 85°, 90°, 95° From 4" up to 26" in <sup>3</sup> ⁄4" ind 150 lbs 12 14 16 11.75" 13.75" 15.5" 11.5" 13" 13.5" 13.5" 15" 15.5" 22.25" 21.25" 23.75"

SUNRISE MEDICAL.

Customer Service: 800.333.4000 www.SunriseMedical.com