



Reclining back with Heavy Duty seat and base with wide seat and 375 lb. weight capacity, provides maximum comfort stability and support.

FEATURES:

- Heavy duty seat and base provides maximum stability and support.
- Fits most style bathtubs.
- Makes getting in and out of the bathtub safe, simple and worry-free.
- Fully waterproof, floating battery hand control with easy to use convex and concave operating buttons.
- Easy release suction cups on base provide maximum stability in the tub.
- Comfortable, machine washable fabric covers.
- Safe. Will not lower unless there is sufficient power to raise the user back out of the tub.
- Simple to understand charging process.
- 375 lb. weight capacity.
- Limited 3 Year Warranty.





FULLY WATERPROOF, FLOATING BATTERY HAND CONTROL WITH EASY TO USE CONVEX AND CONCAVE OPERATING BUTTONS!



SPECIFICATIONS

Item #	Seat Height Range	Seat Width w/Flaps	Recline	Powerlift and Recline	Product Weight	Carton Shipping Weight	Weight Capacity
BL100-DR	3.1 " - 18.1 "	27 "	up to 40°	Yes	28.9 lbs.	39.5 lbs.	375 lbs.

^{*}Important Notice: Item numbers that display an asterisk (*) reflect an item that will be discontinued or package change.