

NEW!



AquaJoy Saver Fixed Back Bathlift

1P BL200-DR 1/ea

Non-reclining back allows for spacious leg room. Heavy Duty seat and base with wide seat and 350 lb. weight capacity, provides maximum stability and support.

FEATURES:

- Non-reclining, sturdy one piece design maximizes legroom.
- Fits most style bathtubs.
- Makes getting in and out of the bathtub safe, simple and worry-free.
- Fully waterproof, floating battery hand control with easy to use convex and concave operating buttons.
- Easy release suction cups on base provides maximum stability in the tub.
- Safe. Will not lower unless there is sufficient power to raise the user back out of the tub.
- 350 lb. weight capacity.
- Limited 3 Year Warranty.



SPECIFICATIONS

1P = 1-Pack Available

Item #	Seat Height	Seat Width w/Flap	Product Weight	Carton Shipping Weight	Weight Capacity
BL200-DR	2.7" - 18.3"	28.5"	24.2 lbs.	32 lbs.	350 lbs.

Important Notice: Item numbers that display an asterisk () reflect an item that will be discontinued or package change.