


## PRO POOL LIFT DIMENSIONS & SPECIFICATIONS



### DIMENSIONS:

Height at deck level:	38" (97cm)
Width at deck level:	27-1/2" (69cm)
Total footprint on deck:	27-1/2" (69cm) wide x 46" (117cm) long
Set back range:	Maximum 14" (35cm) to center of front anchors
Lowering depth:	27" (69cm) to top of seat surface
Product weight:	120 Pounds (54kg)
Shipping weight:	200 Pounds (91kg)
Shipping Dimensions:	48" x 42" x 48" (122 x 107 x 122cm)
Shipping Method:	LTL Freight

### SPECIFICATIONS:

Frame:	Stainless Steel (type 304)
Finish:	White UV Stabilized epoxy power coat
Actuator:	24v Electric Linear Actuator
Power:	Rechargeable 24v battery (120v Charger included with lift)
Controls:	Waterproof push button
Chair & Footrest:	UV Stabilized Polyethylene (LLDPE)
Weight capacity:	450 pounds (204kg) – (400 pounds (181kg) for deck to water levels greater than 6")
Warranty:	Limited lifetime structural warranty