

# BRADSHAWJUNIOR™

## Fully Electric, Height-Adjustable Pediatric Bed



## Safe, Simple Sleep Positioning for Toddlers to Teens

The Bradshaw Junior is the smart choice for providers serving families who need clinical positioning without complex customization. Reimbursable under HCPCS E0329 and ideal for conditions like cerebral palsy, muscular dystrophy, and other non-behavioral needs, it delivers support and safety without bulk and high price tags, making it a go-to option for providers and families alike.

It's easy to deliver, quick to assemble and even easier to maintain, helping your team save time and reduce costs with every install. And with intuitive safety features, an easy-to-clean slat-deck frame and a warm oak finish that blends in at home, families love it too.

## Benefits & Features

### Clinical Functionality

Provides adjustable foot elevation and Reverse Trendelenburg positioning to support respiratory function, circulation, and pressure redistribution - plus an adjustable height range promotes transfer safety for caregivers

### Saves on Delivery & Setup

Ships on a single pallet, can be easily assembled, and mounts onto its own transport stand, making storage and moves quick and efficient; each section weighs under 44 lb.

### Safety & Security

Features a lockbox to prevent tampering with controls & side rails with two-stage locking mechanism

### Fewer Service Calls for Happier Customers

Simple, durable construction means fewer maintenance issues, reduced follow-ups & more satisfied end users



Mattress Included



## Operational Functions



Raise  
Backrest



Knee  
Break



Auto  
Contour



Variable  
Tilt



Variable  
Height

## Specifications

Item #	BSBOAK-STD-FULL
Overall Length	84.6"
Overall Width	36.6"
Platform Width	30.7"
Platform Length	74.8"
Foot Down Tilt	0°-12°
Max User Weight	363 lb.
Platform Height	15.4 – 31.9"
Caster Size	6"



877.224.0946  [www.drivemedical.com](http://www.drivemedical.com)

 99 Seaview Blvd., Port Washington, NY 11050

©2025 Medical Depot, Inc. dba Drive Medical, 99 Seaview Boulevard, Port Washington, NY 11050. All trademarks used in association with the sale of products of Drive Medical are trademarks owned by Medical Depot, Inc. All other trademarks, trade names, service marks and logos referenced herein belong to their respective companies. Refer to product instruction guide for warranty details.