

Bradshaw Junior™ Bed

www.inspiredbydrive.com



Dealer: _____ Marked For: _____ Account #: _____

Date: _____ PO #: _____ Ordered By: _____

Phone: _____ Fax: _____ Email: _____

Ship to Address: _____

BRADSHAWJUNIOR™



Package #	Description	MSRP	HCPCS
BSBOAK-STD-FULL	Bradshaw Bed, Standard, Oak, w/ Lock Box	\$5,500	E0329

Specifications

Weight Capacity	363 lb.
Color	Honey Oak
Overall Width	36.6"
Overall Length	84.6"
Platform Width	30.7"
Platform Length	74.8"
Overall Height Range (min-max)	36.6" – 53.1"
Platform Height Range (min-max)	15.4" - 31.9"
Under bed clearance (to underside of frame)	13.9"
Caster Size/Type	6"/Locking
Motor	Full Electric w/ Hi Lo
Reverse Trendelenburg	0 – 12°
Backrest Section Tilt	0 – 74°
Thigh Section Tilt	0 – 31.5°
Mattress Type	Castellated Foam
Mattress Dimensions	75.8"(L) x 33.6"(W) x 6"(H)