

Thomashilfen is a family owned company in the 3<sup>rd</sup> generation passionate about providing safe and comfortable products for special needs. We are committed to developing and delivering products that help families and their children for more Quality of Life.

Let us **move** you!

“The Defender Reha allows children with special needs 360° of impact protection. The head, trunk and thighs supports, combined with our unique accessories, offer many positioning options. The seat allows for growth while remaining narrow enough for a second car seat as well as another child to fit in the back seat.”

Positioning-options  
**INCLUSIVE**



Cup/snack holder swivels and can easily be installed on either armrest



Smooth up-front harness adjuster



Lateral trunk supports for an upright, stable sitting position

Ultra plush seat pad made with durable premium fabrics and extra foam for maximum comfort is removable for machine washing



9-position, 1-hand height adjustable headrest



Head support holds the child's head safely in place



Chest clip and soft harness covers included



Lateral side supports are removable to adjust seat width for plenty of room for the child to grow with the seat

The no-threaded 5-point harness provides additional support in Booster Mode for children weighing 40 to 120 lbs. Small crotch pad included.

Deep, supportive headrest, backrest and seat base for 360° impact protection

Thomashilfen

**Defender Reha US**

The Positioning 2-in-1 Special Needs CarSeat

## Accessories & Measurements



Footrest, short with footrest adapter and footrest pad



Footrest, long with footrest adapter



Seat wedge, below (+ 10°)



Table = playtray and additional impact shield



Abduction block



Crotch pad, large

Defender Reha AND all accessories exceed US FMVSS 213 and Canadian CMVSS 213 standards



Measurements & weights	Defender Reha US
Seat depth	13.5"
Seat width*	8.5 - 12"
Seat to head**	23.5 - 30.5"
Shoulder height (child)***	10 - 19.5"
Lower leg length	(7") 11 - 17"
Overall dimensions (w x h x d)	18.25 x 26.5 x 20"
Weight of seat	22 lb
<b>HARNESSES MODE</b>	
Weight of child	26.5 - 65 lb
Height of child	27 - 57"
<b>BOOSTER MODE</b>	
Weight of child	40 - 120 lb
Height of child	43.4 - 57"

\* To reach max. seat width remove additional hip guide pads from additional seat pad  
 \*\* To reach max. back height remove additional seat pad and foam pad under seat cushion  
 \*\*\* To reach max. shoulder height remove additional seat pad and foam pad under seat cushion

Product & accessories	item code
Defender Reha US	886
Footrest, short (depth 4.7")* (tiltable and adjustable in angle)	806
Footrest pad (height 4") — (for 806)	834
Footrest, long (depth 7.1")* (tiltable and adjustable in angle)	805
Footrest adapter + Base plate*	821
Seat wedge, below (+ 10°), black	810
Table, black	831
Crotch pad, large, black**	832
Abduction block, black	833

\* The footrest adapter must be ordered to use footrest  
 \*\* We recommend using larger crotch pad only in booster mode.



# DEFENDER REHA US

The Special Needs 2-in-1 Car Seat Solution  
 Helping Families Enjoy Safe Traveling