



A



A

Deluxe Hip-High Chair

Hip-High Chair Steel Construction

Item #	Product Description	UOM
17100	Maple Wood Frame Construction	1/cs

- Designed for post-hip surgery residents.
- Allows for sitting without bending at the hip.
- Armrests make sitting down or pushing up to get out from the chair safe and easy.
- Attractive vinyl seat and backrest are soft, washable, comfortable and water resistant.
- Constructed of solid maple hardwood—not a plywood veneer.
- The footrest has three height adjustments to provide maximum comfort and an ergonomic fit. (Figure A)

specifications

Floor to Seat Height	26"
Width Between Arms	19"
Floor to Arm Height	31"
Floor to Footrest Height (adjustable)	5.5", 7.5", 9.5"
Overall Dimensions	24" (W) x 40" (H) x 25" (D)
Weight	33 lbs.
Carton Shipping Weight	44 lbs.
Weight Limit	300 lbs.
Warranty	1 Year on Upholstery, Limited 2 Year Warranty on Frame

Item #	Product Description	UOM
17100-BV	Steel Brown-Vein Frame Construction	1/cs

- Designed for post-hip surgery residents.
- Attractive, durable brown-vein finish is easy to maintain and less institutional looking. (Figure A)
- Contoured armrests are comfortable and durable.
- Attractive vinyl seat and backrest are soft, washable, comfortable and water resistant.

specifications

Floor to Seat Height	26"
Width Between Arms	17.5"
Floor to Arm Height	32"
Floor to Footrest Height	6"
Overall Dimensions	22" (W) x 40.75" (H) x 29" (D)
Weight	37 lbs.
Carton Shipping Weight	47 lbs.
Weight Limit	300 lbs.
Warranty	1 Year on Upholstery, Limited Lifetime Warranty on Frame