

<b>MAXIMUM SPEED</b>	15 MPH (varies based on riders weight and driving terrain)
<b>TRAVEL RANGE</b>	40 MILES (varies based on riders weight and Driving terrain)
<b>Bluetooth</b>	Yes, Advanced Bluetooth system with sound
<b>MOTOR</b>	500 watt Trans-Axel Motor
<b>WEIGHT CAPACITY</b>	350 LBS
<b>BATTERIES</b>	48 Volt, 20 AH
<b>BRAKES</b>	Front and Rear Drum Brakes
<b>BRAKE LEVERS:</b>	Front and Rear
<b>STORAGE</b>	Rear Box Storage & Under Seat Storage both secure lockable
<b>VARIABLE TWIST THROTTLE</b>	Yes
<b>SEAT SIZE</b>	14 " Wide and 23" Deep
<b>SEAT BASE HEIGHT</b>	From Floor Height 30
<b>SMART CHARGER</b>	Yes
<b>BAG HOOK</b>	Yes
<b>PARKING BRAKE</b>	Yes
<b>Front Storage Area</b>	Two Storage Bins
<b>REVERSE SWITCH</b>	Yes
<b>ELECTRIC HORN</b>	Yes
<b>GUAGES</b>	Yes, Speedometer with Battery Gauge
<b>3 Speed Settings</b>	Low, Medium and High
<b>CLIMBING:</b>	12 degrees
<b>LIGHTING PACKAGE</b>	Turn Signals, Brake Lights, High/Low Beam Headlight, Rear Light
<b>KEYLESS KEY FOB</b>	Yes, Power on/off
<b>ANTI-THEFT ALARM SYSTEM</b>	Yes, Keyless Remote
<b>BATTERY INDICATOR GAUGE</b>	Yes
<b>TIRE SIZE</b>	Front and Rear Tire Size 15"/3"
<b>REAR REFLECTORS</b>	Yes
<b>ALUMINUM RIMS</b>	Yes
<b>SCOOTER LENGTH</b>	70"
<b>SCOOTER WIDTH</b>	28"
<b>SCOOTER HEIGHT</b>	40"
<b>DECK HEIGHT</b>	12"
<b>GROUND CLEARANCE</b>	6.5"
<b>WHEEL BASE</b>	52"
<b>FRONT SUSPENSION FORK</b>	Yes
<b>REAR SUSPENSION</b>	Two Heavy Duty Shocks
<b>SCOOTER WEIGHT</b>	209 LBS
<b>WARRANTY</b>	3 YEAR LIMITED