

Table of contents

Back supports

Cover technology	3
Hardware technology	4 - 5
Back support technology.....	6
General use back supports	7
Positioning back supports	8 - 14
Positioning back supports with laterals	15 - 16

Cushions

Cushion features	19
General use	20 - 21
Skin protection	22
Positioning	22
Skin protection and positioning.....	23
Adjustable skin protection.....	23 - 27
Adjustable skin protection and positioning	27 - 31
Custom	32 - 33
Non coded cushions	33 - 34
ROHO® cushion sizing guide	34

Supports & accessories

BodiLink®	37 - 43
Lower extremity	44 - 48
Upper extremity	49 - 52
Solid seat support	52 - 53
Replacement covers.....	53 - 54
Additional accessories.....	55 - 58

References

Comfort Elements

E2611/E2612
General use back



Features

- Easily folds with wheelchair for transportation and storage
- Adjustable tension strapping system can accommodate for kyphosis and provide lumbar support
- 1” foam provides a comfortable back support
- Back securely installs directly over wheelchair canes
- Steel D-Rings provide superior strength and durability
- Cover options: Comfort-Tek or Stretch Air
- Warranty: 5 years back and hardware, 6 months for cover

Dimensions		Weight capacity	
Width	15” – 21”, 22” – 25”	All sizes	300 lbs
Length	16”		

Comfort Elements part number
ELB-BSXXXXW16L
<i>*Length (L) refers to the actual size dimension from bottom to top edge</i>
XXXX to be replaced XXXX = width in inches

Comfort Unite 1

E2613/E2614
Positioning back



Features

- 2.25" of shell contour.
- Adjustable moldable stays to conform and support the back.
- Keyhole (in the shell) and stud strap (on the cover) attachments for secure fit.
- Lightweight, durable aluminum shell providing pelvic and scapular relief.
- Hardware quick install: easy adjustment for depth, height and angle while

Comfort Unite 2

E2615/E2616
Positioning back, Posterior



Dimensions		Weight capacity	
Width	12" – 24", (2" increments)	≤ 18" wide	250 lbs
Length	12" – 22", (2" increments)	20" wide	300 lbs
Contours	4"	22" wide	350 lbs
		24" wide	400 lbs

Standard sizes							
	12" W	14" W	16" W	18" W	20" W	22" W	24" W
22" L			16 x 22	18 x 22	20 x 22	22 x 22	24 x 22
20" L		14 x 20	16 x 20	18 x 20	20 x 20	22 x 20	24 x 20
18" L	12 x 18	14 x 18	16 x 18	18 x 18	20 x 18	22 x 18	24 x 18
16" L	12 x 16	14 x 16	16 x 16	18 x 16	20 x 16	22 x 16	
14" L	12 x 14	14 x 14	16 x 14	18 x 14	20 x 1		

Back supports with included trunk supports

- Lateral trunk supports (LTS) are constructed with moldable stainless steel and protective foam.
- Maximum distance between two lateral supports is equal to the width of the wheelchair. Minimum distance between two lateral supports is 8” less than the width of the wheelchair.
- Lateral trunk supports have up to 8” of width adjustment (4” per pad) and 1.5” depth adjustment each.
- Lateral trunk supports are available in four different sizes: extra small, small, medium, and large.
- Fixed lateral hardware comes standard.
- Extra lateral hardware links available to provide optimal contact with the user.

Lateral trunk support options

4” D x 3” L (LAT-4D3L-COMFORT-TEK or LAT-4D3L-STRETCH-AIR)

5” D x 4” L (LAT-5D4L-COMFORT-TEK or LAT-5D4L-STRETCH-AIR)

6” D x 4” L (LAT-6D4L-COMFORT-TEK or LAT-6D4L-STRETCH-AIR)

7” D x 5” L (LAT-7D5L-COMFORT-TEK or LAT-7D5L-STRETCH-AIR)

Lateral trunk support option can be selected by adding the (part number) to the Comfort Unite 1 LTS, Comfort Relief LTS or Comfort Embrace LTS part number.

Comfort Unite 1 LTS



Weight capacity	
10” - 18” wide	250 lbs
20” wide	300 lbs
22” wide	350 lbs
24” wide	400 lbs

Comfort Support Pro 7 Series

E2602
General use



Features

- Non-contoured molded foam surface allowing for easier transfers.
- Coccyx relief to assist in offloading bony prominences.
- Kwik strap loops around wheelchair canes to securely attach cushion during transfers.
- Warranty: 6 months cover/ 5 years cushion

Dimensions		Weight capacity	
Seat thickness	3.5"	All sizes	750 lbs

Standard sizes					
	22" W	24" W	26" W	28" W	30" W
18" D	22 x 18	24 x 18	26 x 18	28 x 18	30 x 18
20" D	22 x 18	24 x 20	26 x 20	28 x 20	30 x 20

Support Pro 7 Series part number	
53- WWDD -B	Without QuadraGel
53G- WWDD -B	With QuadraGel
WW & DD to be replaced WW = width & DD = depth in inches	

Comfort Ascent

Comfort Saddle 7 Series



Dimensions		Weight capacity	
Seat thickness	4.5"	All sizes	750 lbs
Inferior thigh support thickness	3.5"		
Abductor	5"		

Standard sizes	22" W	24" W	26" W	28" W	30" W
18" D	22 x 18	24 x 18	26 x 18	28 x 18	30 x 18
20" D	22 x 18	24 x 20	26 x 20	28 x 20	30 x 20

E2608
Skin protection and positioning

Features

- Lateral and medial contours promote postural alignment and stability.
- Optional QuadraGel pack to redistribute pressure for skin protection.
- Coccyx relief to assist in offloading bony prominences.
- Kwik strap loops around wheelchair canes to securely attach cushion during transfers.
- Warranty: 6 months cover/ 5 years cushion

Saddle 7 Series part number	
53S- WWDD -B	Without QuadraGel
53SG- WWDD -B	With QuadraGel
WW & DD to be replaced WW = width & DD = length in inches	

ROHO Single Compartment

E2622 / E2623
Adjustable skin protection



Single Compartment part numbers

Cushion style	Standard cover	Heavy duty cover
High profile	1R WWDDC	1R WWDDHD
Mid profile	1R WWDDMPC	-
Low profile	1R WWDDLPC	-

	Sensor ready standard cover	Sensor ready heavy duty cover
High profile	1R WWDDH -CA-SR	1R WWDDH -CH-SR
Mid profile	1R WWDDM -CA-SR	R WWDDM -CH-SR

WWDD to be replaced | **WW** = width & **DD** = depth in cells

Features

- Adjustable to shape, weight, and changes in an individual’s position over time.
- Hydrostatic distribution of load across the entire contact area.
- Provides immersion and envelopment for pressure redistribution.
- SR: Cushion is Sensor Ready™ for Smart Check, which can be easily attached and detached by means of quick disconnect.
- No weight limit if appropriately sized to the individual.
- Warranty: 6 months cover / 2 years cushion

Smart Check part number

SMARTCHECK-01

Dimensions

References

Weight conversion

Pounds	Kilograms	Pounds	Kilograms
100 lbs	45.4 kg	350 lbs	158.8 kg
150 lbs	68.0 kg	400 lbs	181.4 kg
200 lbs	90.7 kg	450 lbs	204.1 kg
250 lbs	113.4 kg	550 lbs	249.5 kg
300 lbs	136.1 kg	750 lbs	340.2 kg
330 lbs	149.7 kg		

Size conversion

Inches	Centimeters	Inches	Centimeters
10"	25.4 cm	22"	55.9 cm
12"	30.5 cm	24"	61.0 cm
14"	35.6 cm	26"	66.0 cm
16"	40.6 cm	28"	71.1cm
18"	45.7 cm	30"	76.2 cm
20"	50.8 cm		

© 2025 Permobil, Inc. Products may be covered by one or more U.S. and foreign patents and trademarks, including but not limited to: THE ADAPTOR®, AGILITY®, AIRLITE®, CONTOUR SELECT®, DRY FLOATATION®, HIGH PROFILE®, Hybrid Elite SR®, LOW PROFILE®, MINI-MAX®, QUADTRO SELECT®, ROHO®, Sensor Ready™, Smart Check®, Comfort-Tek™, Compass®, Curve®, QUADROGEL®, and Elements®.

Product information is changed as needed; current product information is available at permobil.com.

ROHO, Inc. and The Comfort Companies, LLC are part of Permobil, Inc.

GlideWear™ is a trademark of Tamarack Habilitation Technologies, Inc. BOA® is a trademark of BOA Technology Inc. Zytel® is a trademark of Dupont Polymers, Inc. Vicair®, SmartCells® and Vector® are trademarks of Vicair B.V.

