

Cirrus Plus HD

HCPCS: K0824-1, K0824




FIGURE A

- Foldable, lightweight frame for easy portability (Figure A)
- Steel double cross-brace frame
- Adjustable-tension back upholstery
- Armrests are padded, adjustable and removable
- Adjustable-length controller mount
- Adjustable-tension back upholstery
- Swing-away footrests with heel loops
- Optional elevating leg rests (Item #CPELR) sold separately
- Optional 2" seat depth extension (Model #STDS1818N) sold separately
- Positioning belt included
- Padded back upholstery with storage pouch and cushioned seat
- Programmable controller
- Flat-free, non-marking tires
- Easy free-wheel operation
- Manual wheel locks
- General foam seat cushion included
- Anti-tip wheels

PRODUCTS

Description

CPN22FBA	CirrusPlusPowerChair,22in FlpBkArm,1/cs
----------	---

SPECIFICATIONS

Description

Base Weight	100 lbs. (22"), 105 lbs. (24")
Batteries (pair)	12V x 36AH (22"), 12V x 50AH (24")
Battery Weight (pair)	56 lbs. (22"), 66 lbs. (24")
Brakes	Electromagnetic
Caster Wheels	8" x 2" Solid
Charger	4A Offboard
Climbing Angle	6°
Controller	PG VSi 50A
Dimensions	42.5" x 29.5" (22"), 46" x 31" (24")
Drive Wheels	12.5" x 2.25" Flat-free
Floor to Seat Height	18.5"
Ground Clearance	4"
HCPCS	K0824
Maximum Range	15 miles

SPECIFICATIONS	Description
Motors (2)	24V x 320W x 4600 rpm
Seat Dimensions	22" x 16", 24" x 19.5"
Top Speed	5 mph
Total Weight (excl. hand control)	156 lbs. (22"), 159 lbs. (24")
Turning Radius	35"
Warranty on Batteries	12 Months
Warranty on Electronic Controller/Drive Train Components	24 Months
Warranty on Frame	Lifetime
Weight Capacity	450 lbs
Weight Capacity (13046 13046-HD)	400 lbs. (22"), 450 lbs. (24")

PARTS	Description
CPN8006	Pivot Release Brkt.CPN,1/ea, EA
CPN8014	Lower Arm Frame,CPN,1/ea, EA
CPN8004	Upper Arm Frame,CPN,1/ea, EA
CPN8008R	Front Seat Saddle,Rt.CPN,1/ea, EA
CP8005-20	Adj. Tension Back Upholstery, CPN20, 1ea, EA
CP8005-18	Adj. Tension Back Upholstery, CPN18, 1ea, EA
CP8005-16	Adj. Tension Back Upholstery, CPN16, 1ea, EA
CP8006-20	Padded Back Upholstery Cover, CPN20, 1ea, EA
CP8006-16	Padded Back Upholstery Cover, CPN16, 1ea, EA
CP8006-22	Padded Back Upholstery Cover, CPN22, 1ea, EA
CP8006-18	Padded Back Upholstery Cover, CPN18, 1ea, EA
CP8019-20	Link 20",CirrusPlus1/ea, EA
CPN8005	Rear Seat Saddle,CPN,1/ea, EA
CP8019-16	Link 16",CirrusPlus1/ea, EA
CP5005C	BatteryBasketBackClip,Cir.Plus,1/ea, EA
CP5005CF	BatteryBasket2pcFrontClip1/ea, EA
CP8015-20	Foldable,20" Bat Basket,Cirrus,1/ea, EA
CP8015-18	Foldable,18" Bat Basket,Cirrus,1/ea, EA
CP8015-16	Foldable,16" Bat Basket ,Cirrus,1/ea, EA
CP8035R	Brake,RightCirrusPlus1/ea, EA