

## Cobra GT4

HCPCS:



- New, innovative design for a distinctive, stylish appearance
- One-touch, infinite tiller angle adjustment allows the user to easily position the tiller into the ideal driving position
- Ergonomic delta control handles makes it easier for users with limited strength and dexterity to operate the scooter
- Large, air-filled, low profile tires and full suspension for smoothest outdoor ride quality
- "Blade" style alloy wheel rims
- Efficient, extra bright LED lighting (Tail Lights and Turn Signals)
- Choice of Full size, Sport Style, Captain's Seats in 20" or 22" widths
- Seats are height adjustable, swivel and feature a fold-down, reclining backrest, and auto-style slides allow forward or backward adjustment
- Padded armrests feature easy turn cylinders for angle adjustment
- Positioning belt included
- "High/Low" switch for instant speed control
- Anti-tip wheels for greater safety
- Easy, freewheel mode

**PRODUCTS****Description**

COBRAGT422CS	COBRA GT4 SCOOTER WITH 22" CAPTAIN SEAT
--------------	---

**SPECIFICATIONS****Description**

Base Weight	191 lbs
Batteries (pair)	12V x 75AH
Battery Weight (pair)	94 lbs
Brakes	Electromagnetic
Charger	5A Offboard
Climbing Angle	8°
Controller	Dynamic Rhino II 160A
Dimensions	63" (L) x 28.8" (W)
Floor to Seat Height	27" - 30"
Front Wheels	14" x 14" Pneumatic
Ground Clearance	5.5"
Maximum Range	22 miles
Motor	24V x 950W x 5200 rpm
Rear Wheels	14" x 14" Pneumatic
Seat Dimensions	20" (W) x 20" (D)   22" (W) x 20" (D)
Seat Weight	55 lbs
Top Speed	10 mph

SPECIFICATIONS	Description
Total Weight (excl. hand control)	340 lbs
Turning Radius	82"
Warranty on Battery	12 Month Limited
Warranty on Electronics	24 Months Limited
Warranty on Frame	Limited Lifetime
Weight Capacity	450 lbs

PARTS	Description
29002W7011	AxleKey,Cobra,Prowler,PhoenixHD,Vent,1/e, EA
C31-029-00101	VelcroStripForBattery,Cobra/KingCobra,1e, EA
C26-051-00501	Controller Cover, Cobra/King Cobra, 1/ea, EA
C10-029-00400	Controller, Rhino 2- 160A, Cobra, 1/ea, EA
S10-062-00100	Controller Holder, Cobra, 1/ea, EA
S18-062-00400	Wire Harness Assembly, Cobra, 1/ea, EA
S01-062-00200	Rear Frame, Cobra, 1/ea, EA
C04-062-00200	Shock Absorber, 1,050lbs, Cobra, 1/ea, EA
C09-062-02100	Transaxle Only, Cobra, 1/ea, EA
S09-062-00400	Motor & Transaxle Assembly, Cobra, 1/ea, EA
C09-062-02000	Motor Only, 5200RPM / 950W,Cobra, 1/ea, EA
C09-022-02300	MotorBrush,6.75mmX17mmX18.5mm,Cobra,1/ea, EA
C09-022-02700	Carbon Brush Cap, 7x17mm, Cobra, 1/ea, EA
C09-062-01100	Electromagnetic Brake, Cobra, 1/ea, EA
C09-062-00903	Freewheel Lever, Cobra/King Cobra, 1/ea, EA
S06-062-00700	RearWheelAssembly,8"Rim,Right,Cobra,1/ea, EA
S06-062-00600	RearWheelAssembly,8"Rim,Left,Cobra, 1/ea, EA
C06-022-02600	Brake Pad, 14"Rim,Cobra,1/ea, EA
S08-062-00577	Rear Shroud, Grey, Cobra, 1/ea, EA
C14-062-00201	RearSignalLightRight,Cobra/KingCobra,1ea, EA