



Yes, you can.®


CS3™

Bed featured with assist bar,
traditional bed ends in Biltmore Cherry
and Solace® mattress



It's time to change your definition of value. The CS3 was developed for those who require a bed that is stable throughout the travel range and mobile in the low position.

CS3™



Tool-less Assembly

Fast connect system allows for faster setup and easier alterations.



Hand Control Pendant

Carroll CS beds feature a rugged Hand Control, with bright, universal symbols, allowing residents and staff to easily operate the comfort controls of the bed.



Mobility-On-Demand

The CS3 features the Maxi-Lock floor lock system, securing the bed to the floor in all elevated positions. Mobility is enabled in the low position.



Slat Deck

Exclusive arched slat deck design allows for improved infection control and reduces wear and tear on mattresses.

CS3 Specifications

- 8.5" – 27" Height Range
- 450 lbs. Weight Capacity
- Locking Casters
- I.V. Rod Holders
- Meets UL60601-1
- Meets CAN/CSA 22.2 No. 601.1-M-90
- Meets CAN/CSA 22.2 No. 601.1-2-38:03
- Tool-less Assembly
- Universal Deck Length (76" or 80")
- Four DC Actuators
- Ergonomically Designed, Rugged, Water Resistant Pendant
- Optional Attendant Control Panel
- Fixed Heel Lift

CS3 Accessories

- Hi-Low Lockout
- CPR Release
- 4" Bed Extender

Invacare Continuing Care

www.invacare-ccg.com

info_iccg@invacare.com

USA

One Invacare Way
Elyria, Ohio 44035
USA
800-668-2337

Canada

994 Hargrieve Road
London, Ontario
N6E 1P5 Canada
800-668-2337

Trademarks are identified by the symbol™ and®.
All trademarks are owned or licensed to Invacare Corporation or its subsidiaries unless otherwise noted.

© 2009 Invacare Corporation. Form 08.28.1/09



Yes, you can.®