

## Ultimate Foot Circulator w/ Remote Electrotherapy Pain Relief



A simple, drug-free solution for relaxing stiff muscles and increasing blood circulation to improve and facilitate muscle performance.

The **AccuRelief™ Ultimate Foot Circulator** is an advanced electrical muscle stimulator for home, work or travel. This powerful device delivers prescription strength relief by relaxing stiff muscles, reducing swollen feet and ankles while increasing blood circulation. Electrical muscle stimulation delivers pulses into the muscle tissue to assist the muscle to contract and relax while sending a signal to the brain to help the body's natural contraction ability.

The Ultimate foot circulator features a three channel electrical stimulator with a stepper function, two gel pads and remote control for ease of use. With up to 99 levels of adjustable intensity and 15 therapeutic phases throughout the cycle time, finding a comfortable level of relief is simple and quick.

The AccuRelief Ultimate Foot Circulator is FDA approved for use without a prescription.

### Product Features:

- Uses EMS technology (Electrical Muscle Stimulation) for safe and effective pain relief.
- FDA cleared for over-the-counter sale (no prescription required).
- Pre-programmed with fifteen therapeutic phases adjusting throughout the cycle time for quick, safe and effective relief.
- Powerful with 99 levels of adjustable intensity – increase or decrease stimulation with remote control.
- Comes complete with remote control, four re-usable electrodes, pair of lead wires, two AAA batteries (for remote), AC Adaptor and instruction manual.

### Technical Specifications:

- Power sources: Adaptor (included) or batteries (optional) for foot circulation device.
- Waveform: Symmetrical Biphasic square pulse with zero DC current
- Frequency: 10Hz – 36Hz
- Pulse Width: 250us for foot stimulation; 200us for body stimulation
- Intensity Levels: 0-99 adjustable levels
- Treatment Time: 1-99 min
- Channels: Three, Isolated, (1) for Foot Stimulation, (2) for Body Stimulation
- Size: 16.4" x 4.09" x 16.82"
- Weight: 4.21 lbs
- Warranty: One year warranty form date of purchase
- Replenishment kit available (ACRL-0001)

| Item #    | Description  | Ship Dimensions          | Ship Weight | UPC Code     | Qty |
|-----------|--|--------------------------|-------------|--------------|-----|
| ACRL-5500 | AccuRelief Ultimate Foot Circulator                      | 19.49" x 12.20" x 17.32" | 12 lbs      | 023601255006 | 2   |
| ACRL-0001 | AccuRelief Universal Supply Kit, (4) packs of electrodes | 10.81" x 3.31" x 8.87"   | 6.196 lbs   | 023601200013 | 6   |

[www.accurelief.com](http://www.accurelief.com)

This device is FDA cleared for over-the-counter sale for pain relief associated with sore and aching muscles.

