



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Company: _____

Ship To: _____

Account Number: _____ Date: _____

Address: _____

Quote PO Number: _____

City: _____ State: _____

RTS: _____

Zip: _____ Phone: _____

Marked For: _____

Email/Fax: _____

Select if chair quote has: Additional Ki Products with Different Quote/Order Number **Quote/Order Number:** _____

Catalyst 5Vx K0005 _____ **\$1895** Effective 08/28/2020
 300 lb weight capacity

HCPCS codes provided are based on code verification by the PDAC or our interpretation of the code definitions. Proper use of HCPCS codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

Frame

Seat Width

- 14" _____ N/C
- 15" _____ N/C
- 16" _____ N/C
- 17" _____ N/C
- 18" _____ N/C
- 19" _____ N/C
- 20" E2201 _____ \$500
- 21" E2201 _____ \$500
- 22" E2201 _____ \$500

Seat Depth

- 14" _____ N/C
- 15" _____ \$90
- 16" _____ N/C
- 17" _____ \$90
- 18" _____ N/C
- 19" _____ N/C
- 20" E2203 _____ \$500

Weight Capacity

- 300 lb _____ Std
- 350 lb (Heavy Duty) _____ \$500
Requires minimum 6" caster height

Backrest

Back Post / Height Range / Back Height

Back Post	Price	Height Range	Back Height			
<input type="checkbox"/> Standard	N/C	<input type="checkbox"/> Short	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11.5
		<input type="checkbox"/> Medium	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15.5
		<input type="checkbox"/> Tall	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19.5
<input type="checkbox"/> Angle Adjustable w/ Push Handle	\$185	<input type="checkbox"/> Short	<input type="checkbox"/> 14	<input type="checkbox"/> 15	<input type="checkbox"/> 16	
		<input type="checkbox"/> Medium	<input type="checkbox"/> 16	<input type="checkbox"/> 17	<input type="checkbox"/> 18	
		<input type="checkbox"/> Tall	<input type="checkbox"/> 18	<input type="checkbox"/> 19	<input type="checkbox"/> 20	

Transit

- Transit System _____ \$320
Not available with Removable Seat Pan or 350 lb Weight Capacity
- Omit Transit _____ N/C

Cross Brace

- Standard Cross Braces _____ Std
- Seating System Cross Braces _____ N/C
Allows for easy attachment of seating system mounting hardware. Not available with Seat Sling.

Seat Option

- Removable Aluminum Seat Pan E2231 _____ \$265
Not available with Transit System. Pan depth 1/2" shorter than Seat Depth.

Wheelchair Drive Ready

Chair will ship from Ki Mobility with adaption kit installed for selected Alber system. Order system separately from Alber USA. Select Rear Seat Height (full profile selections only) in drive wheel size on Rear Seat Height chart (page 5). Requires Rear Anti-Tips. Not available with Height Adjustable Flip Back T-Arm.

- e-motion Ready _____ N/C
- twion Ready _____ N/C

Standard Back Post

Select option below only if a Standard Back Post option is selected on chart to the left

- Straight with Push Handle _____ Std
- Straight without Push Handle _____ N/C
- 8° Bend with Push Handle _____ \$25
Not available with Short Back Height Range
- Straight with Fold Down Push Handle _____ \$125
Not available with 350 lb Weight Capacity

Push Handle Options

Not available with Straight Back Post with Fold Down Push Handle or 350 lb Weight Capacity

- Low Profile Push Handle _____ \$65
2" shorter than standard push handles

Backrest continued

Back Upholstery

If selecting both upholstery and Axiom Back, the upholstery will be installed.

- Standard _____ Std
Not available with Short or Medium Back Height Ranges on 20" - 22" Seat Widths
- Tension Adjustable E2611
 - Standard Cover _____ \$315
 - Air Mesh Breathable Cover _____ \$390
- Omit Back Upholstery _____ - \$15

Military Patch

- Attached only if ordered with Tension Adjustable Back Upholstery*
- US Military Service Patch _____ \$90
 - Army Navy Marines Air Force Coast Guard
 - US Purple Heart Patch _____ \$90
 - Purple Heart with Military Service Patch Combo _____ \$145
 - Army Navy Marines Air Force Coast Guard

Footrest / Frame Type

Frame Type

- Swing Away _____ N/C
- Super Low - Swing Away _____ N/C
- Fixed Front _____ N/C

Minimum Front Seat Height with Front Mount Hangers						
Hanger Length	Short		Medium		Long	
Hanger Bend	70°	80/90°	70°	80/90°	70°	80°
Minimum Front Seat Height with Standard Selection	13.5"	13.5"	14"	15"	16"	17"
Minimum Front Seat Height with 3" Extension	14.5"	16"	16.5"	18"	18.5"	20"

Hanger / Front Frame Bend

Frame	Hanger	60°	70°	80°	90°	Retail
Swing Away	Extension Tube Mount	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/> ¹		Std
	Front Mount		<input type="checkbox"/>	<input type="checkbox"/> ¹	<input type="checkbox"/> ¹	N/C
	PRO ELR				<input type="checkbox"/>	\$295
	Omit Hangers				<input type="checkbox"/>	-\$20
Fixed	Extension Tube Mount		<input type="checkbox"/>	<input type="checkbox"/> ¹		N/C

¹ Not available with 7" or 8" casters

Footrest / Extension / Seat to Footrest Length

Select Extension and Seat to Footrest Length in row corresponding to desired Footrest. Footrest Length of at least 2.5" shorter than Front Seat Height recommended. Extension Tube Hangers and Fixed Front Frame measurements are based on 70° Hanger / Front Frame Bend.

Footrest	Extension	Seat to Footrest Length											
Extension Tube Hangers (60°, 70°, 80°)	<input type="checkbox"/> Short	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15			
	<input type="checkbox"/> Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17			
	<input type="checkbox"/> Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5			
Front Mount Hangers ¹	<input type="checkbox"/> Short (70°, 80°, 90°)	Standard Selections									Requires 3" Extension for \$35		
		<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7					<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11
	<input type="checkbox"/> Medium (70°, 80°, 90°)	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5					<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5
		<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8			<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13
<input type="checkbox"/> Long (70°, 80°)	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5			<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	
	<input type="checkbox"/> 5	<input type="checkbox"/> 5.5	<input type="checkbox"/> 6	<input type="checkbox"/> 6.5	<input type="checkbox"/> 7	<input type="checkbox"/> 7.5	<input type="checkbox"/> 8	<input type="checkbox"/> 8.5	<input type="checkbox"/> 9	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	
		<input type="checkbox"/> 9.5	<input type="checkbox"/> 10	<input type="checkbox"/> 10.5	<input type="checkbox"/> 11	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5
PRO ELR Low Pivot Setting	<input type="checkbox"/> Short	12 - 17.5											
	<input type="checkbox"/> Long	15 - 21											
PRO ELR High Pivot Setting	<input type="checkbox"/> Short	10.5 - 15.5											
	<input type="checkbox"/> Long	13 - 19											
Fixed Front Frame (70°, 80°)	<input type="checkbox"/> Short	<input type="checkbox"/> 11.5	<input type="checkbox"/> 12	<input type="checkbox"/> 12.5	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15				
	<input type="checkbox"/> Medium	<input type="checkbox"/> 13	<input type="checkbox"/> 13.5	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17			
	<input type="checkbox"/> Long	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5			

¹ See Minimum Front Seat Height chart above.

Footrest / Frame Type continued

3" Extension

Available only with Front Mount Hangers

- 3" Extension \$35
- 3" Extension (Provided as Parts) \$35


Footplate

- Composite N/C
Not available with Front Mount Hanger or PRO ELR
- Composite Angle Adjustable K0040 \$185
Not available with Front Mount Hanger or 14" Seat Width
- Aluminum Angle Adjustable K0040 \$185
- One Piece Flip Up Angle Adjustable K0040 \$190
Not available with Front Mount Hanger or PRO ELR

Hanger Release

If no selection is made Classic Hanger Release will be used

- Classic Std
- 4-Way - In line N/C
See reference A to the right. Not available with PRO ELR.
- 4-Way - Out N/C
See reference B to the right. Not available with PRO ELR.



Heel Loop

Not available with One Piece Flip Up Angle Adjustable Footplate

- Heel Loops E0951 \$30
- Adjustable Heel Loops E0951 \$50
- Omit Heel Loops N/C

Calf Strap

- Velcro Adjustable K0038 \$40
- Padded Velcro Adjustable K0038 \$68
- Bodypoint Padded AeroMesh K0038 \$102

Residual Limb Support

Residual Limb Support E1020

- Right and Left \$790
- Left \$395 Right \$395
- Omit Left Hanger -\$10 Omit Right Hanger -\$10

Cover

- Smooth Air Mesh Smooth Air Mesh

Pad Size

- | | | | |
|---------------------------------|----------------------------------|---------------------------------|----------------------------------|
| <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 | <input type="checkbox"/> 5 x 8 | <input type="checkbox"/> 9 x 14 |
| <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 | <input type="checkbox"/> 7 x 10 | <input type="checkbox"/> 11 x 10 |
| <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 | <input type="checkbox"/> 9 x 10 | <input type="checkbox"/> 11 x 14 |

Front Seat Height and Casters

Caster Size / Front Frame Type / Front Seat Height

Select Front Seat Height in row corresponding to desired caster size.

Note: Selections in first three Front Seat Height columns use Super Low Front Frame. Super Low Front Frame not available with Fixed Front Frame.

Caster Size	Front Seat Height													
	Uses Super Low Front Frame													
4"	<input type="checkbox"/> 13.5 ¹	<input type="checkbox"/> 14	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15 ²	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20
5"	<input type="checkbox"/> 14.5	<input type="checkbox"/> 15	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	
6"	<input type="checkbox"/> 15.5	<input type="checkbox"/> 16	<input type="checkbox"/> 16.5	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21		
7"	<input type="checkbox"/> 17	<input type="checkbox"/> 17.5		<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21				
8"	<input type="checkbox"/> 17.5	<input type="checkbox"/> 18	<input type="checkbox"/> 18.5	<input type="checkbox"/> 19	<input type="checkbox"/> 19.5	<input type="checkbox"/> 20	<input type="checkbox"/> 20.5	<input type="checkbox"/> 21						

¹ Front Seat Height not available with 4 x 1" Poly Caster

² Super Low Front Frame required when using 4 x 1" Poly Caster

Caster

Caster Type	Price	Caster Size				
		4"	5"	6"	7"	8"
.75" Roller Blade	N/C	<input type="checkbox"/>				
.75" Lighted Roller Blade	\$70	<input type="checkbox"/>				
1" Poly	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
1" Poly Aluminum	\$160	<input type="checkbox"/>	<input type="checkbox"/>			
1" Pneumatic	\$85			<input type="checkbox"/>		<input type="checkbox"/>
1.5" Poly	\$50	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
1.5" Soft Roll Aluminum	\$175	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
2" Poly	\$70			<input type="checkbox"/>		
2" Pneumatic w/ Foam Insert	\$225					<input type="checkbox"/>

Maximum Seat Dump Reference

(Front Seat Height minus Rear Seat Height)

Seat Depth	14"	15"	16"	17"	18"	19"	20"
Maximum Dump Standard Caster Housing	1.5"	2"	2"	2"	2.5"	2.5"	2.5"
Maximum Dump Performance Caster Housing	3"	3"	3"	3.5"	3.5"	4"	4"

Depending on configuration, chairs with Fixed Front Frame Type may achieve up to an additional .75"

Front Seat Height and Casters continued

Fork Stem

- Best Fit _____ Std
Ki Mobility will pick Fork and Stem combination to optimize performance based on Caster Size and Front Seat Height
- Standard Fork Stem _____ N/C
- 3/4" Longer Fork Stem _____ N/C
- 1 1/2" Longer Fork Stem _____ N/C

Caster Housing

- Standard _____ N/C
- Performance _____ \$40

Fork Options

- Frog Legs Ultra-Sport Fork - Black _____ \$425
Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.
 - User Weight Under 80 lbs
 - User Weight 80 - 120 lbs
 - User Weight 120 - 160 lbs
 - User Weight 160 - 220 lbs
 - User Weight 220 - 260 lbs
 - User Weight 260 - 300 lbs
- Frog Legs PHASE 2 Carbon Fork _____ \$425
Lightweight shock absorbing fork. Available only with 4", 5" or 6" Caster Sizes. Not available with Caster Pin Locks, 350 lb Weight Capacity, 2" Poly caster or 2" Pneumatic w/ Foam Insert caster.
- Caster Pin Locks K0073 _____ \$195
Not available with Performance Caster Housing, 1 1/2" Longer Fork Stem or Frog Legs Fork

Wheel Locks

Wheel Lock

- Push to Lock _____ Std
- Push to Lock Flush Mount _____ N/C
- Pull to Lock _____ N/C
- Low Profile Scissor Lock _____ \$65
Not available with Low Poly, Full Poly or Sentinel Solid tires
- Short Thro Scissor Lock _____ \$110
- Omit Wheel Locks _____ N/C

Wheel Lock Options

- Extension Handles E0961 _____ \$70
Available only with Push to Lock or Pull to Lock Wheel Locks
- Grade Aids E0974 _____ \$150
Available only with Push to Lock Wheel Locks and Pneumatic or Pneumatic with Airless Insert tires

Rear Seat Height and Wheels

Wheel Size / Tire Profile / Rear Seat Height

Select Rear Seat Height in row corresponding to desired Wheel Size and Tire Profile. 14" Seat Depth not available with 24", 25" or 26" wheel sizes.

Wheel Size	Tire Profile	Rear Seat Height					
20"	Low	<input type="checkbox"/> 11.5"	<input type="checkbox"/> 12.5"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16.5"
	Full	<input type="checkbox"/> 12"	<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"
22"	Low	<input type="checkbox"/> 12.5"	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17.5"
	Full	<input type="checkbox"/> 13"	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"
24"	Low	<input type="checkbox"/> 13.5"	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18.5"
	Full	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"
25"	Low	<input type="checkbox"/> 14"	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"
26"	Low	<input type="checkbox"/> 14.5"	<input type="checkbox"/> 15.5"	<input type="checkbox"/> 16.5"	<input type="checkbox"/> 17.5"	<input type="checkbox"/> 18.5"	<input type="checkbox"/> 19.5"
	Full	<input type="checkbox"/> 15"	<input type="checkbox"/> 16"	<input type="checkbox"/> 17"	<input type="checkbox"/> 18"	<input type="checkbox"/> 19"	<input type="checkbox"/> 20"

Rear Seat Height and Wheels continued

Rear Wheel		Rear Wheel Size					
		Price	20"	22"	24"	25"	26"
Rear Wheel Type	Assembled Six-Spoke Mag w/ Treaded Tire and Composite Handrim ^{1,2}	-\$25			<input type="checkbox"/>		
	MAXX Performance Spoke	\$200	<input type="checkbox"/>	<input type="checkbox"/>		<input type="checkbox"/>	<input type="checkbox"/>
	MAXX Performance Spoke	N/C			<input type="checkbox"/>		
	Mag	N/C	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>		
	Superlight Spoke ²	\$295			<input type="checkbox"/>	<input type="checkbox"/> ³	
	Spinergy Spox	\$835		<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

¹ No substitutions available. Do not select Rear Tire Type or Handrim. Measures as Full Profile tire.
² Not available with 350 lb Weight Capacity
³ Available only with Superlight and Plastic Coated Handrims

Rear Tire Type

Full Profile

- Street Pneumatic 1 3/8" E2211, E2212 _____ Std
Not available on 25" wheel size
- Full Poly 1 3/8" E2218 _____ N/C
Not available on Superlight wheel or 25" or 26" wheel sizes
- Pneumatic with Airless Insert 1 3/8" E2211, E2213 _____ \$95
Not available on Superlight wheel or 25" wheel size
- Schwalbe HS-127 1 3/8" E2211, E2212 _____ \$125
Available only on 24" wheel size. Not available on Mag wheel.
- Euro Trek Knobby 1 3/8" E2211, E2212 _____ \$95
Available only on 24" wheel size

Low Profile

- Low Poly 1" E2218 _____ \$50
Not available on Superlight wheel, 25" or 26" wheel sizes or 350 lb Weight Cap.
- Sentinel Solid Tire 1" E2218 _____ \$95
Not available on Superlight wheel or 20" and 22" Mag wheels
- High Pressure 1" E2211, E2212 _____ \$95
Not available on Mag wheel
- Schwalbe Marathon Plus Evolution 1" E2211, E2212 _____ \$135
Not available on Mag wheel, 20" or 22" wheel sizes

Axles

Axle

- Quick Release _____ Std
- Quad Quick Release _____ \$175
Not available with Mag or Spinergy wheels

Camber

Combination of Standard and Adjustable Camber cannot exceed 4°. Selection will default to a max of 4° combined if exceeded.

- Standard Camber _____ N/C
Chair will be set up with this camber selection
- 0° 1° 2° 3° 4°

- Adjustable Camber _____ N/C
Extra camber spacers mounted to the inside of the rear frame for later use
- 1° 2° 3° 4°

Spinergy Options

Select spoke color and handrim tab for Spinergy wheels

Spoke Color:

- Black Red White Blue Yellow

Handrim Tab:

Available only with Aluminum Anodized and Plastic Coated Handrims

- Short Tab _____ \$25
- Long Tab _____ N/C

Handrim

Standard Handrims

- Aluminum Anodized _____ Std
- Superlight Aluminum Anodized _____ \$50
Available only on 24" and 25" wheel sizes. Available only on MAXX Performance or Superlight wheels. Not available with 350 lb Weight Capacity.
- Plastic Coated _____ \$100
- Projection (8 Vertical Projections) _____ \$310
Available only with full profile tires in combination with 22"-26" MAXX Performance Spoke, or 24" Mag wheel
- Omit Handrim _____ N/C

Natural Fit Handrims

Not available on 20" wheel size or 22" Mag wheel. Not recommended with Mag wheel.

- Natural Fit LT Regular Grip _____ \$465
- Natural Fit Standard Regular Grip _____ \$445
- Natural Fit LT Super Grip _____ \$515
- Natural Fit Standard Super Grip _____ \$495
- Natural Fit LT Omit Grip _____ \$395
- Natural Fit Standard Omit Grip _____ \$345

Handrim Option


- 3M Micro Grip _____ \$65
Not available on Plastic Coated or Projection Handrims

Axle Position

If no selection is made position 3 will be used. Position 4 is 7/8" forward of the back cane.

Center of Gravity

- 1
- 2
- 3
- 4



Armrests

Armrest	
Height Adjustable T-Arm E0973	
<input type="checkbox"/> Desk (8.5" - 13.5")	\$275
<input type="checkbox"/> Full (8.5" - 13.5")	\$275
Height Adjustable T-Arm Low E0973	
<input type="checkbox"/> Desk (7.5" - 11.5")	\$275
<input type="checkbox"/> Full (7.5" - 11.5")	\$275
Height Adjustable Flip Back T-Arm E0973	
<input type="checkbox"/> Desk (9.5"-13.5")	\$285
<input type="checkbox"/> Full (9.5"-13.5")	\$285
Flip Back Desk (10.5")	
<i>Not available with Seat Depths less than 16" or Angle Adjustable Back Posts</i>	
<input type="checkbox"/> Desk (10.5")	\$50
<input type="checkbox"/> Full (10.5")	\$50
<input type="checkbox"/> Swing Away (10.5" - 12.5") E0973	\$210
<i>Not available with Angle Adjustable Back Posts or Arm Pads</i>	
<input type="checkbox"/> Swing Away with Side Guard (10.5" - 12.5") E0973	\$285
<i>Not available with Angle Adjustable Back Posts or Arm Pads</i>	
<input type="checkbox"/> Omit Armrests	N/C

Arm Pad

Not available with Swing Away Armrests

<input type="checkbox"/> Waterfall Arm Pad	\$25
<i>Available only with Height Adjustable T-Arm and Height Adjustable T-Arm Low</i>	

Accessories

Removable Side Guard	
Not available with T-Arms, Flip Back or Swing Away with Side Guard Armrests	
<input type="checkbox"/> Composite Side Guards K0108	\$145
Aluminum Straight K0108	
<input type="checkbox"/> Black Matte	\$195
<input type="checkbox"/> Match Frame Color	\$245
<i>Not available with Hydrographics</i>	
Fixed Mount Fender Side Guards K0108	
<i>Not available with 24" wheel size in combination with 14" Rear Seat Height</i>	
<input type="checkbox"/> Aluminum	\$150
<input type="checkbox"/> Carbon Fiber	\$475
Aluminum Adjustable Fender with Foam Pad K0108	
<input type="checkbox"/> Black Matte	\$275
<input type="checkbox"/> Match Frame Color	\$325
<i>Not available with Hydrographics</i>	
<input type="checkbox"/> Carbon Fiber Adjustable Fender K0108	\$485
Anti-Tip	
<input type="checkbox"/> Rear Anti-Tips E0971	\$145
Swing-In Anti-Tip E0971	
<i>Not available with 350 lb Weight Capacity</i>	
<input type="checkbox"/> Single - Left	\$150
<i>Not available with left side O₂ Holder</i>	
<input type="checkbox"/> Single - Right	\$150
<i>Not available with right side O₂ Holder</i>	
<input type="checkbox"/> Pair	\$300
<i>Not available with Cane and Crutch Holder or O₂ Holder</i>	
<input type="checkbox"/> Omit Rear Anti-Tips	N/C
Positioning Belt	
<input type="checkbox"/> 1 1/2" Auto Buckle E0978	\$70
<input type="checkbox"/> 1 1/2" Padded Auto Buckle E0978	\$155
<input type="checkbox"/> 2" Airline Buckle E0978	\$110
Bodypoint Evoflex - Pelvic Stabilizer E0978	\$183
Size: <input type="checkbox"/> Small (16-24") <input type="checkbox"/> Medium (19-29") <input type="checkbox"/> Large (25-37")	
Bodypoint Belt Mounting Kit	\$105
Style: <input type="checkbox"/> Bracket <input type="checkbox"/> Band Clamp 1"	
Other Accessories	
<input type="checkbox"/> Tool Kit	\$45
<input type="checkbox"/> Removable Underseat Pouch	\$50
<input type="checkbox"/> Pediatric Ki Mobility Backpack	\$50
<input type="checkbox"/> Ki Mobility Backpack	\$65
<input type="checkbox"/> Luggage Carrier	\$85
<i>Not available with PRO ELR or Residual Limb Support</i>	
<input type="checkbox"/> Black Neoprene Impact Guards	\$90
<input type="checkbox"/> Clear Spoke Guards K0065	\$115
<i>Not available with Mag wheel</i>	
<input type="checkbox"/> Black Spoke Guards K0065	\$115
<i>Not available with 20" wheel size or Mag wheel</i>	
<input type="checkbox"/> Cane and Crutch Holder	\$195
<input type="checkbox"/> IV Pole	\$265
<i>Available only with Standard Straight Back Posts</i>	
Step Tube	\$40
<i>Not available with Transit Wheels or 350 lb Weight Capacity</i>	
<input type="checkbox"/> Right Side <input type="checkbox"/> Left Side	
O ₂ Holder	
O ₂ Holder E2208	\$295
<input type="checkbox"/> Right Side <input type="checkbox"/> Left Side	
Transit Wheels	
<input type="checkbox"/> Transit Wheels K0108	\$145
<i>Not available with Step Tube, Rear Seat Heights less than 14" or 350 lb Weight Capacity</i>	

Color

Gloss Finish	
<input type="checkbox"/> Black	N/C
<input type="checkbox"/> Blue Skies	N/C
<input type="checkbox"/> Bubblegum Pink	N/C
<input type="checkbox"/> Candy Blue	N/C
<input type="checkbox"/> Candy Purple	N/C
<input type="checkbox"/> Candy Red	N/C
<input type="checkbox"/> Copper	N/C
<input type="checkbox"/> Granny Smith	N/C
<input type="checkbox"/> Metallic Orange	N/C
<input type="checkbox"/> Yellow	N/C

Matte Finish	
<input type="checkbox"/> Charcoal Gray	N/C
<input type="checkbox"/> Emerald Green	N/C
<input type="checkbox"/> Iridescent Burgundy	N/C
<input type="checkbox"/> Iridescent Navy	N/C
<input type="checkbox"/> Military Green	N/C
<input type="checkbox"/> Retro Red	N/C
<input type="checkbox"/> Royal Blue	N/C
<input type="checkbox"/> Shadow Black	N/C
<input type="checkbox"/> Smooth White	N/C

Hydrographics	
Hydrographics	\$695
<input type="checkbox"/> AOR1 (Desert Digital Camo)	
<input type="checkbox"/> Pink Camo	
<input type="checkbox"/> Muddy Girl	
<input type="checkbox"/> NEXT Camo	
<input type="checkbox"/> Carbon Fiber ¹	
<input type="checkbox"/> Zebra ¹	
¹ Select Hydrographics Base Color. Zebra not available in Black.	
Hydrographics Base Color	N/C
<input type="checkbox"/> Black	<input type="checkbox"/> Blue
<input type="checkbox"/> Granny Smith	<input type="checkbox"/> Pink
<input type="checkbox"/> Purple	<input type="checkbox"/> Red

Colored Anodizing Package				
<i>Select desired color. If no selection is made, black will be used.</i>				
	Option	Standard Caster Housing ²	Caster Fork ³	Axle Plate
	Retail	\$75	\$75	\$75
Color	Black¹	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Blue	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Bright Green	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Gold	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Orange	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Pewter	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Pink	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Purple	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>
	Red	<input type="checkbox"/>	<input type="checkbox"/>	<input type="checkbox"/>

¹ Black is no charge
² Available only if corresponding option selection was made
³ Not available with 2" wide casters, Frog Legs Ultra-Sport or PHASE 2 Carbon Forks

Catalyst Front Seat Height Matrix																											
4" Caster				5" Caster				6" Caster				7" Caster				8" Caster											
FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork	FSH	Super Low	Stem	Fork								
15 ¹	13.5 ¹	Std	4"	16	14.5	Std	4"	17	15.5	Std	5"	18	17	Std	6"	19	18	Std	6"								
15.5	14			16.5	15.5	3/4		18	16.5	3/4		18.5	17.5			20	18.5	3/4									
16 ²	14.5 ²	17.5		16	1 1/2	18.5		17	1 1/2	19		17.5	20.5	n/a		1 1/2											
16.5	15	3/4		5"	16	15	Std	5"	17	16	Std	6"	19.5	18		3/4	7"	19	17.5	Std	7"						
16.5 ²	15 ²	1 1/2			16.5	15.5			17.5	16.5			18.5	18.5				1 1/2	19.5			18					
17	15.5	3/4			17	15.5	18		17	18.5			17.5	19		19		18.5									
15.5	14.5	Std			5"	17.5	16		3/4	5"	18		16.5	1 1/2		6"		18.5	17.5	Std	7"	20.5	19	Std	7"		
16	15					17.5	16.5				1 1/2		18.5					17	19			18	19.5			18.5	
16.5	15	3/4				18	17				19		17.5					19.5	18.5	20		19					
17	15.5	1 1/2				5"	16.5		15.5		Std		5"	18.5				17.5	1 1/2	6"		20	19	3/4	7"	20.5	n/a
17	16		17				16		19					18	19			18									
17.5	16.5	17.5	16.5				19.5		18.5					20	19			20.5				19	21	n/a			
17	15.5	Std	6"				17.5		16		3/4			6"	18			16.5	Std			7"	20	19		3/4	7"
17.5	16			18			16.5	18.5	17			19			17.5		20.5	n/a									
17.5	16.5	3/4		18.5			17	19.5	18			20			18.5		20.5	19									
18	17	1 1/2		6"			18	17	1 1/2		6"	19.5			18		1 1/2	7"	20.5				19	1 1/2		7"	
18.5	17				18.5		17.5	18.5		17.5		19			18	20			19								
19	17.5	19			18		19	18		19		18			20	19											
18	17	Std			7"		18	16.5	Std	7"		19.5			18.5	3/4	7"		20		19		1 1/2	7"			
18.5	17.5					18.5	17	19.5				18	20		19												
19	17.5	3/4				19	17.5	20				18.5	20.5		19				21	n/a							
19.5	18	1 1/2				7"	18.5	17.5	3/4			7"	20.5		19	1 1/2			7"	21	n/a		1 1/2		7"		
19.5	18.5		19				18	20.5					19	21	n/a												
20	19	19.5	18				20	18.5					20.5	19	21					n/a							
							19.5	18	1 1/2																		
				20			18.5																				
				20.5			19																				

¹ Front Seat Height not available with 4 x 1" Poly Caster
² Not available with 4 x 1" Poly, 4 x 1.5" Poly or 4 x 1.5" Soft Roll Aluminum

Seat heights with "FSH" heading available with Standard Front Frame. Seat heights with "Super Low" heading available with Super Low Front Frame. Actual seat height may vary depending upon configuration.

Axiom Back Accessories continued

Therafin Posture Support

Selection of a modification from the Axiom Back Modifications section required in order to mount posture supports. Not available with Low Backs.

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Small Static Therafit 2-Piece Chest Strap with H & L Adj and 1" Safety Strap	60153T	E0960	\$144
<input type="checkbox"/>	Medium Static Therafit 2-Piece Chest Strap with H & L Adj and 1" Safety Strap	60154T	E0960	\$144
<input type="checkbox"/>	Large Static Therafit 2-Piece Chest Strap with H & L Adj and 1" Safety Strap	60155T	E0960	\$144
<input type="checkbox"/>	Small Dynamic Therafit Chest Strap with Dual D-Ring Adjustment	60110T	E0960	\$113
<input type="checkbox"/>	Medium Dynamic Therafit Chest Strap with Dual D-Ring Adjustment	60111T	E0960	\$113
<input type="checkbox"/>	Large Dynamic Therafit Chest Strap with Dual D-Ring Adjustment	60112T	E0960	\$113
<input type="checkbox"/>	Small Dynamic Therafit Comfortfit Full Vest	30589T	E0960	\$167
<input type="checkbox"/>	Medium Dynamic Therafit Comfortfit Full Vest	30590T	E0960	\$167
<input type="checkbox"/>	Large Dynamic Therafit Comfortfit Full Vest	30591T	E0960	\$167
<input type="checkbox"/>	XLarge Dynamic Therafit Comfortfit Full Vest	30592T	E0960	\$167
<input type="checkbox"/>	XXLarge Dynamic Therafit Comfortfit Full Vest	30593T	E0960	\$167
<input type="checkbox"/>	Small Dynamic Therafit Comfortfit Trim Vest	30618T	E0960	\$169
<input type="checkbox"/>	Medium Dynamic Therafit Comfortfit Trim Vest	30619T	E0960	\$169
<input type="checkbox"/>	Large Dynamic Therafit Comfortfit Trim Vest	30620T	E0960	\$169
<input type="checkbox"/>	XLarge Dynamic Therafit Comfortfit Trim Vest	30621T	E0960	\$169
<input type="checkbox"/>	XXLarge Dynamic Therafit Comfortfit Trim Vest	30622T	E0960	\$169

Axiom Back Modifications

Extended lead times will apply for all modifications

Axiom Back Modifications

Choose modification to add holes to the back shell in order to mount posture supports. Not available with Low Backs.

	Description	Item Code	HCPC	Retail
<input type="checkbox"/>	Modification for 4 Point Harness	AXM4PT	K0108	N/C
<input type="checkbox"/>	Modification for 2 Point Chest Strap	AXM2PT	K0108	N/C
<input type="checkbox"/>	Modification for Shoulder Guides <i>Includes shoulder guides</i>	AXMSG	K0108	\$50

Painted Back Shell Modification

If no selection is made, standard Shadow Black will be used with no additional lead time

	Gloss Finish	Item Code	Retail
<input type="checkbox"/>	Black	AXMBLK	N/C
<input type="checkbox"/>	Blue Skies	AXMSKY	N/C
<input type="checkbox"/>	Bubblegum Pink	AXMPNK	N/C
<input type="checkbox"/>	Candy Blue	AXMCBL	N/C
<input type="checkbox"/>	Candy Purple	AXMPUR	N/C
<input type="checkbox"/>	Candy Red	AXMRED	N/C
<input type="checkbox"/>	Copper	AXMCOP	N/C
<input type="checkbox"/>	Granny Smith	AXMGSM	N/C
<input type="checkbox"/>	Metallic Orange	AXMMOR	N/C
<input type="checkbox"/>	Silver	AXMSIL	N/C
<input type="checkbox"/>	Teal	AXMTEA	N/C
<input type="checkbox"/>	Yellow	AXMYLW	N/C

	Matte Finish	Item Code	Retail
<input type="checkbox"/>	Charcoal Gray	AXMCGM	N/C
<input type="checkbox"/>	Emerald Green	AXMEGM	N/C
<input type="checkbox"/>	Iridescent Burgundy	AXMIBM	N/C
<input type="checkbox"/>	Iridescent Navy	AXMINM	N/C
<input type="checkbox"/>	Military Green	AXMAGM	N/C
<input type="checkbox"/>	Retro Red	AXMRRM	N/C
<input type="checkbox"/>	Royal Blue	AXMRBM	N/C
<input type="checkbox"/>	Smooth White	AXMSWM	N/C