

Introducing the REDESIGNED

Effective January 1, 2024

2024 Evaluation Form



Chill-Out Chairs alternative seating system can help improve the user's positioning, providing the benefit of reduced transfers throughout the day and a comfortable seating experience. It features the redesigned, more defined reclined angle that aids in the reduction and easing of kyphosis of the occupant.

Four-Step Order Process:

1. Select style and size.
2. Select foam package and fabric.
3. Choose accessories and foam warranty.
4. Call our team if you have any questions.

1-800-661-9915 or visit chill-outchair.com

Rock'er

Known for its adjustable tilt and calming rocking motion, the Rock'er can improve the alignment of users of all ages and abilities, while redistributing pressure.



Roll'er




With a rolling wheelbase, the Roll'er can reduce the number of transfers required throughout the day – perfect for home, school or therapy settings.

The Chill-Out Chair comes standard with a multi-purpose ottoman and multi-use modular support wedge, covered in matching fabric.

Sizing and Dimensions:

Recommended User Size Ranges / Chair Dimensions

Chill-Out Chair Size	Height Range (inches)	Height Range (cm)	Weight Range (lbs)	Weight Range (kg)	Overall Length (A)	Overall Height (B)	Overall Width (C)	Seat Width (#1)	Seat Depth (#2)	Seat Height (#3)	Rock'er Weight (lbs) ***	Roll'er Weight (lbs) ***
Small	40" - 52"	102 - 132	30 - 80	14 - 36	29"	33"	19"	12"	12"	26"	21	30
Medium	50" - 62"	127 - 157	55 - 120	25 - 54	34"	36"	21"	14"	14"	31"	26	37
Narrow Large	60" - 68"	152 - 173	70 - 155	32 - 70	40"	41"	22"	14"	16"	38"	33	42
Large	60" - 68"	152 - 173	90 - 175	41 - 79	40"	41"	24"	16"	16"	38"	36	46
Extra-Wide Large	60" - 68"	152 - 173	175 - 300	79 - 136	40"	46"	31"	21"	17"	38"	46	56



Select Your Model (check selection)

	SMALL	MEDIUM	NARROW LARGE	LARGE	EXTRA WIDE LARGE
Rock'er	\$2,150 (COCROCKRSM)	\$2,500 (COCROCKRMD)	\$2,650 (COCROCKRNLG)	\$2,750 (COCROCKRLG)	\$3,050 (COCROCKRXLGW)
Roll'er	\$2,450 (COCROLLRSM)	\$2,800 (COCROLLRMD)	\$2,950 (COCROLLRNLG)	\$3,050 (COCROLLRLG)	\$3,400 (COCROLLRXLGW)


Select Your Foam Package (check selection)

	SMALL	MEDIUM	NARROW LARGE	LARGE	EXTRA WIDE LARGE
Essential (Included with base price)	No Charge (3640121)	No Charge (3640321)	No Charge (3640521)	No Charge (3640421)	Not Available
TEC Foam Package	\$575 (3640125)	\$585 (3640325)	\$660 (3640525)	\$660 (3640425)	\$680 (3640625)

Essential Foam Package	TEC Foam Package (Therapists'/Educators Choice)
Essential Foam Package is made of one density to offer great comfort and support for those 100lbs or less. Included in base price of Chill-Out Chair	Provides soft comfort, firm support and best durability for all user weights. It features softer foam under the knees and beside the headrest to increase comfort and eliminate pressure points. Recommended for heavier users.
Standard One-Year Warranty	Standard Three-Year Warranty

Select Your Fabric Options (check selection) Spirit, Prairie and Seaside fabrics are included with the base price of your chair.

Spirit Series: Buttery soft, easy to wipe vinyl fabrics. Endurepel™ Armour topcoat protection provides an easy clean finish that has anti-bacterial and anti-mildew properties.




ENDUREPEL™
ARMOUR

Prairie Series: Woven fabrics with a soft, dry hand and lots of visual, stylish texture. Treated with Endurepel™ Shield offering anti-stain and moisture repellent technology, combined with easy clean and odor protection properties.



Seaside Series: Soft and comfortable distressed suede like fabrics. Treated with Endurepel™ Shield offering anti-stain and moisture repellent technology, combined with easy clean and odor protection properties.




Premium Fabrics: (check size and fabric)

Escape Series: Patterned and solid, intricately textured fabrics. Durable, anti-microbial and built for indoor and outdoor use. SolarShade fabrics are created with high UV rating, water repellency and anti-microbial technology.



Small	(9640221) (\$170)
Medium	(9640321) (\$215)
Narrow Large	(9640521) (\$260)
Large	(9640421) (\$260)
Extra Wide Large	(9640621) (\$285)



Accessories (check selection)

	SMALL	MEDIUM	NARROW LARGE	LARGE	EXTRA WIDE LARGE
Occiput Headrest	\$340 (2640106)	\$340 (2640306)	\$340 (2640306)	\$340 (2640406)	\$340 (2640406)
Activity Wedge	\$290 (2644031)	\$290 (2644032)	\$290 (2644034)	\$320 (2644033)	\$320 (2644035)
Memory Foam Liner	\$145 (3640190)	\$165 (3640390)	\$190 (3640590)	\$190 (3640490)	\$210 (3640690)
Medbloc Liner	\$290 (2640110)	\$290 (2640310)	\$340 (2640510)	\$340 (2640410)	\$340 (2640610)
3-Piece Foam Positioning Aids	\$270 (2040750)	\$270 (2040751)	\$270 (2040751)	\$270 (2040752)	\$270 (2040753)
5-Point Chest Harness	\$640 (2141071)	\$640 (2141072)	\$640 (2141073)	\$640 (2141073)	\$640 (2141074)


New!



Occiput Headrest

Provides targeted support for the neck and head. Can be adjusted to allow for growth.

Memory Foam Liner
Adds an extra layer of foam for comfort.



Activity Wedge
Encourages upright positioning.

Medbloc Liner
Heat-sealed liner that prevents moisture seepage and protects the hygienic integrity of the foam.



3-Piece Foam Positioning Aids
Increase support and positioning with our foam positioning aids.

Includes the new Multi-Use Modular Support Wedge!



New!

3-Point Pelvic Positioning Belt

Our 3-Point Padded Pelvic Positioning Belt is great for pelvic stabilization and safety, keeping the hips back in the chair.




\$480
(2640041)

Large Accessories Sleeve Pocket

Keep all positioning aids at the ready when not in use with the Accessories Sleeve Pocket.

Ideal for storing positioning pads.

\$150
(2640090)



5-Point Chest Harness

Increase support and positioning while holding the user in the chair securely.

Removable Abduction Pommel

The Removable Abduction Pommel sits between the user's legs to help encourage hip abduction. The pommel is attached to the Chill-Out Chair with strong velcro straps and can be removed as desired.



\$165
(2640048)

Ottoman Support Kit

Our ottoman support kit with multi height ottoman keeps the ottoman off the ground and allows you to roll your Chill-Out Chair with the ottoman in place.

***Not available in XW-LRG**



\$200
(2040090)

Pushbar

Our Pushbar makes moving the Roller Chair from room-to-room effortless, resulting in less transfers in-and-out of the chair.



\$265
(2040055)

Small Accessories Sleeve Pocket

Store your portable electronics, adaptive communication devices, and video players right with your chair.

Always have them at hand!



\$75
(2040053)

Feeding/Activity Tray

Perfect for feeding, play or communication, the Feeding and Activity Tray is the ideal accessory to assist with individual development.

XW-LG measures 35" across



SM, MED, NLG, LG **\$715** (2040051)
XWLG **\$1,455** (2040052)

Rocking Base for Roll'er Chair

Turn your Roll'er model into a Rock'er! Roll the Roll'er model onto the track for back-and-forth tilting motion or to add fixed tilt to the chair.

Please note, the Rocking Base requires our Ottoman Support Kit. **Some assembly required.**

***Not available in Small or Extra Wide Large.**



Medium **\$700** (2040041)
Narrow Large **\$700** (2040042)
Large **\$700** (2040043)

Transport Chair Cover

Bring your Chill-Out Chair wherever you go with the Transport Chair Cover that fits over both the chair and ottoman when placed on the seat. Prevents damage and protects integrity of chair when transporting.

New!



Small-Large **\$155** (6340700)
XWide Large **\$200** (6340702)

Foam Warranty

WARRANTY PACKAGE	ESSENTIAL FOAM	TEC FOAM
ADDITIONAL One Year Warranty	\$125 (WARR_1YR)	\$125 (WARR_1YR)
ADDITIONAL Three Year Warranty	\$300 (WARR_3YR)	\$300 (WARR_3YR)
LIFETIME Warranty	N/A	\$450 (WARR_LIFE)

Customizations



Terms & Conditions

Prices and fabric colors are subject to change without notice. Prepaid or net 30 days O.A.C. All taxes extra. Freight extra. Manufacturer reserves the right to change prices, colors, features, dates, locations and products without notice. Residential delivery charge may apply. Duty and brokerage included where applicable.

CUSTOMER SIGNATURE

Photo Release

I hereby give permission to Freedom Concepts Inc. to use photographs of me or my child for the sole purpose of the company's promotional use, and will not be distributed to third parties. YES NO

I would like to receive regular updates. MAIL EMAIL NO

SIGNATURE OF OCCUPANT/PARENT

User Information

User Name: _____

Date of Birth: MONTH _____ /DAY _____ /YEAR _____

Height: _____ INCHES _____ CM Weight: _____ LBS _____ KG

Shoulder Width: _____ INCHES _____ CM Hip Width: _____ INCHES _____ CM

Diagnosis: _____

Purchaser Name: _____

Phone: _____

Email: _____

Address: _____

Prov/State: _____ City: _____

Postal/Zip Code: _____

Order Information

Order Date: _____ Quote Date: _____

Chill-Out Representative: _____

Payment Terms: _____

CREDIT CARD (CHECK ONE) Visa MasterCard Exp. Date: _____

Card Number: _____

Name On Card: _____

Shipping Address: _____

Prov/State: _____ City: _____

Postal/Zip Code: _____ Federal ID/SS#: _____

Privacy Policy

Freedom Concepts Inc. collects, uses and discloses personal information of the Customer and/or Rider to properly assemble and manufacture the mobility device and may contact the Rider's therapist for additional personal information to assist with this process. The collection, use and disclosure of the personal information of the Customer and/or the Rider is governed by any consent which may, from time to time, be provided, and by the Freedom Concepts Inc. Privacy Policy which is available on our website www.freedomconcepts.com or by contacting our office. By signing below, you hereby warrant that you have the authority to bind the Rider and hereby authorize the collection, use, and disclosure of personal information as herein provided and in accordance with the Freedom Concepts Inc. Privacy Policy.

	Total - Page 1	
	Total - Page 2	
	Total - Page 3	
	Sub-Total	
	Additional Freight	
	Tax 1	
	Tax 2	
	Other	
	Total	