



MK Battery

MK Battery
1631 South Sinclair Street
Anaheim, CA 92806
www.mkbattery.com
800-372-9253

WHEELCHAIR APPLICATION CHART

| Manufacturer | Chair Model | Battery Type | Qty. | Recommended Charger |
|--|---|----------------------|--------|-------------------------------|
| Access Point Medical | AXS5018, AXS5020, AXS-6000, AXS-7000, AXS33F | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | AXS8000 | 12260 | 2 | LS24/4 |
| | AXS-31, AXS42P, AXS32P | 12-12 | 2 | LS24/4 |
| Amigo Mobility | FD (Front Drive), Value Shopper, Value Shopper XL (2 each) | U1 | 1 | LS24/4, M-D-EN0801, M-L-12090 |
| | RD (Smart Shopper, Excite, Power Shopper, Escort), EXT350, HD450 | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | RT Express, MC, MCX, RT Express Jr. | 17-12 | 2 | LS24/4 |
| | Escort | 12260 | 2 | LS24/4 |
| | GT Transport, B-6 | 40-12 | 2 | M-D-EN0801, M-L-12090 |
| | Classic FD, TravelMate RD TravelMate | 12-12 7-12 | 2 2 | LS24/4 LS24/4 |
| Aquatec | Fortuna Bath Lift | 3-12 | 1 | LS12/0.6 |
| Braun | T1100/ T1200F (series 1-3) | 24 | 2 | M-D-EN0801, M-L-12090 |
| | T1100/ T1200F (series 4 & after) | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| Bruno | *PWC-2200 RWD, PWC-2300 FWD, Cub RWD (*Optional Cub 46 only) | U1 (*22NF optional) | 2 | M-D-EN0801, M-L-12090 |
| | Stairlifts (Electra-Ride - all models) | 7-12 | 2 | LS24/4 |
| | Typhoon C3 RWD | 17-12 | 2 | LS24/4 |
| | Cub 35 FWD (U1 optional) | 17-12 or U1 | 2 | LS24/4 |
| Chauffer Mobility | Chauffer Series, Viva Power heavy duty mini models | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Lil Taxi | 17-12 | 2 | LS24/4 |
| | Viva Power 645 | 22NF or 24 | 2 | M-D-EN0801, M-L-12090 |
| Dalton Medical | ePower aluminum, ePower Pediatric | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Primechair, MW, RW | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | ePower Bariatric, Tacahe Rehab PC1450RHB | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Tacahe PC1350, Tacahe Hvy Dty PC1450, Tacahe PC1450, Tacahe PC2450 Tilt | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Tacahe Hvy Dty PC-P710, P720 | 24 | 2 | M-D-EN0801, M-L-12090 |
| | 3 wheel SC-135 HD, 4 wheel SC-S145 HD | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | 3 wheel SC-43 | 17-12 | 2 | LS24/4 |
| | 3 wheel SC-235, 4 wheel SC-245 | 22NF | 2 | M-D-EN0801, M-L-12090 |
| Damaco | Ovation, Electro-Lite, Elite | 12260 or U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | D90, D99 | 12260 or U1 or 17-12 | 2 | LS24/4 |
| | Child's Elite (14 x14) | 17-12 | 2 | LS24/4 |
| Drive Medical Design | Cirrus - (DP116, DP118, DP120), Cirrus Plus EC | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Image, Cirrus Plus, Odyssey, Sunfire & Sunfire Plus, Gladiator GT, Daytona | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Hawk & Hawk Compact, Phantom & Mini Phantom, Geo, Phoenix | 12-12 | 2 | LS24/4 |
| | Sunfire Gladiator | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Denali | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Gladiator | 24 | 2 | M-D-EN0801, M-L-12090 |
| Electric Mobility | Turnabout/Stowaway, Mini base scooter 130 | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Rascal Fold & Go 760, Rascal 115 | 12-12 | 2 | LS24/4 |
| | Rascal AutoGo 550, Rascal AutoGo 555, Rascal Fold & Go 370, Rascal 320 | 17-12 | 2 | LS24/4 |
| | Rascal 500T, Rascal 140T | 12260 | 2 | LS24/4 |
| | Rascal - 300, 312, 314, 318, 326, 330 Tilt, 610T Scooter, 610F Scooter | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Candy Apple | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Rascal 655 | 40-12 | 2 | M-D-EN0801, M-L-12090 |
| | Rover, Sparky, Squire, Rascal 310, Rascal 710 | 24 | 2 | M-D-EN0801, M-L-12090 |
| Etac Balder | Finesse, Junior, Liberty | 22NF | 2 | M-D-EN0801, M-L-12090 |
| Everest & Jennings | Metro Power, Metro Power III, Navigator, Quest, Sprint II, 3N, Tempest, Kid Power | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Hot Wheels, Carrette, Mobie, Everest | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Sabre LTD, Marathon, Marathon LE, 3P, 3W, 3H, 3V, Magnum | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Xcaliber, Sprint, Sprint Plus, MX, Sabre ES, Sabre ES-GT (22NF) | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Xcaliber Power Recliner, MX Power Recliner, Magnum Power Recliner | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Lancer 2000, Lancer Base, Solaire, Solaire Base, Lancer | 24 | 2 | M-D-EN0801, M-L-12090 |
| | Everest Lite | 12-12 | 2 | LS24/4 |
| | Explorer, Vortex | 24 | 2 | M-D-EN0801, M-L-12090 |
| Evermed | ECS, EBS | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | ECW, EBW | 24 | 2 | M-D-EN0801, M-L-12090 |
| Falcon Rehab | Hi-Rider | 27 | 2 | M-D-EN0801, M-L-12090 |
| Frank Mobility | E-Fix (E19 Model) | 7.2-12 | 2 | LS24/4 |
| | E-Fix | 12-12 | 2 | LS24/4 |
| | E-Fix (with high capacity case) | 17-12 | 2 | LS24/4 |
| | E-Fix Basic | 12-12 or 17-12 | 2 | LS24/4 |
| | E-Fix Plus | 22NF | 2 | M-D-EN0801, M-L-12090 |
| Golden Technology | Golden Companion, Golden Companion II, Alero, Companion | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Buzzaround, Buzzaround Lite | 12-12 | 2 | LS24/4 |
| | Alante, Golden Compass | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Alante Jr. | 12260 | 2 | LS24/4 |
| | Alante II, Avenger | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Liteway | 17-12 | 2 | LS24/4 |
| Handicare Wheelchairs (formerly Movingpeople.net) | Raptor, Raptor PS | 34 | 2 | M-D-EN0801, M-L-12090 |
| | Eligo 2 | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Ibis XP | U-1 or 40-12 | 2 | M-D-EN0801, M-L-12090 |
| | Beatle | 40-12 | 2 | M-D-EN0801, M-L-12090 |
| | Alex | 40-12 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Atlantic, Pacific | 24 | 2 | M-D-EN0801, M-L-12090 |
| Hoveround | MPV1-MPV5, HRV 100, Activa LX, Activa DM | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Activa GLX, Activa Forerunner | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Activa Mini, Activa | 12-12 | 2 | LS24/4 |
| | Teknique RWD,FWD,GT | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Teknique HD, XHD | 24 | 2 | M-D-EN0801, M-L-12090 |

| Manufacturer | Chair Model | Battery Type | Qty. | Recommended Charger |
|---|--|---------------------|-------------|-------------------------------|
| Invacare | Pronto R2, Ranger II FWD, MWD | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Xterra GT, Power Tiger | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Power 9000 (16" or wider), Ranger II RWD, R51LXP | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | New Nutron Series: R32LX, R50LX, R51LX, R51, Excel | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | New Nutron Series: R32, Pronto M6, Pronto M71, P9000XD1 | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | 3G Storm Arrow, 3G Storm Ranger X, Storm Arrow FWD, 3G Ranger X RWD | 24 | 2 | M-D-EN0801, M-L-12090 |
| | TDX-4, TDX-5 | 24 | 2 | M-D-EN0801, M-L-12090 |
| | P7E, LX-3, LX-3 plus & LX-4 Scooters, Ranger II, FWD Jr. | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Pronto M50, M51 & M51-CG, Pronto R2 250-series, Pronto M41, Excel 250-series | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Ranger II 250-S FWD/RWD/MWD, TDX Spree, TDX SC | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Zoom-3, Zoom 220 | 12-12 | 2 | LS24/4 |
| | Lynx SX-3, ATM Take Along | 17-12 | 2 | LS24/4 |
| | 3G Storm Torque SP, 3G Storm Torque 3, TDX SR, TDX-3 | 22NF or 24 | 2 | M-D-EN0801, M-L-12090 |
| Karma Medical Products | KS-636, KS-646, KS747, KP25, KS-343, KS-333 | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | KS-848, KS-838 | 24 | 2 | M-D-EN0801, M-L-12090 |
| | KS-242 | 12260 | 2 | LS24/4 |
| | KP10, KP10.2 | 17-12 | 2 | LS24/4 |
| | KS-747.2, KS-737.2, KP45.3 T/TR, KP60 | 22NF | 2 | M-D-EN0801, M-L-12090 |
| Leisure Lift/Pace Saver/ Burke Mobility | KP40, KP45.3 | 40-12 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Scout Midi-Drive: (RF, NP) | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| Levo | Scout Midi-Drive (RF4) | 22NF or 24 | 2 | M-D-EN0801, M-L-12090 |
| | Scout M Series: (M1, M1 PBR, M2) | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Passport (All), Plus III (All), Junior Premier | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Espre (All), Eclipse (All) | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | LCM 25amp, Comfort II | 12260 | 2 | LS24/4 |
| Lifestand | LCM 50amp | 12260 | 4 | LS24/4 |
| | LCM 36amp | 40-12 | 2 | M-D-EN0801, M-L-12090 |
| | Combi, Combi Jr. | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Levo LC | 7-12 | 1 | LS6/1.2 |
| | LCE, LCE-V, LCE-KID/JR. | 2.9-12 | 2 | N/A |
| Medline Industries | LSC Compact Standing Wheelchair | 17-12 | 4 | LS24/4 |
| | Gemini 1 | 17-12 | 2 | LS24/4 |
| | Strider Micro, Strider Mini | 12-12 | 2 | LS24/4 |
| | Strider Maxi | 24 | 2 | M-D-EN0801, M-L-12090 |
| Merits Health Products | Gemini 2, Zenith, Strider Midi | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | MP11N, MP11A, MP11W (Travel Ease), MP11N-FR, MP11A-FR | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | MP11X, MP11U, MP3R (Gemini), MP3W | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Travel-Ease Commuter (P181, P182, P183, 184) | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Travel-Ease Regal P301, P3011 | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | MP3CF, MP3F, MP3C, MP3U | U1 or 40-12 | 2 | M-D-EN0801, M-L-12090 |
| | Pioneer 3 (S131, S1312, S1313), Pioneer 4 (S141, S1412, S1413) | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Pioneer (S132, S133, S135, S142, S143, S145) | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Travel-Ease (P101, P107, P171, P200) | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Travel-Ease Regal (P31862, P310, P318, P328), Mid-Wheel Drive P32E | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Pioneer 9 (S331), Pioneer 10 (S341) | 24 | 2 | M-D-EN0801, M-L-12090 |
| | Travel-Ease Regal (P313, P314, P31312, P31461, P31411, P31362) | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Big Boy Regal P710 & 720 | 24 | 2 | M-D-EN0801, M-L-12090 |
| | Pioneer 5 (S534), Pioneer (S541), Travel-Ease Regal (P321) | 12-12 | 2 | LS24/4 |
| | Pioneer (S538, S548) | 17-12 | 2 | LS24/4 |
| Meyra | Pioneer 1 (S235), Pioneer 2 (S245), Pioneer (S246), Travel-Ease Regal (P120, P320) | **22-12 | 2 | LS24/4 |
| | Pioneer 1 (SP232, SP23), Pioneer 2 (SP242, SP24) | 17-12 | 2 | LS24/4 |
| | Pioneer 3 (S132, SP43), Pioneer 4 (S142, SP44) | U1 or 40-12 | 2 | M-D-EN0801, M-L-12090 |
| | 967, 2472, 2482-Narrow | 40-12 | 2 | M-D-EN0801, M-L-12090 |
| Movingpeople.net (Fortress Scooters) | 1700FS, 1704FS | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | 2000FS Three & Four Wheel | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Winner 3 & Four Wheel, Champior | 22NF or 24 | 2 | M-D-EN0801, M-L-12090 |
| Orthofab/ Lifestyles (Fortress Wheelchairs) | 655,760,760V,760N, Spirit, Victory, 720,770, Kameleor | 24 | 2 | M-D-EN0801, M-L-12090 |
| | 755FS,1000FS, Commuter, Star 3 & Star 4 Scooter | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Titan | 40-12 | 2 | M-D-EN0801, M-L-12090 |
| | Lifestyle Lift | 2,3-12 | 2 | N/A |
| Ortho-Kinetics | Sierra, Colt, Encore, Triumph, MVP, Explorer, Lark 3 & 4 | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Lark XT, Bravo, Bravo Plus, Pony | U1 | 1 or 2 | LS24/4, M-D-EN0801, M-L-12090 |
| Otto Bock | A200, Skippi | U-1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | B600, C1000 | 24FT | 2 | M-D-EN0801, M-L-12090 |
| Permobil | Chairman, Max 90, Super 90, Chairman 2K, Chairman Entra, Chairman HD: | 24 FT | 2 | M-D-EN0801, M-L-12090 |
| | Chairman Basic, Playman/Robc | 24 FT | 2 | M-D-EN0801, M-L-12090 |
| | C300: Corpus & Corpus Jr. , 300PS, 300 PS Jr. | 34 | 2 | M-D-EN0801, M-L-12090 |
| | C350: Corpus & Corpus Jr. | 34 | 2 | M-D-EN0801, M-L-12090 |
| | C400: Corpus, Corpus Jr., Lowrider, 400PS, Stander Jr. | 24 FT | 2 | M-D-EN0801, M-L-12090 |
| | C500: Stander, Corpus, Lowrider, Street Corpus, 500PS | 24 FT | 2 | M-D-EN0801, M-L-12090 |
| | K300RS Jr. | 34 | 2 | M-D-EN0801, M-L-12090 |
| | Hexior, Koala | 40-12 | 2 | M-D-EN0801, M-L-12090 |
| | Trax | 27 or 31 | 2 | M-D-EN0801, M-L-12090 |
| Pillar Technology | Express, Express LX, Blazer Plus | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Special Edition, 4 Wheeler 409 /410, Powerchai | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Deluxe Powerchair | 22NF | 2 | M-D-EN0801, M-L-12090 |
| Pride * Chairs w/ power seats require 22NF | Jazzy (610, 1107, 1103, 1113, 1113 ATS, 1143), Dynamo, Jet 3, Jet 7, Pride LX | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Jazzy (600, 1105, 1115, 1121), Jet 2, Jet 2 HD, Jet 12, PHC 5, Jet 1 w/out Active Trac | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Jazzy (1100, 1104, 1120-2000, 1170XL, 1400, 1420, 1470) | *24 | 2 | M-D-EN0801, M-L-12090 |
| | Blast (650, 850, HD), PHC 1, PHC 10, Jet 1 w/ Active Trac | 24 | 2 | M-D-EN0801, M-L-12090 |
| | Quantum 6000 & 6000XL, 6000Z (group 4), R400C | 24 | 2 | M-D-EN0801, M-L-12090 |
| | Shuttle, Sidekick, Sundancer, Victory, Maxima, Dynamo (discont.), Legend | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Legend, Celebrity, Celebrity XL, Cyclone, Lase | U1 or 40-12 | 2 | M-D-EN0801, M-L-12090 |
| | Legend XL | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Sonic, Go Chair | 12-12 | 2 | LS24/4 |
| | GoGo Ultra X, GoGo SC40E & SC44E | 12-12 | 2 | LS24/4 |
| | GoGo SC40LR & SC44LR, GoGo Elite Traveller Plus | 17-12 | 2 | LS24/4 |
| | Revo, Jazzy Select Traveller, Select, Select 6 & 7, Quantum 610, Select 6 Ultra | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Jazzy 600XL, Jazzy 600, Select 14XL, Select 14, Quantum 600, 6000Z (group 3) | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Rally | 17-12 or U1 | 2 | M-D-EN0801, M-L-12090 |
| | Hurricane, Celebrity XL | 40-12 | 2 | M-D-EN0801, M-L-12090 |

| Manufacturer | Chair Model | Battery Type | Qty. | Recommended Charger |
|------------------------------|--|---------------------|-------------|-------------------------------|
| Pride (cont.) | Hurricane (2005-present) | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Wrangler | 27 | 2 | M-D-EN0801, M-L-12090 |
| | Victory XL, Boxster (discont.) | U1 or 40-12 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| Quickie | P222, P300(All), P200, P210, P320 | 24 | 2 | M-D-EN0801, M-L-12090 |
| | S626, G424, Aspire F10/F11, Freestyle F11, Rhythm, Groove, P22C | 22NF or 24 | 2 | M-D-EN0801, M-L-12090 |
| | S646, S646SE, P222SE | 22NF or 24 | 2 | M-D-EN0801, M-L-12090 |
| | V100, P110, P190, Z500, V121, S525, V521 | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | P110 (14" wide), V121 (14" wide), V521 (14" wide), S-11, Rhapsody, Prelude | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Melody, Guardian Trek 3 & Trek 4, P-11 | 12-12 | 2 | LS24/4 |
| | Freestyle Mini, Zippie Z-Boç | U1 or 50-12 | 2 | M-D-EN0801, M-L-12090 |
| | Aspire M10/M11, Freestyle M11 | U1 or 40-12 | 2 | M-D-EN0801, M-L-12090 |
| Ranger All Seasons | 1x3, 2x3 | U1 | 4 | LS24/4, M-D-EN0801, M-L-12090 |
| | Solo, Safari | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| Redman Wheelchairs | 107SRX | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | 107SR | 22NF or 24FT | 2 | M-D-EN0801, M-L-12090 |
| Shoprider Mobility | Streamer (888W, 888WA), Sovereign (888-3, 888-4), Sunrunner 3 & 4 (incl. Deluxe, | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Wizz (888WNL, 888WNLB), 6Runner 10 (888WNLM), Trooper | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | 6Runner 10 Deluxe | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | FPC (PHFW-1118, PHFW-1120), Jetstream L (888WAL), Jetstream M (888WAM) | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Streamer 888WS, Sprinter 889-3 XL | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | HD (888WHD), 6Runner 14 (888WNLL, 888WNLLHD), Sprinter XL3 Deluxe | 24 | 2 | M-D-EN0801, M-L-12090 |
| | Sprinter, 889-3, 889-4, Sprinter XL4 (8889B-4) | U1 or 40-12 | 2 | M-D-EN0801, M-L-12090 |
| | TE889DX2-4, 889DX4-4, Sprinter Deluxe | 22NF or 40-12 | 2 | M-D-EN0801, M-L-12090 |
| | Flagship, Sprinter Jumbo XL, Sprinter XL4 & XL3, Sprinter XL4 Deluxe | U1, 22NF or 24 | 2 | M-D-EN0801, M-L-12090 |
| | TE999 | 40-12 | 2 | M-D-EN0801, M-L-12090 |
| | Scootie & Scootie Jr., Traveler, Smartie, Jimmie (UL8PB, UL8PBS), XtraLite Jiffy | 12-12 | 2 | LS24/4 |
| | Start, Dasher 3 & 4, Xtralite 3 & 4, E-Z Chopper, Sunrunner S | 12-12 | 2 | LS24/4 |
| | Sunrunner, Snazzy, 777E, Dasher 9, Start E | 17-12 | 2 | LS24/4 |
| | Echo | 10-12 | 2 | LS24/4 |
| | Suiter | Destination | 17-12 | 2 |
| Suntech | Sterling, Scoota, Regent 3 & 4 (new series) | U1 or 22NF | 2 | M-D-EN0801, M-L-12090 |
| | Regent 3 & 4 (early series) | 40-12 | 2 | M-D-EN0801, M-L-12090 |
| | Scoota Bug, Indigo 3 & 4 | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| Teftec | OmegaTrac, BetaTrac, AlphaTrac | 24 FT | 2 | M-D-EN0801, M-L-12090 |
| Theradyne | EV1782, EV1852, EV1882 | 40-12 | 2 | M-D-EN0801, M-L-12090 |
| | EV1868, EV1872 | 40-12 | 2 | M-D-EN0801, M-L-12090 |
| | EV1776, EV1866 (pediatric models) | 12260 | 2 | LS24/4 |
| | EV1864 | 12V1.9 | 2 | N/A |
| | Rover | 24 | 2 | M-D-EN0801, M-L-12090 |
| Tuffcare | Escort, Limo, Challenger | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |
| | Challenger Extra Wide | 24 | 2 | M-D-EN0801, M-L-12090 |
| 21st Century | Bounder, Bounder Plus, Big Bounder | 27 | 2 | M-D-EN0801, M-L-12090 |
| Tzora | Classic, Elite, Lite | 7-12 or 12-12 | 2 | LS24/4 |
| | Titan | 12260 | 2 | LS24/4 |
| Wheelcare | Superlight Scooter | 17-12 | 2 | LS24/4 |
| | Sportster Vehicle | 24 | 2 | M-D-EN0801, M-L-12090 |
| | Sportster Stretch | 27 | 2 | M-D-EN0801, M-L-12090 |
| Wheelchairs of Kansas | BCWPRAD- Power Advantage | 22NF | 2 | M-D-EN0801, M-L-12090 |
| | BCWPR- Power Chair, PEV/Overlander | 27 | 2 | M-D-EN0801, M-L-12090 |
| | Prairie Cruiser | 24 | 2 | M-D-EN0801, M-L-12090 |
| Win Med | Tourist Xpress | 12-12 | 2 | LS24/4 |
| | Tango Elite | U1 | 2 | LS24/4, M-D-EN0801, M-L-12090 |

Check our website at www.mkbattery.com for periodic updates to this chart.

© Copyright 2009 MK Battery