

Cruiser Features

HCPCS:E1236

CRUISERS

- **30° Fixed Tilt**
Promotes upper body positioning for individuals with limited upper body strength or trunk stability
- **Self Tensioning Seating Surface**
Frame is uniquely designed to open wider with the weight of the occupant – this design keeps the fabric taut, and provides a flat/stable seating surface for optimal positioning
- **85° – 95° Seat-to-Back Angle Adjustment**
Enhanced positioning for individuals with limited ability to bend at the waist
- **One-Piece Adjustable Push Handle**
For easy maneuvering
- **Under Seat Storage Basket & Side Bags**
Provides extra storage for quick trips to the park
- **Cruiser EI - Size 10**
Early Intervention (EI) wheelchair designed to address the positioning needs of children under 50 lbs. (suggested age range from 6 months to 4 years)
- **Frame Folds Compactly**
For easy storage and transport – fold and go!
- **2 Upholstery Options Available**
Padded Cordura for added comfort or Textilene for breathability
- **In-Depth Positioning Options**
Provide optimal positioning for increased therapeutic benefits – all day seating
- **Crash Tested Transit Models**
Adhere to national safety standards, and provide safe transport in bus or van



Cruiser 18 (Cordura)

- 2-Piece Push Handle
- 8" Front and 11" Rear Tires
- Silver Frame
- Adult Style Cordura Upholstery
- Headrest Extension
- H-Harness with Padded Cover

Standard Solid Tires

CX10/12/14/16: 7.5" x 2" Front, 11.5" x 2.5" Quick Release Rear Wheels
CX18: 8" Front, 11" Rear

Positioning options and accessories sold separately.

Cruiser Specifications

Measurement	Chair Model	Seat Width*	Seat Depth**	Seat Back Height ^	Seat to Footplate*	Seat to Floor*	Seat-to-Back Angle	Seat Angle	Headrest Extension	Shoulder Strap Height ◊	Weight Capacity Regular/Transit	Weight of Chair Δ
in.	CX 10	10	6 - 11	21 C / 20 T	5 - 21	23	85° - 90° - 95°	30°	6	11 - 19	75 lbs / 66 lbs	27 lbs
cm	CX 10	25.5	15.2 - 27.9	53 C / 51 T	13 - 53.4	58			15	28 - 48	34 kg / 30 kg	12 kg
in.	CX 12	12	8 - 13	22 C / 20 T	5 - 25	24	85° - 90° - 95°	30°	8	12 - 20	75 lbs / 66 lbs	27 lbs
cm	CX 12	30.5	20 - 33	56 C / 51 T	13 - 63.5	61			20	30.5 - 51	34 kg / 30 kg	12 kg
in.	CX 14	14	10 - 15	25 C / 23 T	4 - 22	22	85° - 90° - 95°	30°	8	15 - 23	100 lbs / 100 lbs	28 lbs
cm	CX 14	35.5	25.4 - 38	64 C / 58 T	10 - 56	56			20	38 - 58	45.5 kg / 45.5 kg	13 kg
in.	CX 16	16	11 - 16	27 C / 24 T	4 - 22	23	85° - 90° - 95°	30°	8	17 - 25	170 lbs / 170 lbs	32 lbs
cm	CX 16	40.5	27.9 - 40.6	71 C / 61 T	10 - 57	58.4			20	43 - 63.5	77 kg / 77 kg	14 kg
in.	CX 18	18	14 - 21	25 C / 25 T	8 - 23	26	85° - 90° - 95°	30°	10	16 - 22	250 lbs / 170 lbs	30 lbs
cm	CX 18	45.5	35.6 - 53.3	64 C / 64 T	20.3 - 58	66			25.4	41 - 56	144 kg / 77 kg	14 kg

in. - Inches, cm - Centimeters

* Plus or minus .5 inch (1.27 cm) due to fabric. ** 2 inch cut-out for minimum seat depth. Available as an option. ◊ - Models with Textilene Fabric, Shoulder Strap Height is 2 inches shorter.

^ Seat Back Height: (C) Cordura Upholstery / (T) Textilene Upholstery. Δ - Weight of chair is without front rigging.

Measurement	Chair Model	Overall Width	Overall Height	Overall Length	Push Handle Adjustments °	Folded Length	Folded Width	Folded Height
in.	CX 10	22	37	48	28 - 43	47	15	17
cm	CX 10	56	94	122	71 - 109	119	38	43
in.	CX 12	23.5	38	44	29 - 46	48	15	17
cm	CX 12	60	97	112	74 - 117	122	38	43
in.	CX 14	25.5	37	49	28 - 45	47	16	17
cm	CX 14	65	94	124.5	71 - 114	119	41	43
in.	CX 16	27.5	40	52	31 - 48	53	16	17
cm	CX 16	70	102	132	78 - 122	135	41	43
in.	CX 18	28.5	38	55	N/A	50	15	17
cm	CX 18	72.4	97	140	N/A	127	38	43

° Adjustable Push Handle is not available on 18 inch models

Use only Convoid options and parts on Convoid products. Convoid parts are not interchangeable with other manufacturer's products. Replace any worn parts immediately.

Note: Convoid maintains a policy of continual product improvement; therefore, we reserve the right to make changes without notice.

Limited Lifetime Warranty

Frame and X-Braces.....Lifetime of original retail buyer
Fabric and webbing..... One year
Other components.....One year

This warranty does not cover normal wear and tear, or damage caused by accident or misuse.



Upholstery Styles & Colors

Removable, Washable Fabric

Textilene



Cordura



Textilene



Cordura