

DASH3

FLIGHTWEIGHT
BY TRAVEL BUGGY



Less
weight
More
value



Upgraded
comfort
seating



Smoother,
quieter
power delivery



OUR LIGHTEST, SMARTEST DASH YET

Meet the
DASH3

The DASH has always been a fan favourite, and now it's better than ever. Introducing the DASH 3: **lighter, smarter, and ready to go where life takes you.**

KEY FEATURES

4 LBS LIGHTER. BUILT STRONGER.

Featuring magnesium components for improved durability and performance.

ENHANCED PRECISION & CONTROL

Experience smoother, more responsive control with the upgraded joystick and CPU system.



A FRESH NEW LOOK

Finished in a sleek dark grey metallic frame, the DASH 3 looks as sharp as it performs.



ROLLS LIKE A SUITCASE



POWER YOU CAN FEEL, NOT HEAR

New tech motors deliver quieter, smoother performance, while magnesium rear rims reduce vibration.

BUILT-IN STORAGE



WIDER, CONTOURED ARMRESTS



STRONGER, EXTENDABLE FOOTPLATE



YOUR DRINK, WITHIN REACH



SPECIFICATIONS

MODEL SKU	16MNA0060100
FRAME COLOUR	Dark Grey Metallic
NET WEIGHT	42 lbs
GROSS WEIGHT	46 lbs
WEIGHT CAPACITY	265 lbs
DRIVING RANGE	15 km (9.5 mi)
INCLINE CAPABILITY	10°
TURNING RADIUS	24"
MAXIMUM SPEED	6 km/h (3.7 mph)
MOTOR POWER	2 x 170 W Brushless
UNFOLDED DIMENSIONS	38"L x 24"W x 37"H
FOLDED DIMENSIONS	30"L x 24"W x 13.5"H
FRONT / REAR WHEEL DIAMETER	8" / 12"

SEATING DIMENSIONS

SEAT SIZE	17.5"W x 17"D x 2"H
BACKREST HEIGHT	19"
BETWEEN ARMRESTS	18.5"
SEAT TO FLOOR HEIGHT	20"

BATTERY

BATTERY TYPE	Li-Ion 24 V / 10 Ah
BATTERY INCLUDED	One (1)

WARRANTY

FRAME	5 Years
MOTOR	2 Years
CONTROL SYSTEM	1 Year
BATTERY	1 Year
PARTS	3 Months