

# 2023 Evaluation Form

# DCP12

CREATING A CYCLE OF mobility




Inseam Range: 13" - 21"

The DCP12 allows kids ages 4-8 to get out and enjoy adaptive biking with their friends. This bike provides excellent stability and fits through standard doorways, making it a popular choice for homes, schools, and therapy clinics.

DCP12 MODEL SPECIFICATIONS	
Rider Age	4-8
Max Load (lbs)	125
Bike Weight (lbs)	45
Frame Height	41"
Total Length	54"
Bike Width	29"
Step-Through Height	14"

Measurements Include Rear Steering Handle and High Back Seating

Sizing Guide: Inseam 13" - 21"



## Tips For Choosing an Adaptive Bike

### 1. Inseam Measurement

The measurement we use to ensure each rider gets the correct size of adaptive bicycle is called an inseam. To find an inseam, measure the distance between the bottom of the rider's heel to the top of the rider's inner thigh. Make sure the rider's inseam measurement matches the inseam range for the adaptive bicycle.

### 2. Know Your Accessories

We offer many different accessories to help people of all abilities experience freedom. Make sure to look through all of our available accessories to make sure the adaptive bike will work for the rider. More information on all of our accessories is available at: [www.freedomconcepts.com/product-lines/adaptive-bikes/accessories/](http://www.freedomconcepts.com/product-lines/adaptive-bikes/accessories/)

### 3. We Do Custom

If none of our current models of adaptive bicycles work for you, we are able to custom build a bicycle that will. Please contact us at 1-800-661-9915 for more information.

## RIDER INFORMATION

Rider's Name: \_\_\_\_\_ Male Female Height: \_\_\_\_\_ Inches \_\_\_\_\_ cm

Date of Birth: \_\_\_\_\_ Weight: \_\_\_\_\_ Lbs \_\_\_\_\_ Kgs



Inseam: \_\_\_\_\_




Diagnosis: \_\_\_\_\_




# DCP12

- Inseam Range: 13" - 21" - SHADED ITEMS INCLUDED IN BASE MODEL PRICE





MODEL		FRAME ADD-ONS - Select one (Optional)	
			
<p>The standard DCP12 model accommodates an inseam range of 22"-29".</p>		<p>The Wide Glide Frame provides the rider with more stability. This model is typically selected for riders with more control over the steering. This model also fits through standard doorways with a width of 33" and length of 43".</p>	
<input checked="" type="checkbox"/> DCP12 (DCP12) <span style="float: right;">\$4,075</span>		<b>Wide Glide Frame</b> <span style="float: right;">\$215</span> <small>*Adds 4" to width of frame (WIDE FRAME)</small>	<b>Wide Glide Extended Frame</b> <span style="float: right;">\$330</span> <small>(WGEXT)</small>

FRAME COLORS- Select one			
Standard Options		Additional Options	
			
<b>Candy Blue</b> (CANDY BLUE) (STD)		<b>Arctic Teal</b> (TEAL) <span style="float: right;">\$110</span>	<b>Midnight</b> (BLACK) <span style="float: right;">\$110</span>
<b>Hot Rod Purple</b> (HOT ROD PURPLE) (STD)		<b>Blaze Orange</b> (ORANGE) <span style="float: right;">\$110</span>	<b>Snow</b> (WHITE) <span style="float: right;">\$110</span>
<b>Electric Green</b> (GREEN) (STD)		<b>Radiant Pink</b> (PINK) <span style="float: right;">\$110</span>	<b>Fossil</b> (FOSSIL) <span style="float: right;">\$110</span>
<b>Cherry Red</b> (RED) (STD)		<b>Chrome</b> (CHROME) <span style="float: right;">\$110</span>	
<b>Custom Color</b> (CUSTOM COLOR) <span style="float: right;">\$260</span> Please indicate custom color note here:		<b>Color combination</b> (COLORCOMBO) <span style="float: right;">\$125</span>	
		Please indicate primary frame color:	
		Please indicate secondary accent color:	

DECAL OPTIONS - Select one		NAMEPLATE (Optional)	BELL & SQUEAKY HORN (Select one Bell & one Horn)			
<b>Racer</b> (STD) (6401232) 		<b>Custom Nameplate</b> <b>Nameplate</b> (NAMEPLATE SP) <span style="float: right;">\$95</span> Notes: 	<b>Bells</b> (STD)	<b>Squeaky Horns (Optional)</b> (STD)		
<b>Hibiscus</b> (STD) (6400250) 			<b>Red</b> (5015014) <b>Blue</b> (5015011) <b>Purple</b> (5015013) <b>Gold</b> (5015012) <b>Black</b> (5015010)	<b>Dog</b> (5015063) <b>Frog</b> (5015065) <b>Cat</b> (5015061) <b>Parrot</b> (5015070) <b>Orca</b> (5015067)	<b>Dino</b> (5015062) <b>Penguin</b> (5015068) <b>Unicorn</b> (5015071) <b>Tiger</b> (5015069) <b>Shark</b> (5015072)	<b>Turtle</b> (5015074) <b>Giraffe</b> (5015075)

TIRE SELECTION - Select one		
		
<b>Smooth Solid Tires</b> (STD) <small>(SMOOTH SOLID TIRES)</small> Ideal for indoor and outdoor use, with limited speeds and/or with caregiver present		
	<b>Knobby Pneumatic</b> <span style="float: right;">\$110</span> <small>(2112252)</small> Ideal for more independent riders, wanting to ride at faster speeds	
		<b>Smooth Pneumatic</b> <span style="float: right;">\$110</span> <small>(2112251)</small> Ideal for more independent riders, wanting to ride at faster speeds


## HANDLEBAR OPTIONS - Select one

			
<b>Hoop Handlebars</b> (STD) (2111001)	<b>High Rise Handlebars</b> \$60 (2111005)	<b>Enclosed Hoop Handlebars</b> \$70 (2111003)	<b>Wide Handlebars</b> \$60 (2111007)

## STEERING OPTIONS - Select setup

<b>Unlocked Steering</b> (STD) (STEERUNLOCKED) <i>This option allows the rider to control the steering.</i>	<b>Locked Steering</b> (N/C) (STEERLOCKED) <i>This option takes steering control away from the rider so the caregiver can have full control of the bike with the Rear-Steering Handle.</i>
<small>*Please note you can change the steering from unlocked to locked by using the Freedom Concepts multi-tool that comes with every bike.</small>	

## REAR STEERING


<b>Rear Steering Handle</b> \$685 (REAR STEER)


## REAR STEERING ADD-ONS

	
<b>Bicycle Handlebars</b> \$180 (2125006) <small>*Requires Rear Steering Handle</small>	<b>Hoop Handlebar</b> \$180 (2125005) <small>*Requires Rear Steering Handle</small>

## BRAKE OPTIONS AND POSITIONS

Select First Brake Position (Included)	Select Second Brake Position (Optional)
<b>*Please select only one option</b> (STD)	<b>*Please select only one option - must differ from first brake position</b>
<b>Left Handlebar</b> (2113216) *Not available with High Rise Handlebars	<b>Left Handlebar</b> (2113212) *Not available with High Rise Handlebars \$110
<b>Right Handlebar</b> (2113217) *Not available with High Rise Handlebars	<b>Right Handlebar</b> (2113213) *Not available with High Rise Handlebars \$110
<b>Rear Steering Handle</b> (2113215) *Rear Steering Handle Required	<b>Rear Steering Handle</b> (2113211) *Requires Rear Steering Handle \$110
<b>T-bar</b> (2113214)	<b>T-bar</b> (2113210) \$110

## PROPULSION - Please fill out one column



<b>Direct Drive</b> (STD) (DIRECT DRIVE) <i>When the bike is moving, the pedals are moving.</i>	<b>Free Wheel</b> \$80 (FREE WHEEL) <i>Allows legs to rest while the bike coasts forward.</i>	<b>On/Off Pin Reciprocator</b> \$235 (ON-OFF) <i>Alternate from Free Wheel mode or Direct Drive by engaging the On/Off Pin Reciprocator.</i>
<b>Pedal Power:</b> Extra Easy (EXTRA_EASY) Easy (EASY) Average (AVERAGE) Challenging (CHALLENGING)	<b>Pedal Power:</b> Easy (EASY) - 3" crank only Average (AVERAGE) Challenging (CHALLENGING)	<b>Pedal Power:</b> Easy (EASY) Average (AVERAGE) 

# Seating System



Please Fill Out Either The Child Set Up or The Adult Set Up  
\*Adult Seats are only available with Wide Glide frames\*

Child Seat	12.5" width, 12.5" depth
Child High Back	13" width, 18" height
Child Low Back	13" width, 7.5" height
Adult Seat	16.75" width, 14" depth
Adult Low Back	18" width, 10" height
Adult High Back	18" width, 18" height



Please note that all seat bases come with a seat belt and all highback backrests come with a chest strap



## CHILD SET UP





	
<b>Seat Base</b>	
<b>Child Seat</b> (STD) (2114060)	
<b>Child Narrow Seat</b> \$110 (6114014)	
<b>Child Pommel Seat</b> \$110 (6114010)	
<b>Child Narrow Pommel Seat</b> \$215 (6114016)	
<b>Four-Point Neoprene Hip Stabilizer Belt</b> \$200 (2125273) *OPTIONAL	
<b>Back Rest</b>	
<b>Child Low Back</b> (STD) (2114009)	
<b>Child High Back</b> \$110 (2114010)	

## ADULT SET UP

	
<b>Seat Base</b>	
<b>Adult Seat</b> \$110 (6114015)	
<b>Adult Narrow Seat</b> \$110 (6114017)	
<b>Adult Pommel Seat</b> \$110 (6114018)	
<b>Adult Narrow Pommel Seat</b> \$215 (6114020)	
<b>Four-Point Neoprene Hip Stabilizer Belt</b> \$200 (2125274) *OPTIONAL	
<b>Back Rest</b>	
<b>Adult Low Back</b> \$60 (2114019)	
<b>Adult High Back</b> \$110 (2114020)	







## SEATING ACCESSORIES (Optional)

<b>Neoprene Chest Harness</b> 	<b>10 Degree Fixed Tilt</b> 	<b>Laterals</b> ( <i>Width and height adjustable, comes with an adjustable track system</i> ) 	<b>Flip-Up Abductor</b> 
<b>Medium</b> \$120 (6314025)  <b>Large</b> \$120 (6314026)	<b>10 Degree Fixed Tilt</b> \$195 (TILT SEAT)  <b>*This option moves the bottom bracket position forward 1".</b>	<b>Pair of Standard Small Laterals</b> \$710 <i>Small Laterals Chest Width: 8½" - 11¼"</i> (2114070)  <b>Pair of Standard Large Laterals</b> \$710 <i>Large Laterals Chest Width: 12"-15"</i> (2114072) Chest Width: _____ Seat-to-Armpit: _____  <b>Laterals Add-On:</b> <b>Narrow Laterals</b> \$145 (2125033) <i>*Shifts chest width range 2" narrower</i>	<b>Flip-Up Abductor</b> \$360 (2020205)  <i>This flip-up abductor assists with leg/knee separation during pedaling and is mounted securely to the underside of the seat.</i>



## HEADRESTS (Optional)

		
<b>Adjustable Backrest Extension</b> \$210 (Backrest Ext)  <i>The backrest extension extends the seat backing height an additional 7" and can be removed for ease of transport.</i>  <i>Tip: A Low Back Seat and the Adjustable Backrest Extension can be combined to create a height-adjustable backrest.</i>	<b>Adjustable Headrest with Occiput Support Cushion</b> \$390 (HEADREST PLUS) <i>*Requires Highback Seat</i>  <i>The adjustable headrest with an occiput support cushion provides additional neck support for the rider. The cushion comes with velcro and can be removed or adjusted as required.</i>	<b>Head and Neck Support System</b> <i>*Requires Highback Seat</i>  <b>*Please call 1-800-661-9915 for sizing and part number. In some cases we are able to adapt our seatback to accommodate common wheelchair headrests.</b>



## PEDAL OPTIONS - Select one

			
<i>Mountain Pedals are recommended for riders that have more control and high tone muscle control.</i>	<i>The Neoprene adjustable footplate is a quick and easy way to maintain proper positioning and pedal pressure. Recommended for children with low muscle tone.</i>	<i>Molded ABS Footplates are designed to provide precise positioning for children with more spasticity. They also help with leg alignment and knee flexion.</i>	<i>Molded ABS Footplates with Ratchet Straps helps to secure the rider's foot in a specific position. Recommended for children with plantar flexion.</i>
<b>Mountain Pedals with Toe Cage</b> (STD) \$211 (2112081)	<b>Neoprene Adjustable</b> \$375 SML (2012031) *4" x 7" MED (2012032) *4½" x 10" LRG (2012033) *5" x 11½"	<b>Molded ABS</b> \$425 SML (2012042) *3" x 7" MED (2012043) *3" x 10¾" LRG (2012044) *3½" x 11½"	<b>Molded ABS with Ratcheting Straps</b> \$545 SML (2012052) *3" x 7" MED (2012053) *3" x 10¾" LRG (2012054) *3½" x 11½"


## PEDAL ADD-ONS - (Optional)



	
<b>1/4 Inch Footplate Spacers</b> \$5 each (2412050) Please add how many spacers are needed below. (Maximum 8 per side) _____ left + _____ right = _____ total	<b>Pulley String</b> \$35 <i>Front pulley assists with plantar flexion and rear pulley helps self-level the pedals.</i> <b>Front</b> (2115020) <b>Rear</b> (2115025)




## CHAINGUARD - Select one

	
<b>Standard Chainguard</b> (STD) \$230 (STD CHAINGUARD)  <b>*The Chainguard color will match the frame.</b>	<b>ABS Enclosed</b> \$230 (ABS ENCLOSED)  <b>*Recommended for multi-user environments.</b>








CRANK ARMS - Select one	
<b>3" Alloy Crank (STD)</b> (5012252) *Recommended	 <b>Crank Arm</b>
<b>4" Alloy Crank (N/C)</b> (5012254)	
<b>Crank Arm Length Shortener (optional)</b> *Reduces length by 1/2"	
<b>Left Side (2012049)</b>	<b>\$80</b>
<b>Right Side (2012050)</b>	<b>\$80</b>
<b>Both Sides (2012047)</b>	<b>\$160</b>

FRONT FORK - Select one	
	
<b>Rigid Chrome Fork (STD)</b> (5011012)	<b>Adjustable Tilt Fork (2111261)</b>
	<b>\$700</b>
<i>The adjustable tilt fork allows therapists and caregivers to add or reduce the amount of tilt (up to roughly 10 degrees) this helps to improve the rider's overall comfort and positioning.</i>	

ACCESSORIES - (Optional)				
				
<b>Velcro Hand-Eye Coordination Gloves</b>	<b>Indoor Stationary Trainer</b>	<b>Leg Braces</b>	<b>Rear Platform</b>	<b>Adjustable Sun Shade Canopy</b>
<b>\$95</b>	<b>\$560</b>	*Only available with molded footplates. *Requires Rear Pulley.	<b>\$420</b>	<b>\$365</b>
<b>XS (6315031)</b> 5 1/2" height x 1 3/4" wrist	<b>SML (6315033)</b> 6" height x 2" wrist	<b>Left Side Only (LEG BRACE L)</b>	<b>(REAR PLATFORM)</b>	<b>(RAIN &amp; SUN CANOPY)</b>
<b>MED (6315034)</b> 7" height x 2 1/4" wrist		<b>Right Side Only (LEG BRACE R)</b>		*Requires the Adjustable Headrest.
		<b>Pair (LEG BRACE PAIR)</b>		
		<b>\$585</b>		
		<b>\$585</b>		
		<b>\$1,170</b>		
		Inseam: _____ Calf Circumference: _____ *Measurement with AFO's*		

ACCESSORIES - (Optional)				
				
<b>Spoke Covers (3)</b>	<b>Adjustable Tablet/Phone Holder</b>	<b>Cane/Crutch Holder</b>	<b>Padded Adduction Strap</b>	<b>Oxygen Tank Holder</b>
<b>\$120</b>	<b>\$210</b>	<b>\$320</b>	<b>\$125</b>	<b>\$370</b>
<b>(2112240)</b>	<b>Tablet (5025070)</b> <b>Phone (5025080)</b>	<b>(2125060)</b>	<b>(6320060)</b>	<b>(2415070)</b>
				<b>*Fits Size D/E Tanks</b>

ACCESSORIES - (Optional)				
				
<b>Communication Tray</b>	<b>Simulated Sheepskin Cover</b>	<b>Mirror</b>	<b>Tool Pouch</b>	<b>Bike Cover</b>
<b>\$360</b>	<b>Child Seat (6114070)</b>	<b>\$65</b>	<b>\$35</b>	<b>\$150</b>
<b>(2125281)</b>	<b>\$120</b>	<b>(5015040)</b>	<b>(6300022)</b>	<b>(6315600)</b>
	<b>Adult Seat (6114080)</b>			
	<b>\$140</b>			

**CUSTOMER INFORMATION**

Dealer/Vendor: \_\_\_\_\_  
Contact Person: \_\_\_\_\_  
Address: \_\_\_\_\_  
City: \_\_\_\_\_ Prov./State: \_\_\_\_\_  
Postal/Zip Code: \_\_\_\_\_ Ordered by: \_\_\_\_\_  
Phone: \_\_\_\_\_ Fax: \_\_\_\_\_  
Email: \_\_\_\_\_

**PRIVACY POLICY**

Freedom Concepts Inc. collects, uses and discloses personal information of the Customer and/or Rider to properly assemble and manufacture the mobility device and may contact the Rider's therapist for additional personal information to assist with this process. The collection, use and disclosure of the personal information of the Customer and/or the Rider is governed by any consent which may, from time to time, be provided, and by the Freedom Concepts Inc. Privacy Policy which is available on our website www.freedomconcepts.com or by contacting our office. By signing below, you hereby warrant that you have the authority to bind the Rider and hereby authorize the collection, use, and disclosure of personal information as herein provided and in accordance with the Freedom Concepts Inc. Privacy Policy.

**ORDER INFORMATION**

Purchase Order: \_\_\_\_\_  
Order Date: \_\_\_\_\_ Quote Date: \_\_\_\_\_  
FCI Representative: \_\_\_\_\_  
Payment Terms: \_\_\_\_\_  
CREDIT CARD (CHECK ONE)  Visa  MasterCard Exp. Date: \_\_\_\_\_  
Card Number: \_\_\_\_\_  
Name On Card: \_\_\_\_\_  
Signature: \_\_\_\_\_

**TERMS AND CONDITIONS**

Prices are subject to change without notice. Prepaid or net 30 days O.A.C. All taxes extra. Freight extra. Manufacturer reserves the right to change prices, colors, features, dates, locations and products without notice. Residential delivery charge may apply. Duty and brokerage included where applicable.

**SHIPPING INSTRUCTIONS**

Name: \_\_\_\_\_  
Address: \_\_\_\_\_  
City: \_\_\_\_\_  
Province/State: \_\_\_\_\_  
Postal/Zip Code: \_\_\_\_\_ Phone: \_\_\_\_\_  
Federal ID/SS#: \_\_\_\_\_

**PHOTO RELEASE**

I hereby give permission to Freedom Concepts Inc. to use photographs of me or my child for the sole purpose of the company's promotional use, and will not be distributed to third parties.  YES  NO

I would like to receive regular updates.  SNAIL MAIL  EMAIL  NO

\_\_\_\_\_  
CUSTOMER SIGNATURE

\_\_\_\_\_  
SIGNATURE OF RIDER/PARENT

Questions?  
1-800-661-9915

**NOTES:**

[Empty rectangular box for notes]