



Diamond Deluxe Aluminum Rollator/Transport Chair

745R Deluxe, Cherry Red, 1/cs

745B Deluxe, Marine Blue, 1/cs

FEATURES:

- Transport chair and rollator in one.
 - Easily converts from a rollator to a transport chair
 - Comes with Tool Free adjustable footrests
 - Tool free removable front and rear wheels.
 - Deluxe, durable backrest.
 - Comes with removable carry pouch.
 - Large 8" casters are used for indoor & outdoor use.
 - Comfortable mesh seat.
 - Handles are height adjustable.
 - Soft, comfortable grips.
 - Easy to use loop locks.
 - Limited Lifetime Warranty.
- Lifetime Warranty on brake cable.**



RETAIL PACKAGED!

USE AS A TRANSPORT CHAIR... OR AS A ROLLATOR!



SPECIFICATIONS

Item #	Handle Height Highest	Handle Height Lowest	Unit Width	Unit Length	Seat (depth)	Seat (width)	Seat (height)	Product Weight	Carton Shipping Weight	Weight Limit
DIAMOND DELUXE	37.5"	31.5"	24"	25.5"	9"	18"	21"	18 lbs.	23 lbs.	300 lbs..

drive MEDICAL DESIGN & MANUFACTURING
 rollators
 To order call toll free: 877.224.0946