

Cirrus Plus EC

HCPCS: K0822-1, K0822



FIGURE A

- Foldable, lightweight frame for easy portability (Figure A)
- Attractive steel cross-brace frame
- Armrests are padded, adjustable and removable
- Adjustable-length controller mount
- Adjustable-tension back upholstery
- Swing-away footrests with heel loops
- Optional elevating leg rests (Item #CPELR) sold separately
- Positioning belt included
- Padded back upholstery with storage pouch and cushioned seat
- Programmable controller
- Flat-free, non-marking tires
- Easy free-wheel operation
- Manual wheel locks
- General foam seat cushion included
- Anti-tip wheels

PRODUCTS	Description
CPN16FBA	Cirrus Plus EC, CS
CPN18FBA	Cirrus Plus EC, CS
CPN20FBA	Cirrus Plus EC, CS

SPECIFICATIONS	Description
Alarm mute	Yes
Base Weight	90 lbs.
Batteries (pair)	12V x 36AH (U1)
Battery Weight (pair)	56 lbs.
Brakes	Electromagnetic
Caster Wheels	8" x 2" Solid
Charger	4A Offboard
Climbing Angle	6°
Controller	PG VSi 50A
Dimensions	42.5" (L) x 24" (W) 42.5" (L) x 26" (W) 42.5" (L) x 28" (W)
Drive Wheels	12.5" x 2.25" Flat-free
Floor to Seat Height	18.5"
Ground Clearance	4"
HCPCS	K0822

SPECIFICATIONS	Description
Low pressure alarm	Yes
Mattress Base Material	Zippered vinyl
Mattress Dimensions	80" (L) x 60" (W) x 10" (H)
Mattress Maximum Support Capacity	1000 lbs.
Maximum Range	15 miles
Memory function	Yes
Motors (2)	24V x 320W x 4600 rpm
Power failure alarm	Yes
Pump Adjustable cycle time	10, 15, 20, 25 minutes
Pump Airflow	12-LPM
Pump Auto firm	60 mmHg, Auto return in 30 mins
Pump Power	120 VAC, 60 Hz
Pump Dimensions	14.2" (W) x 9.8" (H) x 4.5" (L)
Pump Function lock	Auto lock in 30 seconds
Pump Power Cord Length	15 Ft
Pump Seat function	Pressure preset + 10 mmHg
Pump Pressure Range	20 mmHg to 60 mmHg
Pump UL	60601-1, 1st Edition-2003
Seat Dimensions	16" (W) x 16" (D) 18" (W) x 16" (D) 20" (W) x 16" (D)
Static selectable with auto-return to alternating mode	Yes
System (alternating) alarm	Yes
Top Speed	5 mph
Total Weight (excl. hand control)	146 lbs.
Turning Radius	30" 32" 33"
Warranty on Batteries	12 Months
Warranty on Drive Train Components	24 Months
Warranty on Electronic Controller	24 Months
Warranty on Frame	Lifetime
Warranty on Pump	Limited 18 Months
Warranty on Mattress	Limited 18 Months
Weight Capacity	300 lbs
Weight Capacity (13046 13046-HD)	300 lbs.

PARTS	Description
STDS871	Hand Grip,Sentra,1/ea, EA
CPN6205	Controller Arm Brkt.CPN,1/ea, EA
CP6200T	JoystickTip,PG Controller,1/ea, EA
SI69001	ContrllrScrewSet,P&G,VSI,1/ea, EA
CP8002-16	Seat Upholstery, CPN16, 1/ea, EA
CP6205P	Controller Brket,CirrusPlus,1/ea, EA
CP8002-18	Seat Upholstery, CPN18, 1/ea, EA
CP8002-20	Seat Upholstery, CPN20, 1/ea, EA
CP8002-22	Seat Upholstery, CPN22, 1/ea, EA
CP8001-16	Seat Cushion, CPN16, 1/ea, EA
CP8001-18	Seat Cushion, CPN18, 1/ea, EA
CP8001-20	Seat Cushion, CPN20, 1/ea, EA
CP8001-22	Seat Cushion, CPN22, 1/ea, EA
CP6039	CalfSupportStrap,16"Cir.Pwr.1/ea, EA
CP8039	CalfSupportStrap,18"Cir.Pwr.1/ea, EA
CP2039	CalfSupportStrap,20"Cir.Pwr.1/ea, EA
CPNSFR	Footrest Rt.CPN,1/ea, EA

PARTS	Description
CPNSFL	Footrest Lt.CPN,1/ea, EA
STDS1026	Front Bumper for Ft.Rst/LegRest,1/ea, EA
CP8020-16	CrossBar16"CirrusPlus1/ea, EA