

Phoenix HD 3

HCPCS: K0801-1, K0801



FIGURE A



FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Quick and easy disassembly
- Includes interchangeable color panels in Red and Blue (Figure A)
- Armrests are padded and adjustable
- Quick connect batteries
- Ergonomic throttle control
- Height-adjustable swivel seat with two-tone upholstery
- Flat-free, non-marking tires
- Large, plastic carry basket
- Easy-to-adjust tiller
- Anti-tip wheels
- Headlight
- Swivel seat with fold-down backrest
- Easy free wheel operation

PRODUCTS

	Description
PHOENIXHD3	Phoenix HD 3, EA
PHOENIXHD3-20	Phoenix HD 3, EA

SPECIFICATIONS

	Description
Base Weight	63 lbs.
Batteries (pair)	12V x 20AH
Battery Weight (pair)	31 lbs.
Brakes	Electromagnetic
Charger	1.5A Offboard
Climbing Angle	6°
Dimensions	39" (L) x 22.5" (W)
Drive System	24(V) x 350(W) x 3600 rpm
Floor to Seat Height	20" - 22"
Front Wheel	9" x 3" Flat-free
Ground Clearance	3.75"
HCPCS	K0801
Heaviest Piece Weight	32 lbs. (Front Section)
Maximum Range	15 miles
Motor	24V x 350W x 3600 rpm
Rear Wheels	9" x 3" Flat-free
Seat Dimensions	17.5" (W) x 17" (D)

SPECIFICATIONS	Description
Seat Type	Folding
Seat Weight	21 lbs.
Top Speed	4 mph
Total Weight (excl. hand control)	115 lbs.
Turning Radius	43"
Warranty on Batteries	12 Months
Warranty on Drive Train Components	24 Months
Warranty on Electronic Controller	24 Months
Warranty on Frame	Lifetime
Weight Capacity	350 lbs
Weight Capacity (13046 13046-HD)	350 lbs.
Wheel Base Dimensions (L x W)	29.1" x 19.1"

ACCESSORIES	Description
SF8050	Sun Shade for Scooter,1/ea, EA
SF8020	Rear Basket Set, LType, 1/ea, EA
SF8010	Oxygen Holder Set, L Type,1/ea, EA
SF8030	Crutch Holder Set, SType, 1/ea, EA

PARTS	Description
29002W7011	AxleKey,Cobra,Prowler,PhoenixHD,Vent,1/e, EA
S31003	Grip,Phoenix,Odyssey,Day-3/4,1/ea, EA
S31006	Ignition Switch,Phoenix,Day-3/4,1/ea, EA
S31022	Front LampAssy,Phoenix,1/ea, EA
S31209BK	TillerCoverBlackPhoenix,1ea, EA
S31122	Headlight button Phoenix, Daytona,1/ea, EA
S31001	Speed Dial Assembly,Phoenix,1/ea, EA
S31012	Horn button, Phoenix, Vent, 1/ea, EA
S31208B	Battery Indicator,LED,Phoenix,1/ea, EA
S31147	Charger Port,Phoenix, Vent, 1/ea, EA
S31036-3HD	Fork,Phoenix-3HD,1/ea, EA
S350127	Front Wheel Cover,Phoenix3/HD3,1/ea, EA
S31142-3HD	Front Tire,Phoenix-3HD,1/ea, EA
S31011	Tow Handle,Phoenix,1/ea, EA
S350169715	Seat Post Knob, Phoenix / PhoenixHD, 1ea, EA
S350109	Tiller Swivel Assm.Phn3/4/HD3/HD4,1/ea, EA
S31028	Tiller Adj.Knob, Phoenix,1/ea, EA
S31021	Bearing Ball Holder Assembly,1/ea, EA
S350169716	T Handle Frame Lock,Phnx3/4/HD3/HD4,1/ea, EA
S350169717	Seat Adjustment Pin,Phnx3/4/HD3/HD4,1/ea, EA