

Super Light, Folding Transport Chair with Carry Bag and Flip-Back Arms

HCPCS: E1038



FIGURE A

- Padded flip-back armrests
- Fold-down footrests
- Weighs only 19 lbs
- Easy to push or transport
- Aluminum frame with lightweight nylon upholstery
- Composite wheels are lightweight and maintenance free
- Folds like a conventional folding chair (Figure A)
- Comes with carry bag (item #835) with shoulder strap (Figure B)
- Standard rear wheel locks and seat belt

PRODUCTS

PRODUCTS	Description
SL18	Super Light, Folding Transport Chair with Carry Bag and Flip-Back Arms, CS

SPECIFICATIONS

SPECIFICATIONS	Description
Back Height	17"
Carton Shipping Weight	24 lbs.
Overall Dimensions	23" (W) (open) x 37" (H)
Seat Dimensions	18" (W) x 17" (D)
Seat To Floor Height	20"
Seat Width	18"
Warranty	Limited Lifetime
Weight	19 lbs.
Weight Capacity	250 lbs

PARTS

PARTS	Description
SL182A17	Caster Lock Knob/Bolt,SL18,1/ea, EA
SL182A00	Front Caster, 6"SL18,1/ea, EA
SL182A02	Front Axle Bolt,SL18,1/ea, EA
SL182A01	Rear Caster,8"SL18,1/ea, EA
SL182A03	Rear Axle Bolt,SL18,1/ea, EA
SL182A04	Footrest,SL18,1/ea, EA
SL182A05	Frame Dustcap,SL18,1/ea, EA
SL182A06	Frame Spacer,SL18,1/ea, EA
SL182A07	Arm Pad,SL18,1/ea, EA
SL182A08	Seat Uphol.SL18,1/ea, EA
SL182A09	Rt.Armrest Mntng Brkt.SL18,1/ea, EA

PARTS	Description
SL182A0L	Left Brake,SL18,1/ea, EA
SL182A0R	Right Brake,SL18,1/ea, EA
SL182A10	Seat Belt,SL18,1/ea, EA
SL182A11	Folding Hinge Assm.SL18,1/ea, EA
SL182A12	Support Bar,w/out Hardware,SL18,1/ea, EA
SL182A13	Caster Fork,SL18,1/ea, EA
SL182A14	Fork Bearing Set,SL18,1/ea, EA
SL182A15	Support Chain,SL18,1/ea, EA
SL182A18	Back Uphol.SL18,1/ea, EA