

Ventura 3

HCPCS: K0807-1,K0807



FIGURE A




FIGURE B




FIGURE C



FIGURE D



FIGURE E

- Easily disassembles for transport
- 400 lbs. Weight Capacity
- Includes interchangeable color panels in Red and Blue
- Armrests are padded and adjustable
- Floorboard with 23" of foot space
- Ergonomic Programmable throttle control
- Easy-to-adjust delta tiller
- Height-adjustable swivel seat with fold-down backrest can be adjusted forward or backward
- Flat-free, non-marking tires and anti-tip wheels
- Rear suspension
- Headlight and large, plastic carry basket
- Easy freewheel operation

PRODUCTS	Description
VENTURA318FS	Ventura 3, EA
SPECIFICATIONS	Description
Base Weight	88 lbs.
Batteries (pair)	12V x 33AH
Battery Weight (pair)	28 lbs.
Brakes	Electromagnetic
Charger	4A Offboard
Climbing Angle	7.5°
Dimensions	48" (L) x 25" (W)
Floor to Seat Height	21.5" - 23.5"
Front Wheels	10" x 3" Flat-free
Ground Clearance	2.5"
HCPCS	K0807
Heaviest Piece Weight	60 lbs. (Front Section)
Maximum Range	16 miles
Motor	24V x 350W x 5100 rpm
Rear Wheels	10" x 3" Flat-free
Seat Dimensions	18" (W) x 17.5" (D)
Seat Type	Folding
Seat Weight	26 lbs.
Top Speed	6 mph

SPECIFICATIONS	Description
Total Weight (excl. hand control)	142 lbs.
Turning Radius	53"
Warranty on Batteries	12 Months
Warranty on Electronic Controller/Drive Train Components	24 Months
Warranty on Frame	Lifetime
Weight Capacity	400 lbs
Weight Capacity (13046 13046-HD)	400 lbs.
ACCESSORIES	Description
SF8020	Rear Basket Set, LType, 1/ea, EA
SF8050	Sun Shade for Scooter,1/ea, EA
SF8010	Oxygen Holder Set, L Type,1/ea, EA
SF8030	Crutch Holder Set, SType, 1/ea, EA
PARTS	Description
LRM402208	12AH BatteryCaseW/Batteries-Spitfire1310, EA
VENT-73	Seat Plate, Folding Seat Ventura 3&4, EA
VENT-74	Seat Plate, Capt. Seat Ventura 3&4, 1/ea, EA
VENT-65	Seat Slide Set, Ventura 3 & 4, 2/pkg, PKG
SEC2010N	HeadRestNew18"SunfireEC,1ea, EA
GT8072W32	Hazard Button,GladGT,Ventura,1/ea, EA
GT8072W44	WarningButtonCoverGladGT,Vent,1/ea, EA
VENT-76	Turn Signal Relay, Ventura 3 & 4, 1/ea, EA
S41160	Ignition Switch,Ventura,Odyssey,1/ea, EA
VENT-53	Mirror,Right, Ventura,Panther,1/ea, EA
VENT-54	Mirror,Left,Ventura,Panther,1/ea, EA
GT8072W35	FrontSignalLeftGladGT,Ventura,1/ea, EA
GT8072W36	FrontSignalRightGladGT,Ventura,1/ea, EA
GT8072W53	DirectionBulbGladGT,Ventura,1/ea, EA
VENT-55	Basket Holder, Ventura, 1/ea, EA
VENT-61	Front Decoration, Red, Ventura3&4, 1/ea, EA
VENT-62	Front Decoration, Blue, Ventura3&4, 1/ea, EA
VENT-48	FrontForkAssemblyVentura3,1/ea, EA
VENT-45	FrontShroudVentura3,1/ea, EA
VENT-50	FrontWheelAxle,Ventura 3,1/ea, EA