

## Ventura 4

HCPCS: K0807-1,K0807



FIGURE A




FIGURE B



FIGURE C



FIGURE D



FIGURE E

- Easily disassembles for transport
- 400 lbs. Weight Capacity
- Includes interchangeable color panels in Red and Blue (Figure A)
- Armrests are padded and adjustable
- Ergonomic throttle control
- Height-adjustable swivel seat with fold-down backrest
- Flat-free, non-marking tires
- Large, plastic carry basket
- Seat position can be adjusted forward or backward
- Front and rear suspension
- Programmable controller
- Easy-to-adjust tiller
- Anti-tip wheels
- Headlight
- Delta tiller
- Easy free wheel operation

### PRODUCTS

### Description

VENTURA418FS

Ventura 4, EA

### SPECIFICATIONS

### Description

Base Weight

99 lbs.

Batteries (pair)

12V x 33AH

Battery Weight (pair)

28 lbs.

Brakes

Electromagnetic

Charger

4A Offboard

Climbing Angle

7.5°

Dimensions

48" (L) x 25" (W)

Floor to Seat Height

22" - 24"

Front Wheels

10" x 3" Flat-free

Ground Clearance

2.5"

HCPCS

K0807

Heaviest Piece Weight

70 lbs.

Maximum Range

16 miles

Motor

24V x 350W x 5100 rpm

Rear Wheels

10" x 3" Flat-free

| SPECIFICATIONS                                           | Description         |
|----------------------------------------------------------|---------------------|
| Seat Dimensions                                          | 18" (W) x 17.5" (D) |
| Seat Type                                                | Folding             |
| Seat Weight                                              | 26 lbs.             |
| Top Speed                                                | 6 mph               |
| Total Weight (excl. hand control)                        | 153 lbs.            |
| Turning Radius                                           | 54"                 |
| Warranty on Batteries                                    | 12 Months           |
| Warranty on Electronic Controller/Drive Train Components | 24 Months           |
| Warranty on Frame                                        | Lifetime            |
| Weight Capacity                                          | 400 lbs             |
| Weight Capacity (13046   13046-HD)                       | 400 lbs.            |

| ACCESSORIES | Description                        |
|-------------|------------------------------------|
| SF8030      | Crutch Holder Set, SType, 1/ea, EA |
| SF8010      | Oxygen Holder Set, L Type,1/ea, EA |
| SF8050      | Sun Shade for Scooter,1/ea, EA     |
| SF8020      | Rear Basket Set, LType, 1/ea, EA   |

| PARTS     | Description                                  |
|-----------|----------------------------------------------|
| VENT-86   | Shroud Knob, Ventura, 1/ea, EA               |
| LRM412123 | 21AH Battery,POV, 1/ea, EA                   |
| VENT-15   | ShockAbsorberVentura3&4, EA                  |
| VENT-14   | FrontBumperAssemblyVentura4,1/ea, EA         |
| VENT-44   | Head Light Assembly, Ventura 4, 1/ea, EA     |
| VENT-17   | FrontFrameVentura4,1/ea, EA                  |
| VENT-16   | Anti-Tipper Wheel, Ventura 3 & 4,1/ea, EA    |
| VENT-11   | MiddleShroudVentura3&4,1/ea, EA              |
| VENT-09   | BatteryStrapsVentura,CPN,1/set, PR           |
| VENT-10   | FrontShroudVentura4,1/ea, EA                 |
| VENT-06   | FrontTillerCoverVentura3&4,1/ea, EA          |
| VENT-07   | RearTillerCoverVentura3&4,1/ea, EA           |
| VENT-08   | ControllerVentura3&4,1/ea, EA                |
| VENT-05   | FacepltW/Switches18/20CapSeatVentura3&4, EA  |
| VENT-02   | TillerSwivelAssemblyVentura3&4,1/ea, EA      |
| VENT-04   | FacepltW/Switches18FshSeatVentura3&4/1ea, EA |
| VENT-03   | HandlebarAssemblyVentura3&4,1/ea, EA         |
| VENT-01   | Forward/ReverseAssembly,Ventura3&4,1/ea, EA  |
| VENT-42   | SeatAssemblyVentura3&4,18"CaptSeat,1/e, EA   |
| VENT-43   | SeatAssemblyVentura3&4,20"CaptSeat,1/e, EA   |